

Gas Chromatograph

GC-8A Series



GC-8A Series

Dedicated Single Detector Gas Chromatograph

GC that meets increasing analytical requirements...

GC that can analyze instantaneously...

GC that is tough and useful...

The GC-8A is a gas chromatograph that makes these dreams come true.



GC-8AIT



GC-8AIF

Simple design with emphasis on functionality

Tough and strong diecast frame

Reliability backed by years of experience

Self-diagnostic monitor / thermometer for easy maintenance

Accommodates a range of stable, high-sensitivity detectors (TCD, FID, ECD, FPD)

Easy-to-use digital settings

Comprehensive accessories

Contents

P 3	Column Oven	P 9	GC-8AIE
P 4	Sample Injection Ports	P 10	GC-8APFp
P 5	Detectors	P 11	Accessories
P 7	GC-8AIT		
P 8	GC-8AIF		

Column Oven

Large-capacity, cubic oven

Almost cubic in shape, the column oven opens entirely at the front, making changing columns simple and easy. The column oven's door opens and closes by a simple one-touch operation. The lightweight door moves extremely smoothly.



Fig. 1 Column Oven

Two-stage overhear prevention for protection of columns

If the temperature exceeds the programmed (preset) temperature by 30°C, the heater will be automatically turned off. If the temperature reaches about 420°C, the heater will be automatically turned off.

Temperature control offers excellent stability and response

Excellent temperature control of the column oven is essential for a gas chromatograph. Fig. 2 shows the response and stability when the column oven's temperature setting is changed from 100°C to 101°C. It shows that the oven responds rapidly to a changed temperature setting and is able to immediately achieve stable temperature control.




Fig. 2 Temperature Stability



Fig. 3 Operation of Column Oven Door

Rapid heating and cooling characteristics

Rapid heating and cooling characteristics of the column oven are also important for a gas chromatograph. Fig. 4 shows that heating from 100°C to 350°C requires just 12 minutes, while cooling takes less than one minute and re-stabilizing at the initial temperature takes about three minutes.




Fig. 4 Heating and Cooling Curve

Sample Injection Ports

Easy-to-use on-column sample injection ports

The standard injection ports are on-column type. Since sample solution is directly injected into the head of a glass column, the decomposition and adsorption of sample components are minimized.




Fig. 6 On-Column Sample Injection Port

Capillary column adapter

Adding the optional CLH-800 capillary column holder allows simple mounting of the capillary column for both split and splitless sample injection.



Fig. 7 Quartz Capillary Column Mounted with CLH-800

Detectors

Thermal Conductivity Detector (GC-8AIT, GC-8APT)

Highly stable and sensitive filaments

The rhenium-tungsten filaments, each of which has a resistance of 100 Ω , ensure the highest sensitivity and stability.

Constant current system

The bridge current is kept constant even when the detector temperature is changed. It is not necessary therefore to readjust the current after the detector temperature is stabilized.

Filament protection

If some air should flow into the TCD cell, for example, at the time of column replacement, the bridge current will automatically be turned off and the warning lamp will be lit.




Fig. 8 TCD Cell (GC-8APT)

Flame Ionization Detector (GC-8AIF, GC-8APF)

Cylindrical collector

Combustion gases from the hydrogen flame pass from the nozzle straight through the cylindrical collector to be discharged externally. This reduces contamination of the detector to extremely low levels.

On-detector system

The dead space is eliminated at the connection of the column and the detector. Excellent repeatability is provided even in analysis of trace components.

External igniter

The flame igniter is mounted external to the detector to significantly reduce wire damage due to corrosive combustion components and to make replacement of the igniter filament cheap and simple.




Fig. 9 FID Cell

Electron Capture Detector (GC-8AIE)

Linear ECD with high sensitivity and wide dynamic range

The constant-current ECD maintains a constant current flow in the detector. It offers excellent performance of a linear dynamic range of 10,000 and 0.2 pg (Y-BHC) minimum limit of detection, and greatly restricts peak abnormalities such as negative peaks.

Switchable detector current

The standing current flowing in the detector can be switched in three stages for convenient analysis of samples across a wide concentration range.

Comprehensive contamination countermeasures

The system features a programmed heating system to prevent detector contamination (column oven heating starts only after the detector reaches the set temperature) and a purge system to maintain carrier gas purging and keep the detector in standby status even when the detector is not being used for analysis.

Purge flow line reduces stabilization time

The purge flow line allows purging of the interior of the detector by carrier gas when the detector remains mounted in the gas chromatograph but is not used.

This prevents oxidation by oxygen in the air and ensures stable operation the next time the detector is used.

Flame Photometric Detector (GC-8APFp)

Flame

The FPD cell is designed so that the flame will not be extinguished by a solvent overload.

FID monitor electrode

The FID monitor electrode in the detector permits extraction of the FID signal.

*Requires the FID monitor kit (Special Accessory: P/N 221-23950-91)

Forced cooling system

To ensure a low noise level and long service life, the photomultiplier tube is kept cool by a blower when the detector is maintained at a high temperature.

Complete safety

The photomultiplier tube is protected against damage from exposure to high light levels by a sensing circuit which cuts power to the tube if anomalous conditions are detected.




Fig. 10 ECD Cell




Fig. 11 Chromatogram of Trihalomethane




Fig. 12 FPD Cell




Fig. 13 Analysis of Trace Thiophene in Benzene (S mode)




Fig. 14 Analysis of Organophosphorus Pesticides (P-mode)

GC-8AIT



Column Oven

Temperature range	: Room temperature to 399°C
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^{\circ}\text{C}$
Temperature	: When the temperatures of the column oven and the injection port/detector have reached the set points, the respective READY lamps are lit.
Readout (Monitor)	
Heating speed	: From room temperature to 350°C in 13 minutes.
Cooling speed	: From 350°C to 100°C in 3.5 minutes.
Column	: Stainless steel: 6m x 2 Glass: 5.4m x 2
Overheat Protection	: Two-stage protection: (1) If the temperature exceeds the programmed (preset) temperature by 50°C, the heater will be automatically turned off. (2) If the temperature reaches about 420°C, the heater will be automatically turned off.

Temperature of Injection Port/Detector

Temperature range	: Room temperature to 400°C (10°C steps)
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^{\circ}\text{C}$
Temperature	: When the temperatures of the column oven and the injection port/detector have reached the set points, the respective READY lamps are lit.
Readout (Monitor)	
Overheat Protection	: If the temperature reaches about 420°C, the heater will be automatically turned off.
Injection Port	: On-column injection ports

Flow Control

Carrier gas	: Two pressure regulators Two column inlet pressure gauges
-------------	---

Thermal Conductivity Detector (TCD)

Detector	: Differential type of semi-diffusion flow system Four rhenium-tungsten filaments (100 Ω each)
Maximum temperature	: 400°C
Bridge current	: Constant current system (OFF, 60, 70, 80, 90, 100, 120, 140, 160, 180, 200mA)
Filament protection circuit	: Incorporated

Others

Dimensions	: 440W x 405D x 435Hmm
Weight	: 26.5kg
Power Requirements	: AC100/115V or 220V as ordered. 1500VA max. 50/60Hz.
Others	: (1) Column oven heating starts after the injection port / detector reaches the set temperature. (2) Cannot be used with a preparative system.

GC-8AIF



Column Oven

Temperature range	: Room temperature to 399°C
Temperature control:	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^{\circ}\text{C}$
Heating speed	: From room temperature to 350°C in 13 minutes.
Cooling speed	: From 350°C to 100°C in 3.5 minutes.
Column	: Stainless steel: 6m x 2 Glass: 5.4m x 2
Overheat Protection	: Two-stage protection: (1) If the temperature exceeds the programmed (preset) temperature by 30°C, the heater will be automatically turned off. (2) If the temperature reaches about 420°C, the heater will be automatically turned off.

Temperature of Injection Port/Detector

Temperature range	: Room temperature to 400°C (10°C steps)
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^{\circ}\text{C}$
Overheat Protection	: If the temperature reaches about 420°C, the heater will be automatically turned off.
Injection Port	: On-column injection ports

Temperature Readout

The three temperatures (injection port/detector, column oven, and auxiliary) are selectively indicated on the pyrometer, the selection being made by pushbuttons.

Flow Control

Carrier gas	: Two pressure regulators Two column inlet pressure gauges
Hydrogen	: Two pressure regulators Two pressure gauges
Air	: One pressure regulator One pressure gauge

Flame Ionization Detector (FID)


Detector	: Dual cell with cylindrical collectors
Maximum temperature	: 400°C
Detection limit	: 3×10^{-12} g/s for diphenyl
Dynamic range	: 10^7
Jet	: Mode of quartz

Electrometer

Sensitivity	: 1×10^{-12} A/mV
Noise	: Below 10^{-14} A
Drift	: Below 2×10^{-14} A/h (Constant room temperature)
Input attenuation	: 1, 10, 10^2 , 10^3 , 10^4
Output attenuation	: 1 to 1024 (binary step) and ∞
Linearity	: Wider than 10^5 (fixed input attenuation)
Background suppression	: $\pm 6.4 \times 10^{-10}$ A (NORM mode)
Signal output	: 0 to 1mV for recorder 0 to 1V for data processor

Others

Dimensions	: 440W x 405D x 570Hmm
Weight	: 34.5kg
Power Requirements	: AC100/115V or 220V as ordered. 1500VA max. 50/60Hz.



GC-8AIE



Column Oven

Temperature range	: Room temperature to 399°C
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^{\circ}\text{C}$
Temperature	: When the temperatures of the column oven and the injection port/detector have reached the set points, the respective READY lamps are lit.
Readout (Monitor)	
Heating speed	: From room temperature to 350°C in 13 minutes.
Cooling speed	: From 350°C to 100°C in 3.5 minutes.
Column	: Stainless steel: 12m x 1 Glass: 5.4m x 1
Overheat Protection	: Two-stage protection: (1) If the temperature exceeds the programmed (preset) temperature by 50°C, the heater will be automatically turned off. (2) If the temperature reaches about 420°C, the heater will be automatically turned off.

Temperature of Injection Port/Detector

Temperature range	: Room temperature to 350°C (10°C steps)
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^{\circ}\text{C}$
Temperature	: When the temperatures of the column oven and the injection port/detector have reached the set points, the respective READY lamps are lit.
Readout (Monitor)	
Overheat Protection	: If the temperature reaches about 350°C, the heater will be automatically turned off.
Injection Port	: On-column injection ports

Flow Control

Carrier gas	One pressure regulator One column inlet pressure gauge
-------------	---

Electron Capture Detector (ECD)

Detector	: Coaxial type cell and purge flow line
Radiation source	: ^{63}Ni 370MBq
Detection limit	: 0.2pg/s for γ -BHC
Dynamic range	: 10^4
Maximum temperature	: 350°C
Standing current	: 0.5, 1.0, 2.0mA

Others

Dimensions	: 440W x 405D x 435Hmm
Weight	: 26.0kg
Power Requirements	: AC100/115V or 220V as ordered. 1500VA max. 50/60Hz.
Other	: Column oven heating starts after the injection port / detector reaches the set temperature.

GC-8APFp



Column Oven

Temperature range	: -100°C to 400°C (The optional cryogenic work attachment and cryogenic valve unit are necessary for sub-ambient temperature operation.)
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^\circ\text{C}$
Heating speed	: From room temperature to 350°C in 13 minutes.
Cooling speed	: From 350°C to 100°C in 3.5 minutes.
Column	: Stainless steel: 12m x 1 Glass: 5.4m x 1
Overheat Protection	: Two-stage protection: (1) If the temperature exceeds the programmed (preset) temperature by 30°C, the heater will be automatically turned off. (2) If the temperature reaches about 420°C, the heater will be automatically turned off.

Temperature Programmer/Timer

Initial temperature	: -100°C to 399°C (1°C steps)
Final temperature	: 0°C to 400°C (10°C steps)
Programming rate	: 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 16, 20, 32°C/min.

Temperature of Injection Port/Detector

Temperature range	: Room temperature to 400°C (10°C steps)
Temperature control	: Proportional-integration type of zero cross switching system
Temperature control accuracy	: $\pm 0.1^\circ\text{C}$
Overheat Protection	: If the temperature reaches about 420°C, the heater will be automatically turned off.
Injection Port	: On-column injection ports

Temperature Readout

The three temperatures (column oven, injection port/detector, and auxiliary) are selectively indicated on the pyrometer, the selection being made by pushbuttons.

Flow Control

Carrier gas	: One primary pressure regulator One primary column inlet pressure gauge One differential flow controller One column inlet pressure gauge
Hydrogen	: One pressure regulator One pressure gauge
Air	: One pressure regulator One pressure gauge

Flame Photometric Detector (FPD)

Selective detection of sulfur and phosphorus (by filter interchange)	
Optional filter	: 394nm (S), 526nm (P)
Detection limit	: 5×10^{-11} gS/s for S in thiophene 1.4×10^{-12} gP/s for P in DDVP or parathion
Power supply	: -700V, stabilized
Maximum temperature	: 350°C (air cooling)

Electrometer

Sensitivity	: 1×10^{-10} A/mV
Noise	: Below 10^{-12} A
Drift	: Below 2×10^{-12} A/h (Constant room temperature)
Input attenuation	: 1, 10, 10^2 , 10^3
Output attenuation	: 1 to 1024 (binary step) and ∞
Linearity	: Wider than 10^5 (fixed input attenuation)
Background suppression	: $\pm 6.4 \times 10^{-8}$ A
Signal output	: 0 to 1mV for recorder 0 to 1V for data processor

Others

Dimensions	: 440W x 405D x 650Hmm (Electric control: 440W x 330D x 135Hmm)
Weight	: 35.0kg
Power Requirements	: AC100/115V or 220V as ordered. 1500VA max. 50/60Hz.
Special accessories	: FID monitor and electrometer for monitor (P/N 221-23950-91) Filter for S (P/N 221-00892-01) Filter for P (P/N 221-00897-01)

Accessories

Hydrogen Generator

Compact hydrogen generator by ion exchange membrane water electrolysis. High-purity hydrogen is formed through electrolysis of pure water. The compact design fits perfectly beside a GC-8A gas chromatograph. Requires deionized water with a specific resistance of 5 MΩcm min.

Model	OPGU-2200S
P/N	221-45131-91 (115V) 221-45131-93 (230V)
Max. delivery flow rate	225 mL/min
Max. delivery pressure	20 to 400kPa



Air Compressors

These supply clean (oil mist free) compressed air to support combustion in FID and FPD. Output: 1/4HP
Delivery pressure: 400±50kPa
Maximum delivery flow rate: 20L/min. (NTP)

Model	Oil-less	Oil-less, silent
P/N	221-72380	670-12138
Weight	16kg	44kg
Noise	—	58db



Oil-less type

PPR Series Pressure Regulator for High Purity Gases

Dedicated GC gas pressure regulator that provides a high-purity, high-accuracy supply of carrier gas or hydrogen from the gas cylinder. The air-purge valve on the primary side prevents ingress of air when the cylinder is changed. This ensures rapid baseline recovery after cylinder replacement and extends the column life by restricting oxidative degradation of the solid phase.

Model		P/N
PPR-N ₂	Nitrogen and Argon (Blue)	221-35999-01
PPR-H ₂	Hydrogen (Red)	221-35999-02
PPR-He	Helium (Yellow)	221-35999-03



Accessory/Supply Set (P/N 221-38650-90)

Contains all the items required to operate the gas chromatograph. Major contents: a carrier gas supply pipe, a gas filter, leak check fluid, tools, a mirror for flame checking, a microliter syringe, column tags, quartz beads, color markers, a power outlet, a case, etc.

The separate FID accessories set (P/N 221-38651-90) is also required for a gas chromatograph with FID.



Economy Set of Microsyringes

A range of analytical syringes is available.

- Economy type 10 μL (6 per pack)
- Standard microsyringe 25 to 100 μL
- Gas-sealed microsyringe 25 to 100 μL
- Plunger-guide microsyringe 5 to 25 μL
- Plunger-in-needle microsyringe 0.5 to 5 μL
- Large-capacity gas-tight syringe 1 to 100 mL



Set of Devices for Packing Columns (P/N 221-37387-91)

Contains all the items required to fill or re-pack the column. It allows the column to be easily repacked. Includes detailed instructions and a convenient carrying case.

Major contents: funnels, aspirator, vibrator (electric engraver), trap bottles, column tags, silica wool, column packing, rubber tubes, connectors, etc.



Columns (Tube for packed column)

Length (m)	Stainless Steel Columns	Glass Columns	
	3mm ID, 4mm OD	3.2mm ID, 5mm OD	2.6mm ID, 5mm OD
0.5	201-48705-05	221-22949-05	221-22950-05
1.0	201-48705-10	221-22949-10	221-22950-10
1.5	201-48705-15	221-22949-15	221-22950-15
2.0	201-48705-20	221-22949-20	221-22950-20
2.5	201-48705-25	221-22949-25	221-22950-25
3.0	201-48705-30	221-22949-30	221-22950-30

* Columns exceeding 3 m also available. Consult your Shimadzu representative for details.



Shimadzu Corporation
www.shimadzu.com/an/

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®".
Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

For Research Use Only. Not for use in diagnostic procedures.
The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.