
# MESA SERIES

PRECISION INDUSTRIAL SCALE

## Instruction Manual

Manual de Instrucciones

Bedienungshandbuch

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■ MZ-102

■ MZ-205

## Índice de materias

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# Table of Contents

<b>Introduction</b> .....	1
<b>Important Information</b> .....	1
<b>Installation &amp; Operating Instructions</b> .....	1
AC Adapter	
Battery Operation	
<b>Key Functions / Description</b> .....	2
ON/OFF Key	
CAL/PRINT Key	
“F” Key	
TARE/ZERO Key	
LIGHT Key	
<b>Calibration</b> .....	2
<b>Application Programs</b> .....	3
Weight Unit	
Counting	
Totaling	
<b>Features</b> .....	4
Digital Calibration	
Power Up Test	
Overload	
Auto Off	
RS-232	
<b>Care &amp; Maintenance</b> .....	4
<b>Menu Settings</b> .....	5
<b>Scale Operating Menu</b> .....	5
<b>Interface Description</b> .....	6
<b>Product Specifications</b> .....	7

---

## Introduction

---

Thank you for purchasing a **Denver Instrument Mesa Series Scale**. Please read all of the operating instructions to ensure optimum performance. The unit is pre-calibrated at the factory and is ready to use with the supplied AC adapter.

---

## Important Information

---

- Allow 15 to 20 minutes warm up time before starting calibration routine in order to give internal components sufficient time to stabilize.
- Avoid lengthy exposure to extreme temperatures. Your scale works best when operated within the temperature range listed within the specifications.
- Place scale platform on solid, level surface.
- Keep environment clean, avoid dust, dirt & liquids from getting inside of scale platform and indicator, which can adversely effect and damage the scale.
- Vibrations, air currents and electromagnetic interference may effect the performance of the scale.
- Although the scale is very durable, please prevent excessive overloading and dropping objects on to the platform that may cause permanent damage to the weight cells.

---

## Installation & Operating Instructions

---

Remove and inspect all contents from package.

### Wall Bracket:

If you will be mounting the indicator with the provided bracket we suggest that you install this first. Firmly attach the bracket to the wall using appropriate fastening screws (not included) for your wall type and the weight of the indicator. Slip the indicator into the bracket and attach using the screws included. Note: Do not remove the cable from the base or indicator.

### AC Adapter:


After setting-up the base and the indicator in the desired location plug (included) AC adapter into main outlet (AC adapter volt age dependent upon local region / only use AC adapter provided with scale. Use of another type of adapter may permanently damage the scale and void warranty). Turn unit on by pressing the "ON/OFF" key. Self diagnostic test will automatically run displaying "8.8.8.8." and then to weigh mode with selected unit (Example: "0.000kg"). Scale is ready for use. For additional features and operation please see "Key Functions / Description" and "Features".

### Battery Operation:

One 9 volt Alkaline or 9 volt NICAD rechargeable battery is required (not included with scale). To install battery, open battery cover located on the back of the indicator. Connect battery to the leads accordingly and install within holder and replace cover. It is recommended to remove battery when scale is not in use for extended periods. Also replace low or dead batteries.

Note: Back-lit display is "off" unless activated with "LIGHT" key.

## Key Functions / Description:

### 1. ON/OFF:

Press to turn scale ON or OFF.

### 2. TARE/ZERO:

Pushing TARE/ZERO key “zero’s” display, pushing and holding TARE/ZERO key for more than 3 seconds accesses application menu.

### 3. “F” Key:

Function key for applications (Unit toggling, Counting & Totaling ) and menu scrolling.

### 4. CAL/PRINT:

Press and “hold” CAL/PRINT key (for approximately 3 seconds) and scale will automatically enter calibration mode (displaying one of three weight values). Press and “release” CAL/PRINT key to print displayed information when utilizing RS-232 interface.

### 5. LIGHT :

Turns on backlight during battery operation (Auto Off: One Minute).

## Calibration:

- ① **Weight Calibration:** with the pan cleared and zero stability indicated in display, press and hold the “CAL/PRINT” key for more than 3 seconds to enter the calibration mode.
- ② Depending upon your model a default weight will be displayed in “kg” (note; you have 3 weight values to choose from).
- ③ In order to select additional cal-weight options, press and release the “F” key and the display will toggle between selectable values based on model. (Please see Product Specifications)
- ④ After selection of the desired cal-weight option, press the “CAL/PRINT” key and release to confirm, and the value will flash to begin calibration.
- ⑤ Place corresponding weight onto center of pan. If weight is accepted, the display will stop flashing and show “selected kg” value in display.

- 6 Please remove weight and zero scale.
- 7 Place cal-weight back on to pan and the result should be within +/- 1d.
- 8 Calibration complete.

## Calibration Notes:

- The actual calibration weight needs to be within +/- 5% of the selected weight.
- If the scale is not zeroed before calibration, error "E-01" will be displayed.
- If calibration weight is positively "above" the 5% tolerance, error "E-03" will be displayed.
- If calibration weight is negatively "below" the 5% tolerance, error "E-03" will be displayed.
- Always allow a minimum of thirty seconds for the scale to stabilize. During prolonged periods of use repeat calibration procedure at regular intervals for optimum accuracy.

**Note:** Please do not open or attempt to "internally" calibrate the scale. This will void warranty and may damage scale .

---

## Application Programs:

---

### Weight Units: (Default Application)

- 1 To select "units" press and hold "TARE/ZERO" key for more than 3 seconds to access application menu ("UNIT" will be displayed).
- 2 With "F" key you can scroll through the options (UNIT, COUNT, & TOTAL).
- 3 Select "UNIT" by pressing "CAL/PRINT" key (default unit "grams" or "kg" will be displayed).
- 4 With the "F" key you can scroll through weigh unit options (g, lbs, ozt & kg) depending upon model.
- 5 Select desired weigh unit by pressing "CAL/PRINT" (selected option will be displayed and balance will return to weighing mode).
- 6 Pressing the "F" key will enable you to toggle between selected weight unit and default "kg".
- 7 Selected weigh unit will be retained in memory after balance is turned off.

### Counting:

- 1 To select "Counting" application press and hold "TARE/ZERO" key for more than 3 seconds to access application menu ("UNIT" will be displayed).
- 2 With "F" key you can scroll through the options (UNIT, COUNT, & TOTAL).
- 3 Select "COUNT" by pressing the "CAL/PRINT" key.
- 4 Display will automatically exit menu and return to weighing mode, with a "CNT" displayed in the left hand side of the display & the selected weigh unit displayed (signifying to the user that they are in an application mode).
- 5 Place empty container or box on scale and tare. (If applicable)
- 6 Press "F" key to display reference sample/default "1" (quickly pressing and releasing the "F" key will advance the number 1,2,3, etc... / Holding the "F" key will advance the numbers faster 10, 20, 50, up to 99).
- 7 After selection of the desired number, place corresponding samples on scale and press "CAL/PRINT" to confirm (display will flash and then show "pcs" with corresponding piece count).
- 8 Pressing the "F" key will allow the user to toggle between gross weight, net weight, average piece weight & pcs.

- 9 To exit “Counting” press and hold “TARE/ZERO” key for more than 3 seconds to access application menu and select “UNIT” (default setting).

### **Totaling:**

- 1 To select the application “Totaling” press and hold “TARE/ZERO” key for more than 3 seconds to access application menu.
- 2 With the “F” key you can scroll through the menu options (UNIT, COUNT, & TOTAL).
- 3 Select “TOTAL” by pressing the “CAL/PRINT” Key.
- 4 Display will automatically return to weighing mode with “TTL” displayed (signifying to the end user that they are in application mode).
- 5 Place 1st object on scale and press “F” (value is stored) / remove object and display will return to zero.
- 6 Place 2nd object on scale and press “F” (value is added to 1st) and the “TTL” symbol will flash (indicating total weight stored) / remove object and display will return to zero (“TTL” stops flashing).
- 7 Place 3rd object on scale and press “F” (value will be added and “TTL” will flash) / remove object and display will return to zero (“TTL” stops flashing) process can be repeated up to 99 times.
- 8 Holding down the “F” key for 3 seconds will clear all memory.
- 9 To exit “Totaling”, press and hold “TARE/ZERO” Key for more than 3 seconds to access application menu and select “Unit” (default setting).

---

## **Features:**

---

### **Digital Calibration:**

Allows for convenient user calibration from the keypad with three selectable weight options.

### **Power Up Test:**

When the unit is initially turned on all segments of the display will appear. The segments will be displayed for approx. 3 seconds. After self-test scale will return to zero.

### **Overload:**

When a load or force is applied to the scale platform that exceeds the unit’s capacity “H” will appear in the display indicating an overload has occurred. Immediately removing the object or weight from the unit will allow normal operation to return. If unit does not return to normal operation, damage has occurred.

### **Auto Off:**

During battery operation (in order to conserve battery consumption) the scale will automatically turn off after approximately 4 minutes of inactivity / default setting. (Please see “Menu” options for additional settings).

### **RS-232:**

Standard feature included with scale. (Please see Interface Description)

---

## **Care and Maintenance:**

---

- Avoid lengthy exposure to extreme heat or cold. Please allow the scale to acclimate to room temperature for at least 15 to 20 minutes if exposed to extreme temperatures. (Please see Product Specifications)
- Keep product clean. Remove pan from base and clean accordingly. Do not “wash down” scale or indicator with water or cleaning liquids.
- Please use supplied packaging when transporting scale.
- Remember to remove battery when not in use.

---

**Menu Settings:**

---

**Set up Menu / Configuring Scale**

(Factory set parameters are marked "0")

- ① Turn off scale.
- ② Turn on while holding the TARE/ZERO key.
- ③ Page or level value "1" will be displayed.
- ④ To navigate within a menu level press "F", when you reach the last menu item, the first item is shown again.
- ⑤ To select the next page or level within a group press and release "CAL/PRINT" key
- ⑥ To return to the previous level (or page) press "TARE/ZERO" Key.
- ⑦ To confirm selected parameters press the "CAL/PRINT" key.
- ⑧ "0" indicates the selected menu setting.
- ⑨ To store selected codes and exit menu press and hold "TARE/ZERO" key for 3 seconds or until display returns to weigh mode.

---

**Scale Operating Menu:**

---

**Unit Toggling**

- 141 Unlocked "0"
- 142 Locked

**No. Of Stop Bits**

- 231 1 "0"
- 232 2

**Auto Zero**

- 151 On "0"
- 152 Off

**Manual / Auto Print Mode**

- 311 Manual without stability
- 312 Manual after stability "0"
- 313 Auto without stability
- 314 Auto at stability

**Calibration Setting**

- 171 Calibration on "0"
- 172 Calibration off

**Disable Keys**

(Disables All Keys Except for "TARE/ZERO" Key)

- 411 On
- 412 Off "0"

**Baud Rate**

- 211 1200 Baud "0"
- 212 2400 Baud
- 213 4800 Baud
- 214 9600 Baud

**Battery ON/OFF Mode**

- 511 Turns off after 4 minutes of in-activity "0"
- 512 Always on
- 513 Turns off after 10 minutes of in-activity

**Parity**

- 221 Mark
- 222 Space
- 223 Odd "0"
- 224 Even

**Factory Settings**

- 911 Reset
- 912 Do not reset "0"

## Interface Description

- **Type of Interface:** *Serial*
- **Operating Mode:** *Full Duplex*
- **Transmission Rates:** *1200, 2400, 4800, 9600 Baud*
- **Parity:** *Mark, Space, Odd, Even*
- **Data Output Format:** *16 Characters*
- **Factory settings:** *(Please See Scale Operating Menu)*
- **Pin Assignment:** *Pin2: Tx, Pin3: Rx, Pin5: gnd*

<b>Commands</b>	<b>ASCII</b>	<b>HEX</b>	<b>Remarks</b>
Print	"p" ** CR LF	70 ** 0D 0A	
Tare/Zero	"t" ** CR LF	74 ** 0D 0A	
Change Unit (Can not be used in totaling mode)	"u" ** CR LF	75 ** 0D 0A	
"F" Function Key	"f" ** CR LF	66 ** 0D 0A	
"Unit" Mode	"a" ** CR LF	61 ** 0D 0A	
"Counting" Mode	"b" ** CR LF	62 ** 0D 0A	
"Totaling" Mode	"c" ** CR LF	63 ** 0D 0A	

### **Counting Mode Commands**

Set Number of Pieces	"s" nnn CR LF	73 nnn 0D 0A	= number of pcs
Net Weight	"w" ** CR LF	77 ** 0D 0A	
Gross Weight	"x" ** CR LF	78 ** 0D 0A	
Average Weight	"y" ** CR LF	79 ** 0D 0A	
Pcs	"z" ** CR LF	7a ** 0D 0A	

### **Totaling Mode Commands**

Clear Memory	"1" ** CR LF	6c ** 0D 0A	
--------------	--------------	-------------	--

\*\* : space

### **Output Format**

Position:	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
	-	*	D	D	D	D	D	D	D	D	*	U	U	U	CR	LF

D: Digit or Letter U: Unit CR: Return LF: LineFeed \* : Space

**Product Specifications:**

<b>Model No.</b>	<b>MZ-102</b>	<b>MZ-205</b>
<b>Capacity</b>	100 kg, 220 lbs, 3215 ozt	200 kg, 440 lbs, 6430 ozt
<b>Readability</b>	0.02 kg, 0.05 lbs, 0.5 ozt	0.05 kg, 0.1 lbs, 1 ozt
<b>Weighing Modes</b>	kilograms, pounds, troy ounces	
<b>Tare Range</b>	To capacity by subtraction	
<b>Calibration</b>	5 kg, 20 kg, 50 kg	
<b>Maximum Overload</b>	200 %	
<b>Overload Indication</b>	"H" reading on LCD with audio tone	
<b>Underload Indication</b>	Displays negative readings	
<b>Auto Off</b>	Menu selectable: 4 or 10 minutes no activity / Always on	
<b>Operating Temperature</b>	50 - 104 degrees F (10 - 40 degrees C)	
<b>Low Battery Indication</b>	"LO" on LCD with audio tone	
<b>Power Requirements</b>	Included AC adapter or optional 9 volt battery	
<b>Pan Type</b>	Stainless steel	
<b>Pan Size</b>	403 mm x 520 mm (15.9" x 20.5" inches)	

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**Denver Instrument GmbH