# Mikro-Dismembrator S

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Operating Manual

Mode d'emploi

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## Operating Manual

## **Introductory Note**

This Operating Manual for the homogenizer Mikro-Dismembrator S refers to the design and equipment features at the time given with the Rev.-No. in the bottom line.

Sartorius Stedim Biotech GmbH reserves the right to modify the equipment and this operating manual without prior notice

The information included herein has been carefully compiled. However, this document is not necessarily complete. Continuous development of the device and its equipment can result in incomplete or even erroneous information herein.

Therefore, if you find missing, misleading or wrong information please send us a note.

## **EG Statement of Conformity**

With the declaration of conformity, which is attached or which accompanies the equipment, Sartorius Stedim Biotech GmbH confirms the conformance of the Mikro-Dismembrator S with the specified guidelines (see Point 6.5 EG Statement of Conformity).

## **Guide through this manual**

This manual is organized systematically. Depending on your knowledge about the equipment and its use you can start from the beginning or go directly to the chapter of your interest:

Section 1: Overview on design and function of the equipment and accessories; short reference of possible applications

Section 2: Information on delivery, warranty conditions and installation

Section 3: Operating information, mounting and use of the device

Section 4: Summary of technical data and ordering information

Section 5: Application examples, tips and tricks for successful use

Section 6: Additional specifications and material data.

The Mikro-Dismembrator S is just an example of laboratory equipment products from Sartorius Stedim Biotech GmbH.

Our range of equipment includes homogenizers, shakers and incubator shakers, centrifuges, ultrapure water systems, air samplers and fermenters/bioreactors. You will find further information on our homepage [www.sartorius-stedim.com](http://www.sartorius-stedim.com).

Sartorius Stedim Biotech GmbH  
August-Spindler-Straße 11  
37079 Göttingen, Germany  
Telephone +49.551.308.3118  
Fax +49.551.308.3918  
[info@sartorius-stedim.com](mailto:info@sartorius-stedim.com)  
[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

**Contents****6****Supplement  
21****Introductionary Note****EG Statement of Conformity****Guide through this manual**

|          |  |           |
|----------|--|-----------|
| <b>1</b> | <b>Design and Function</b>   | <b>4</b>  |
| 1.1      | Features and Intended Use  | 4         |
| 1.2      | Equipment and Functional Principles  | 5         |
| 1.2.1    | Basic Device   | 5         |
| 1.2.2    | Shaking Flasks and Special Containers  | 6         |
| 1.2.3    | Grinding Balls   | 7         |
| <b>2</b> | <b>Installation and Start-up</b>   | <b>8</b>  |
| 2.1      | Delivery Information   | 8         |
| 2.1.1    | Unpacking and Checking Completeness  | 8         |
| 2.1.2    | Warranty Regulations   | 8         |
| 2.1.3    | Service  | 8         |
| 2.2      | Installation of the Mikro-Dismembrator S   | 9         |
| <b>3</b> | <b>Operating Information</b>   | <b>9</b>  |
| 3.1      | Preparing the Samples, overview  | 9         |
| 3.2      | Mounting of Shaking Flasks and Containers  | 10        |
| 3.2.1    | Preparing and Mounting of Shaking Flasks   | 10        |
| 3.2.1.1  | Filling the shaking flasks   | 10        |
| 3.2.1.2  | Special Notes  | 10        |
| 3.2.2    | Mounting the Holder for Shaking Flasks   | 11        |
| 3.2.3    | Using Disposable Test Tubes  | 11        |
| 3.3      | Homogenisation of samples  | 12        |
| 3.3.1    | Switching-on the Mikro-Dismembrator S  | 12        |
| 3.3.2    | Setting the Shaking Frequency  | 12        |
| 3.3.3    | Setting the Operation Time   | 12        |
| 3.3.4    | Starting the homogenization  | 13        |
| 3.4      | Maintenance and Cleaning   | 14        |
| 3.4.1    | Cleaning   | 14        |
| 3.4.2    | Maintenance  | 14        |
| <b>4</b> | <b>Technical Data, Ordering Information</b>  | <b>15</b> |
| 4.1      | Technical Data   | 15        |
| 4.1.1    | Drive Unit   | 15        |
| 4.1.2    | Housing  | 15        |
| 4.1.3    | Ordering Information   | 15        |
| 4.2      | Shaking Containers and Accessories   | 16        |
| 4.2.1    | Standard Shaking Flasks  | 16        |
| 4.2.2    | Holder, special Shaking Containers and Special Accessories   | 17        |
| 4.2.3    | Grinding Balls and Process Materials   | 17        |
| <b>5</b> | <b>Applications</b>  | <b>18</b> |
| 5.1      | Application Information Overview   | 18        |
| 5.2      | Use of the Mikro-Dismembrator S in Trace Analysis and Research of Contaminations   | 19        |
| 5.2.1    | Trace Element Analysis via Homogenization by "Brittle Fracture Technique"  | 19        |
| 5.2.2    | Trace - Analysis of Contaminations by Toxic Substances using the Brittle Fracture Technique                              | 20        |
| 5.2.3    | Classification of the Elements, which could be Detected in Human Tissue after Homogenization with a Mikro-Dismembrator S | 20        |

|       |  |    |
|-------|--|----|
| 6.1   | Recommended Combinations of Shaking Flasks and Grinding Balls        | 21 |
| 6.2   | Basic Analysis of Materials for Shaking Flasks and Grinding Balls    | 21 |
| 6.2.1 | Materials and Characteristics  | 21 |
| 6.2.2 | Manufacturer's Data for Grinding Balls made of Glass ("Dragonit 30") | 22 |
| 6.3   | Bibliographical Data   | 22 |
| 6.4   | Safety Instructions Homogenizers                                     | 23 |
| 6.5   | EG Statement of Conformity   | 23 |
| 6.6   | Declaration of Decontamination                                       | 23 |
| 6.7   | Information and Instructions on Disposal and Repairs                 | 23 |
| 6.8   | Dimension sheet Mikro-Dismembrator S                                 | 23 |
| 6.9   | Prospectus Mikro-Dismembrator  | 23 |

# 1 Design and Function

## 1.1 Features and Intended Use

The Mikro-Dismembrator S is a laboratory scale, vibration grinding mill. It can be used for efficient fine grinding, homogenizing and mixing of small sample volumes, for instance. It can be used for any application which requires efficient and reproducible grinding and mixing of a sample, is suitable for the Mikro-Dismembrator S. The Mikro-Dismembrator S allows shaking frequencies of up to 3.000 <sup>1</sup>/min.

Applications include contamination-free preparation of samples for trace element analysis, the disintegration, homogenizing and effective mixing of organic and inorganic materials, or the grinding of samples for isolation of DNA, RNA and proteins, for instance. Dry or wet milling and even treatment of frozen samples is possible. Depending on the material characteristics, the final fineness of the sample can be about 1 µm. Samples disrupted with the Mikro-Dismembrator S may also be used for isolating DNA, RNA or proteins.

The Mikro-Dismembrators S succeeds the Mikro-Dismembrator II, which has proved its worth for many years in daily laboratory work. For this model a broad range of applications has been reported which are transferable to the Mikro-Dismembrator S. Examples for applications include, e.g. <sup>1-1)</sup>.

- Grinding of teeth, bones, frozen skin and similar materials for tissue analysis (histology, toxicology, forensic medicine).

Disruption of parts of plants (leaves, needle-leaves, small fruits, seeds or grains) for trace element analysis in environmental research).

- Disintegration of soft tissue (even when frozen) for histology.
- Disruption of microorganisms, yeast cells or tissue cell cultures (in microbiology, biotechnology, molecular biology, proteome and genome research).
- Grinding of pigments (dye industry).
- Grinding of other tough materials (leather, etc., even when frozen).

Available accessories are shaking flasks and containers of different size and material (stainless steel or PTFE, for instance) and holders for disposable test tubes. Samples will be ground (disintegrated) using grinding balls of different size and material, such as PTFE-coated steel balls, stainless steel balls, balls made of high purity quartz agate (Brazilian agate) or tungsten carbide or glass beads together with the sample, for instance. Effective mixing of samples is possible if they are treated without addition of balls or beads.

---

<sup>1-1)</sup> See additional information on applications of the Mikro-Dismembrator S, practical experiences and typical results in section 5 of this manual. More information is available on request.

## 1.2 Equipment and Functional Principles

### 1.2.1 Basic Device

The Mikro-Dismembrator S can be operated with a selectable shaking frequency for a definable time. Continuous operation is also possible. An eccentric wheel converts the rotation of the brushless drive into reciprocation of the connector with the shaking flask.

- At a given, constant shaking amplitude (16 mm) the shaking frequency can be varied. For the Mikro-Dismembrator S the range is within 100 ... 2.600  $\frac{1}{\text{min}}$ . in continuous mode and within 100 ... 3.000  $\frac{1}{\text{min}}$ . for fixed operation times (timer mode). The shaking frequency is controlled. Therefore, the set value is kept constant during extended runs and with varying loads, ensuring precise, reproducible experimental conditions. The user can select short-term (timer mode) or continuous operation.
- The integrated timer permits adjustable operating times in the range of 00:01 ... 98:59 [Min:Sec], continuous operation is possible with time setting of 99:59 [Min:Sec].

For disintegration, the sample is filled into the shaking flask together with the grinding balls and process materials, if necessary. Due to the reciprocal movement the components are accelerated and will bounce against the inner wall of the flask and against each other. This will disrupt the sample and, at the same time, provide a homogeneous mixing of the sample.

The impact effect depends on the specific weight and skin hardness of the grinding balls, for instance, and changes with the material characteristics. Increasing shaking frequency and higher specific weight of the grinding ball both enhance the input. The effectiveness of the disruption and homogenization of the sample can be controlled by adding certain process materials, such as quartz sand, chemical process materials or dispersing agents. Dry samples, which are treated without grinding balls, will be intensively mixed.

Fig. 1: Mikro-Dismembrator S

- (1) Shaking flask
- (2) Holder for flask and accessories
- (3) Covering with connector
- (4) Key pad
- (5) Digital display
- (6) Clamp for locating the cover

### 1.2.2 Shaking Flasks and Special Containers

The shaking flasks are made of stainless steel or high-purity PTFE, having different volumes:

1. Shaking flasks made of high purity PTFE, with volumes of 3, 5, 7 and 20 ml, available with caps or screw plugs.
  2. Shaking flasks made of stainless steel 1.4301<sup>1-1)</sup>, with volumes of 3, 5 and 7 ml.
- For details refer to Fig. 2 in section 3 of this manual.
  - A detailed of available shaking flasks is included in section 4, "Technical Data, Ordering Information".

Other shaking flasks, devices and corresponding accessories for treatment of samples with a Mikro-Dismembrator S are:

1. Special holder for 20 ml shaking flasks, cat.-no. 8531897. This holder will be used for the 20 ml PTFE flask and for the special container for disposable test tubes.
  2. Container for disposable test tubes, cat.-no. 8531889. It can hold 3 disposable test tubes. For use only with holder cat.-no. 8531897. Suitable test tubes are disposable test tubes made of plastic for sample volumes of about 2.2 ml, for instance.
  3. Container for 4 cryo test tubes, cat.-no. 8531960. For use only with holder cat.-no. 8531897.
  4. Container for 4 reaction vials 2.0 ml with standing rim cat.-no. 8531944 (such as Bio-Rad Platelia) Holder cat.-no. 8531897 required for using this container.
- For details refer to fig. 4 in section 3 of this manual.
  - More detailed information is included in section 4, "Technical Data, Ordering Information".

The flasks and containers are designed according to the requirements of the assays which follow the sample treatment in the homogenizer. Such requirements are the volume and design of the sample, as required for spectroscopic or radiographic assays (standard volumes of 5 ml) or the volumes for trace element analysis of biological material (20 ml), for instance.

The homogenization method and the selection of suitable accessories for the Mikro-Dismembrator S is determined by the assay method applied after the sample treatment, by the characteristics of the raw sample (consistency, abrasiveness) and the characteristics of the treated sample (consistency, fineness or homogeneity).

Usually the user has to find out empirically which combination of shaking flasks, grinding balls, process materials (see information in the next chapter) and operating parameters, will best suit the requirements.

The operating information below and the application examples in section 5 and the literature references in the supplement will give some examples. Information for specific applications will be available on request.

### 1.2.3

---

<sup>1-1)</sup> According to German Steel Standards

## Grinding Balls

The characteristics of the grinding balls, such as their size, abrasiveness and (specific) weight, will essentially determine the final fineness, homogeneity or compactness of the treated sample and the time required for disintegration. Furthermore suitable grinding balls are selected according to the material characteristics of the raw sample (hardness, compactness or toughness, chemical composition, etc.), the type of assay to be applied afterwards and the required fineness, homogeneity or compactness of treated sample. Process materials like silica sand, dispersants or buffers can be added to improve disintegration and homogenization of the sample.

Available grinding balls and process materials are, for instance:

1. Glass beads of  $\emptyset$  in the range of about 0.10 ... 1 mm.
  2. Balls made of quartz agate (Brazilian agate), with a diameter of about 10 mm at a weight of 2 g.
  3. PTFE balls (PTFE coated steel balls with steel core), with a diameter of about 10 mm at a weight 2 g.
  4. Stainless steel balls,  $\emptyset$  about 3 ... 10 mm, weight range within 0.11 ... 4.1 g.
  5. Balls made of tungsten carbide,  $\emptyset$  about 1 ... 9 mm, weight range within 0.07 ... 5.63 g.
- Please refer to more detailed information about the available grinding balls and process materials in section 4. Specific data (material characteristics, chemical composition, etc.) are included in the supplement.

The grinding balls and process materials must not affect the subsequent analysis of the samples. For instance, rubbed-off material of the grinding balls, of the shaking flasks or of the process materials, respectively, must not distort the results of trace element analysis. Therefore some general recommendations for use of grinding balls and process materials should be heeded:

- When soft materials are to be treated and prepared for trace element analysis, we recommend to use PTFE balls with steel core (made of high-purity PTFE) or balls made of quartz agate (Brazilian agate). The latter show low abrasion and can be used for abrasive samples.
- For most of the common applications, like disintegration of semi-solid or solid samples, the stainless steel balls may suffice.
- Tungsten carbide balls, which offer a high specific weight, can be recommended for fast and efficient disruption of solid (or highly abrasive) samples.

- Except for trace element analysis, glass beads or silica sand can be added to the samples to improve the treatment. Furthermore, Tris-buffer can be added to the shaking flask.
- Tough materials (rubber, plastics, etc.) can be frozen in liquid nitrogen to optimize disintegration.
- Likewise biological materials, such as tissue samples at yeast cells, can be quickly disrupted in the Mikro-Dismembrator S after they have been frozen in liquid nitrogen.
- Disintegrated particles tend to conglomerate due to their adhesion interaction. This can be prevented by adding a dispersion material or liquid to the sample and wet milling.

## 2 Installation and Start-up

### 2.1 Delivery Information

#### 2.1.1 Unpacking and Checking Completeness

- The Mikro-Dismembrator S is delivered after a thorough functional test.

Should, nevertheless, some part of the delivery be missing or show any damage, these hints help you to settle your claims.

1. After receipt carefully remove the device and the accessories from their package. Check whether the equipment and of the set of accessories is complete, according to your order. Carefully check all parts for damages caused by transport.

Prior to commissioning, set the voltage selector switch (fig. 6, pos. 2) to the line voltage in your laboratory. The voltage selector switch is located on the rear side of the device. If the Micro-Dismembrator S cannot be set to your line voltage, do not commission the device.

Connecting conditions are:

| Unit                 | Cat.-no. | Power supply                               |
|----------------------|----------|--|
| Mikro-Dismembrator S | 8531609  | 230 V, 50 Hz, 0.8 A<br>115 V, 60 Hz, 1.6 A |

2. After installation at the working place and assembly and connection of all parts, prepare and enable operation as shown in section 3. Check functioning of the device using all possible operation modes. In case of a wrong power version, of malfunctioning of the unit or if parts are damaged or are missing you immediately should return your claim note to your dealer, or to your representative of Sartorius Stedim Biotech GmbH.

- Transit damages should be reported to the respective carrier, railway or postal authority and you should demand for a damage report.

#### 2.1.2 Warranty Regulations

- Sartorius Stedim Biotech GmbH offers a warranty against malfunction of the equipment, defective material and workmanship according to the "General Terms and Conditions of Business" and unless other terms were agreed upon in writing. Date of reference is the date of delivery.
- The warranty will lapse if any modifications or unauthorized repairs are performed on the unit by the owner or third parties.

- Warranty is not granted for consumables and disposables and defective parts due to wear and tear. Nor is warranty granted for malfunctions or defects of the equipment caused by incorrect operation, improper handling, or normal wear and tear.

With aggressive media you must check the resistance of all parts.

- Warranty claims must be confirmed with a respective delivery note or receipt. They expire if the owner or non-authorized third parties modify or repair the device.

#### 2.1.3 Service

- If malfunctions or damages occur, please contact your dealer or your local Sartorius AG representative or directly:

Sartorius AG  
Servicezentrum Nord  
Weender Landstraße 94-108  
37075 Göttingen, Germany  
Telephone +49.551.308.3729 / 3740 / 3741  
Telefax +49.551.308.3730

- You may return defective equipment to Sartorius AG. Sartorius AG or its duly employed representative will perform any repairs according to the actual Conditions for Inspection and Repair, which are available on request.

Those apparatus sent to us for repair have to be in hygienic perfect condition, clean and well packed. If parts were polluted by media and media components, they have to be cleaned, decontaminated, disinfected or sterilised, according to the valid safety rules of this particular application, e.g. for chemical or biological safety.

The sender has to prove that he observed all safety guidelines. For this purpose you have to fill in the enclosed decontamination declaration and attach it to the apparatus. No apparatus will be repaired without declaration of decontamination or description of the measures taken.

Transport damages or subsequent cleaning or disinfection, if required, will be borne by the sender.

## 2.2 Installation of the Mikro-Dismembrator S

The Mikro-Dismembrator S can be operated on standard laboratory benches.

Please note:

- Shaking movements cause the device to vibrate. The attenuation by suspension of the shake drive and the support bases cannot completely prevent the vibrations from being transmitted to the assembly base.
  - Particularly in the case of very high vibration frequencies or self-resonance of the housing and the assembly base, the Mikro-Dismembrator S may travel across the table.
1. Place the Mikro-Dismembrator S on a solid surface, able to hold the weight of the equipment and to compensate any vibrations. Make sure that other equipment in the vicinity cannot be impaired by the operation of the Mikro-Dismembrator S.
    - Dimensions: W x H x L = 297 x 259 x 205 mm
    - Weight: 19 kg
  2. To improve stable standing of the Mikro-Dismembrator S during operation you may place the unit on to the rubber pad delivered with unit.

The reciprocating action of the connector with the shaking flask transmits vibrations to the unit and to the laboratory bench. These vibrations cannot be fully compensated by the construction of the drive system and the rubber stands. Especially at high shaking frequencies or at self resonance of either the homogenizer or of the laboratory bench, resonance effects may force the homogenizer to move over the working place.

## 3 Operating Information

### 3.1 Preparing the Samples, overview

1. Determine a sample volume depending on the requirements of the subsequent assay and the method of processing the sample (grinding, rough disruption, homogenization or mixing, etc.).
  2. Select a suitable shaking flask and accessories required (grinding balls, etc.).
    - For grinding of soft material you can use PTFE or steel flasks and balls made of PTFE or quartz agate.
      - PTFE balls offer low abrasion resistance and cannot be used for abrasive samples.
      - Balls made of quartz agate have a high abrasion resistance.
      - For the preparation for the trace element analysis, PTFE-covered balls or balls made of quartz agate (Brazilian agate) are recommended.
3. Efficient grinding and mixing requires a sufficient dead volume. Consider their size of the grinding balls and the corresponding volume. Recommended percentage portion of a sample volume is 30 ... 40 % of the filling volume.
    - Optimum percentage volumes of sample material, the grinding balls and the dead volume depend on the sample material. This has to be determined experimentally.
  4. If required, add the appropriate buffer solution for wet grinding or silica sand or a dispersing agent. These process materials can speed-up the homogenization or prevent coagulation of the treated sample.
  5. For tough substances (leather, rubber, plastics, etc.), you can freeze-dry the sample material or freeze it in the shaking flask. For this you may dip the sample and/or shaking flask into liquid nitrogen, for instance ("brittle fracture technique", BFT).
    - For grinding of frozen or semi-solid or solid substances you can use stainless steel flasks and steel balls. Tungsten carbide balls have a high specific weight and are very resistant against abrasion (this feature is recommended for homogenization of abrasive materials).
    - Tungsten carbide balls are not recommended for the isolation of RNA
    - For suitable samples (fresh tissues, etc.) you can use disposable test tubes and glass beads. Glass beads show average abrasion. They are, however, unsuitable for treating frozen samples or for preparing samples for trace element analysis.
- For specific or more detailed information refer to the examples for applications in chapter 5 as well as to material data in chapter 6.

## 3.2 Mounting of Shaking Flasks and Containers

### 3.2.1 Preparing and Mounting of Shaking Flasks

- This information refers to the standard shaking flasks made of PTFE or of stainless steel.

#### 3.2.1.1 Filling the shaking flasks


Fig. 2: Standard shaking flask

Upper fig.: PTFE-flasks with screw cap

Lower fig.: Stainless steel flasks

- (1) Cap of shaking flasks (screw cap for PTFE-flasks)
- (2) Washer (only with steel flasks)
- (3) Beaker for sample

Detailed information about the dimensions  $D_i$ ,  $H_i$  and  $H_a$  of the shaking flask are given in the supplement.

1. Remove the cap of the shaking flask selected. Check the washer, if applicable, it should not be damaged.
2. Fill in the sample. Add grinding balls and process materials, as required.
  - For selecting grinding balls please refer to section 1, section 4, the application examples in section 5 or the material data in the supplement.
    - For wet milling add a suitable buffer, dilutant or dispersing additive.
    - You may add glass beads or silica sand, for instance. Note that glass beads or silica sand cannot be used for trace element analysis. Due to their abrasiveness and contents of various trace elements of varying composition, rubbed-off material may interfere with the assay.
  - Allow for the dead volume required.
    - When using screw-cap PTFE containers, screw on the screw-cap.
    - When using steel containers, place the seal on the bottom part and mount the cap.
    - Assemble PTFE containers and their caps.
3. Place the shaking flask into the holder and start the operation.

#### 3.2.1.2 Special Notes

- When the sample is to be refrigerated prior to the disintegration place the shaking flask with the sample for a few seconds to minutes into liquid nitrogen, for instance. Then mount it into the Mikro-Dismembrator S and start operation.

Risk of injury! Deep-frozen flasks must not be handled with bare hands. Use protective gloves and pliers.

### 3.2.2 Mounting the Holder for Shaking Flasks


Fig. 3: Mounting the holder

- (1) Holder for shaking flasks
- (2) Socket head screw
- (3) Cover of connector
- (4) Connector
- (5) Shaking flask
- (6) Locating hole
- (7) Locating disc
- (8) Head clamp bolt
- (9) Locking nut or thumb screw

- For PTFE- and stainless steel flasks of 3 ... 7 ml you can use the holder delivered with the homogenizer.
  - For the PTFE flask of 20 ml, for the container for disposable test tubes, for the reaction vials with standing rim and for the cryo tubes, you will need to use the holder cat.-no. 8531897.
1. Unscrew socket head screw of the connector (4) and insert the holder for shaking flasks (1). Lock holder with socket head screw (2).
  2. Unscrew the locking nut or thumb screw (9) of the holder (1) and turn the head clamp (8) bolt anticlockwise until the flask (5) or container can be inserted.
  3. Insert the shaking flask into the holder, with its bottom part directed to the locating hole of the holder. Turn head clamp bolt clockwise to locate the flask. Tighten the bolt hand-tight.
  4. Lock the head clamp screw with the locking nut or thumb screw

### 3.2.3

### Using Disposable Test Tubes


Fig. 4: Shaking flask with disposable test tubes

- (1) Disposable test tubes
- (2) Receptacle for test tubes
- (3) Cap

- For disposable test tubes, the shaking container cat.-no. 8531889 will be used together with the special holder, cat.-no. 8531897
1. Mount special holder into the connector as shown for the holder of standard flasks, see 3.2.2.
  2. Fill sample and process materials into the disposable test tubes. Close test tubes. Refrigerate, if necessary.
  3. Place test tubes into the container. Close container with cap.
  4. Mount container into the holder and start the unit, see 3.2.2.
- If the container for cryo test tubes cat.-no. 8532001 or cat.-no. 8531960 or for reaction vials 2.0 ml with standing rim cat.-no. 8531944 are to be applied you can proceed in the same way.

### 3.3 Homogenisation of samples

#### 3.3.1 Switching-on the Mikro-Dismembrator S

Fig. 5: Operating panel of the Mikro-Dismembrator S with key pad

1. Switch-on the Mikro-Dismembrator S with the mains switch on the rear of the unit. The homogenizer starts with the operation mode, selected when the unit was switched-off the last time, i.e. "Timer Mode" or "Frequency Mode". The modes can be distinguished by the data shown on the display. The "Timer-Mode" displays the time [TIME 00:00 Min:Sec], the "Frequency-Mode" displays [RPM X.XXX]. You can switch-over to the other mode with either the key [  $\uparrow$  ] or [  $\downarrow$  ].
2. Enter the operation time and shaking frequency via the numeric keys. See fig. 5 and the following paragraphs.
3. Close the cover and the clamps which fix the cover (see fig. 1, pos. 6).

The Mikro-Dismembrator S can not be started, when the cover is open. A cut-off switch opens when the cover is open and switches-off the unit, when running, or prevents it from being started, respectively.

The safety switch must not be disabled!

#### 3.3.2

#### Setting the Shaking Frequency

1. Switch into "Frequency-Mode" with key [  $\uparrow$  ] or [  $\downarrow$  ] if the unit is in "Timer-Mode".
2. Enter the shaking frequency required via the numerical keys. The value range is:

Mikro-Dismembrator S

100...2.600 <sup>1</sup>/min  
continuous mode

100...3.000 <sup>1</sup>/min  
fixed operation times in  
timer mode

3. Confirm the entry with the [ENTER] key. Wrong inputs can be easily corrected. Press deletion key [ C ] to change the actual entry. The cursor starts with the last (right hand) digit of the display, ready for a new entry. Confirm new entry with the [ENTER] key.

#### 3.3.3 Setting the Operation Time

1. Switch-over to "Timer-Mode" with either key [  $\uparrow$  ] or [  $\downarrow$  ] if the Mikro-Dismembrator S is in "Frequency-Mode".
  2. Enter the operation time in [Min:Sec] via the numeric key pad. The timer value range is 00:01 ... 98:59 [Min:Sec]. Press the [ENTER] key to confirm the value time. For a quick entry of the time, you may confirm after input of [Min].
  3. Wrong inputs can be easily corrected. Press deletion key [ C ] to change the actual entry. The cursor starts with the first (left hand) digit of the display, ready for a new entry of [Min]. Confirm new entry with the [ENTER] key.
- The Mikro-Dismembrator S remains inactive at an input of 00:00 [Min:Sec]. Illogical inputs will be ignored, such as 00:70 [Min:Sec], for instance. The input of 99:59 [Min:Sec] switches-over to "Timer-Mode" for continuous operation.

#### 3.3.4

### Starting the homogenization

1. Close the cover of the unit. It provides noise reduction and protects the user from being injured by shaking flasks, which may be flung-off the holder if not mounted properly, especially at high shaking frequencies.
2. Close the cover and the clamps which fix the cover (see fig. 1, pos. 6).

Only use the Mikro-Dismembrator S with the cover closed. Opening the cover activates a safety switch, which interrupts the operation.

The safety switch must not be disabled!

3. Press the [START] key to start operation with the preset values. The unit runs for the preset time or continuously. Pressing the [STOP] key will interrupt the operation. New values will only be accepted after pressing the [ENTER] key.
4. You will see the remaining running time on the display counting down to zero, when the unit is switched to "Timer-Mode". When switched to the "Frequency-Mode" you will see the actual shaking frequency on the display. You can switch-over between the two modes with either key [ ↑ ] or [ ↓ ].
5. The Mikro-Dismembrator S starts when you press the [ENTER] key. For changing an active shaking frequency or time value, press the [STOP] key, enter the new value, confirm with the [ENTER] key and press the [START] key to restart the unit.

### 3.4 Maintenance and Cleaning

#### 3.4.1 Cleaning

1. Dirty surfaces of the unit, the cover or dirty holders, shaking flasks, containers or grinding balls can be cleaned with usual household detergents or alcohol.

Avoid aggressive cleaning agents.  
For example, chloric agents can cause corrosion

- After cleaning, carefully flush with (demineralized) water. Residues of cleaning agents may interfere with the analysis of samples. Especially in trace element analysis it is recommended to consider blank tests. These tests should include any accessory and process material together with a pure substance or agent or a well-known test material.
- 2. Shaking flasks and grinding balls are subject to wear and tear, depending on the characteristics of the sample material. Check for the degree of abrasion and replace the parts regularly. Note the ordering information on spare parts and consumables in section 4.

#### 3.4.2

### Maintenance

- No maintenance measurements are required for the Mikro-Dismembrator S. If the unit does not start or if malfunctions occur, you can do the following:

Fig. 6: Mains switch, voltage selector and fuses

- (1) Mains switch
- (2) Voltage selector
- (3) Fuses
- (4) Mains socket
- (5) Type plate

1. Check the mains supply, the voltage selector and the proper connection of the mains cable. Also check the fuses (see fuse module on the rear of the module).
2. Replace defective fuses.
3. Check setting of the operation modes and values. Check whether the cover is closed properly.

The Mikro-Dismembrator S cannot start when the cover is open.

4. If all settings are correct and proper handling is ensured and the malfunctioning continuous, please contact your service representative of Sartorius AG.

## 4 Technical Data, Ordering Information

### 4.1 Technical Data

#### 4.1.1 Drive Unit

|                          |   |
|--------------------------|---|
| Shaking amplitude        | 16 mm (constant)  |
| Shaking frequency        | Mikro-Dismembrator S<br>100 ... 2.600 1/min (Continuous-Mode)<br>100 ... 3.000 1/min (Timer-Mode) |
| Timer                    | 00:01...98:59 [Min:Sek]<br>continuous operation when set to 99:59 [Min:Sec]                       |
| Interference suppression | Class N   |

#### 4.1.2 Housing

|                      |   |
|----------------------|---|
| Dimensions W x H x D | 297 x 259 x 205 mm  |
| Weight               | 19 kg   |
| Material             | sheet steel cover with acid-proof coat of lacquer                         |
| Mains supply         | Separate cable, about 3 m,<br>for grounded socket, with 3-pin female plug |

#### 4.1.3 Ordering Information

|         |   |
|---------|---|
| 8531609 | Mikro-Dismembrator S<br>- Version 230 / 115 V, 50 – 60 Hz, switchable<br>- standard flask holder for shaking flasks of 3 ... 7 ml |
|---------|---|

## 4.2 Shaking Containers and Accessories

### 4.2.1 Standard Shaking Flasks


Fig. 7: Standard shaking flask for grinding of samples

Upper fig.: PTFE-flasks with screw cap

Lower fig.: Stainless steel flasks

Components:

- (1) Cap of shaking flask (screw cap for PTFE-flasks)
- (2) Washer (only steel flasks)
- (3) Beaker for sample

| Cat.-no. | Material               | Dimensions     |                |                |             |
|----------|------------------------|----------------|----------------|----------------|-------------|
|          |                        | H <sub>a</sub> | H <sub>i</sub> | D <sub>i</sub> | volume [ml] |
| 8531943  | PTFE                   | 37             | 27             | 13.5           | 3           |
| 8531935  | PTFE                   | 42             | 32             | 16             | 5           |
| 8531927  | PTFE                   | 45             | 35             | 17.5           | 7           |
| 8531951  | PTFE                   | 60             | 50             | 25             | 20          |
| 8531838  | PTFE                   | 37             | 27             | 13.5           | 3           |
| 8531846  | PTFE                   | 42             | 32             | 16             | 5           |
| 8531854  | PTFE                   | 45             | 35             | 17.5           | 7           |
| 8531862  | PTFE                   | 60             | 50             | 25             | 20          |
| 8531803  | stainless steel 1.4301 | 37             | 27             | 13.5           | 3           |
| 8531811  | stainless steel 1.4301 | 38             | 32             | 16             | 5           |
| 8531820  | stainless steel 1.4301 | 41             | 35             | 17.5           | 7           |

#### 4.2.2 Holder, special Shaking Containers and Special Accessories

| Cat.-no. | Item, features   |
|----------|--|
| 8531897  | Holder for 20 ml shaking flasks and for container cat.-no. 8531889   |
| 8531889  | Container for disposable test tubes<br>- container with borings of 10.8 x 37 mm<br>- for 3 disposable test tubes of 2.5 ml volume)<br>- requires holder cat.-no. 8531897 |
| 8531960  | Container for 4 cryo test tubes<br>- requires holder cat.-no. 8531897  |
| 8531944  | Container for 4 reaction vials 2,0 ml with standing rim<br>- such as Bio-Rad Platelia tests, for instance<br>- requires holder cat.-no. 8531897                          |
| 8532001  | Container for 4 cryo test tubes (Nalgene 5011 – 0012)  |
| 38240343 | Holder for 3 ... 7 ml shaking flasks   |

#### 4.2.3 Grinding Balls and Process Materials

| Cat.-no. | Item, features  |
|----------|---|
| 8541400  | Glass beads, Ø about 0,10 ... 0,11 mm, bottle of about 570 ml                             |
| 8541507  | Glass beads, Ø about 0,17 ... 0,18 mm, bottle of about 570 ml                             |
| 8541604  | Glass beads, Ø about 0,25 ... 0,30 mm, bottle of about 570 ml                             |
| 8541701  | Glass beads, Ø about 0,40 ... 0,60 mm, bottle of about 570 ml                             |
| 8541809  | Glass beads, Ø about 1 mm, bottle of about 570 ml   |
| 8547505  | Balls made of quartz agate (Brazilian agate),<br>Ø about 10 mm, weight 2 g, set of 10 pcs |
| 8547602  | PTFE balls with steel core, Ø about 10 mm, weight 2 g, set of 2 pcs                       |
| 8546606  | Chromium steel balls, Ø about 3 mm, weight 0.11 g, set of 100 pcs                         |
| 8546703  | Chromium steel balls, Ø about 5 mm, weight 0.51 g, set of 100 pcs                         |
| 8546916  | Chromium steel balls, Ø about 9 mm, weight 3 g, set of 10 pcs                             |
| 8546800  | Chromium steel balls, Ø about 10 mm, weight 4.1 g, set of 10 pcs                          |
| 8547009  | Tungsten carbide ball, Ø about 1 mm, weight 0.07 g, 1 piece                               |
| 8547106  | Tungsten carbide ball, Ø about 3 mm, weight 0.21 g, 1 piece                               |
| 8547203  | Tungsten carbide ball, Ø about 5 mm, weight 0.65 g, 1 piece                               |
| 8547408  | Tungsten carbide ball, Ø about 7 mm, weight 2.6 g, 1 piece                                |
| 8547300  | Tungsten carbide ball, Ø about 9 mm, weight 5.63 g, 1 piece                               |

## 5 Applications

### 5.1 Application Information Overview

This section gives an overview of the application and practical experiences known for the precursor of the Mikro-Dismembrator S, the Mikro-Dismembrator II. The operation principles and results can be converted to the Mikro-Dismembrator S. The examples shall only provide an overview on the possible universal use of this device and, however, do not claim completeness. More detailed information, other examples, application reports or literature information are available on request.

Sartorius Stedim Biotech GmbH intends to continuously add more application information with time. For this we would like to ask you for your kind support and provide us with additional information on your applications and results, if possible. Please contact our responsible product manager for this (see address given in the Introduction to this manual, see page 2 "Guide through this manual").

| Field of application                           | Analysis   | Examples   |
|--|--|--|
| Medicine                                       | assay of inorganic material from animals / human beings                      | grinding of teeth, bones, hairs, etc.  |
|  | tissue analysis  | grinding of soft tissue (kidney, liver, collagenous organs);<br>after deep-freezing or freeze-drying   |
|  | allergenic research  | grinding of wool, feathers, chemical fibres of clothing, etc.  |
|  | trace element analysis e   | sample preparation according to Lyengar  |
| Biology, biotechnology, environmental research | analysis for general applications  | grinding of solid, tough or soft materials;<br>after deep-freezing or freeze-drying                    |
|  | botanic assays, analysis for environm. protection                            | grinding of parts of plants (leaves, needle-leaves, small fruits, seeds, grain,...)                    |
|  | analysis for biotechn. research  | disruption / disintegration of microorganisms, yeast cells or tissue cells                             |
| Chemistry                                      | chemical analysis, in general  | preparation of samples for infrared-spectroscopy, radiographic assays, etc. (RF, Debye-Scherrer, etc.) |
|  | assays in dye industrie  | grinding and homogenizing of pigments for preparation of test dyes or dye compositions                 |
|  | assays for material science in mining or nuclear research and industry, etc. | grinding of oxides;<br>finest grinding of ore, coal, etc.  |

## 5.2 Use of the Mikro-Dismembrator S in Trace Analysis and Research of Contaminations

### 5.2.1 Trace Element Analysis via Homogenization by "Brittle Fracture Technique"

The use of the Mikro-Dismembrator II <sup>5-2)</sup> in trace element analysis and in other applications has proven to be an optimum method to prepare samples for analysis. Applications include grinding and homogenization of soft tissues (after deep-freeze or freeze-drying), medium solid and solid organic and inorganic material from biological source (teeth, bones, hairs, nails, etc.). With the brittle fracture technique (BFT) any kind of biological material, even if temperature-sensitive, can be homogenized without any problems. The BFT - technique uses the brittleness of the material under the temperature of liquid nitrogen (-196 deg C).

1. Fill test sample together with PTFE grinding balls into a PTFE shaking flask. Immerse shaking flask into liquid nitrogen for a few seconds to minutes.
2. Immediately afterwards mount shaking flask into the holder of the homogenizer and start the device. A running time of about 1-2 minutes may suffice for optimum grinding (required time will be according to the necessary fineness for the assay).

The BFT technique allows for fast homogenizing of representative samples of bovine liver, for instance. Assay methods without preceding homogenization, often show large deviations of single analysis, resulting in broad distribution (straggling) of the results. Preceding homogenization provides highly reproducible results with low straggling. See fig. 8 for details.


Fig. 8: Distribution profile of the some element in trace element analysis of bovine liver (6 assays per substance)

part A: distribution profile without homogenizing

part B: distribution profile after homogenizing

<sup>5-2)</sup> Reported for the precursor of the Mikro-Dismembrators S, the Mikro-Dismembrator II.

### 5.2.2 Trace - Analysis of Contaminations by Toxic Substances using the Brittle Fracture Technique

The test sample is covered in the closed PTFE flask. It will not come into contact with the liquid nitrogen. The high-purity PTFE of the shaking flask and the grinding balls ensure a very low contamination of the sample with the fluor ingredients of the PTFE. The table in section 5.2.3 below gives an overview of trace elements of human tissue which can be detected after homogenization of the sample with a Mikro-Dismembrator.

### 5.2.3 Classification of the Elements, which could be Detected in Human Tissue after Homogenization with a Mikro-Dismembrator S

|                |  |
|----------------|--|
| Basic elements | C, Ca, H, N, O, P, S   |
| Elektrolytes   | Ca, Cl, K, Mg, Na (as well as $\text{HCO}_3^-$ , $\text{SO}_4^{2-}$ , $\text{HPO}_4^{2-}$ )  |
| Trace elements | <ul style="list-style-type: none"><li>1 Main ingredients<ul style="list-style-type: none"><li>1.1 Biological important<ul style="list-style-type: none"><li>1.1.1 Group 1 Co, Cr, Cu, Fe, I, Mn, Mo, Se, Zn</li><li>1.1.2 Group 2 As, Ni</li></ul></li><li>1.2 Clinical important (Cr), Cu, Fe, I, (Se), Zu</li></ul></li><li>2 Toxic<ul style="list-style-type: none"><li>2.1 Potentially toxic elements (As, Be, Cd, Hg, Pb, Se, Tl)</li><li>2.2 Most important environmental contaminants (Cd, Hg, Pb)</li><li>2.3 Industrial contaminants like Be, Cr, Mn, Ni, Sb, Si</li></ul></li><li>3 Others Ag, (Al, Av, B, Ba, Br, Ce, Cs, Ga, Ge, Li, Nb, Pt)<br/>rare-earths like Rb, Sc, Sr, Te, Ti, Zr</li></ul> |

## 6 Supplement

### 6.1 Recommended Combinations of Shaking Flasks and Grinding Balls

| Shaking flasks       | Stainless steel | PTFE |
|----------------------|-----------------|------|
| Grinding balls       |                 |      |
| PTFE with steel core | ■               | □    |
| quartz agate         | --              | ■    |
| Chromium steel       | ■               | □    |
| Tungsten carbide     | ■               | □    |

■ = optimum combination    □ = possible combination

### 6.2 Basic Analysis of Materials for Shaking Flasks and Grinding Balls

#### 6.2.1 Materials and Characteristics

| Material         | Chemical composition, special notes   |
|------------------|---|
| PTFE             | high purity PTFE, very low resistant to abrasion; suitable for trace element analysis; not recommended for abrasive samples   |
| stainless steel  | steel classif. 1.4301; corrosion resistant, X5CrNi 1810 (acc. to DIN 17440), medium resistant to abrasion; composition: Fe = 68%, C = 0.07%, Si = 1.0%, Mn = 2.0%, Cr = 17-20%, Ni = 8.5-10%  |
| Chromium steel   | specific weight 7.85 g/cm <sup>3</sup> , medium resistant to abrasion; composition: Fe = 86%, C = 1.5%, Cr = 12%, trace ingredients of W, Mo or V   |
| Tungsten carbide | Specific weight 14.75 g/cm <sup>3</sup> , well suited for abrasive samples; high resistant to abrasion; composition: C = 5.7 - 5.8%, Co = 5.6 - 6%, Ti+Ta = 0.7 - 1%, Fe < 0.1%   |
| quartz agate     | Brazilian agate of spec. weight: 2,6 g/cm <sup>3</sup> , well suited for IR-, UR, and X-ray fluorescent assays and trace element analysis; highly resistant to abrasion; composition: SiO <sub>2</sub> = 99.91%, Al <sub>2</sub> O <sub>3</sub> = 0.02%, Na <sub>2</sub> O = 0.02%, Fe <sub>2</sub> O <sub>3</sub> = 0.01%, K <sub>2</sub> O = 0.01%, MnO = 0.01%, CaO = 0.01%, MgO = 0.01% |

## 6.2.2 Manufacturer's Data for Grinding Balls made of Glass ("Dragonit 30")

| Features                |                               | Parameter   | Unit               |
|-------------------------|-------------------------------|-------------|--------------------|
| Chemical composition    | SiO <sub>2</sub>              | 59          | %                  |
|                         | PbO                           | 24          | %                  |
|                         | K <sub>2</sub> O              | 9           | %                  |
|                         | Na <sub>2</sub> O             | 4           | %                  |
|                         | B <sub>2</sub> O <sub>3</sub> | 2           | %                  |
|                         | + trace elements              |             |                    |
| density                 |                               | approx. 2,9 | g/cm <sup>3</sup>  |
| powder density          | dep. on size                  | 1,2-1,9     | kg/l               |
| transformation temp. Tg | 1013 Poise                    | 424         | deg C              |
| deformation point       | 107.6 Poise                   | 632         | deg C              |
| process point           | 104 Poise                     | 975         | deg C              |
| refractiv index         |                               | 1.556       | n <sub>D</sub>     |
| wide angle              |                               | 58:65       | n                  |
| hardness (Mohs)         |                               | 5           |                    |
| „ (Vickers)             | 200 P.                        | 475         | kp/mm <sup>2</sup> |
| „ (Knoop)               | 200 P.                        | 460         | kp                 |
| „ (Rockwell)            | acc. to DIN 50 103            | 47          |                    |
| chem. consistency       | acc. to                       |             |                    |
| acid class              | 12 116                        | 3           |                    |
| hydrol. class           | 12 111                        | 4           |                    |
| alcali resistance       | 52 322                        | 3           |                    |
| elasticity modulus      |                               | 6200        | kp/mm <sup>2</sup> |
| shear modulus           | 20 deg C                      | 2560        | kp/mm <sup>2</sup> |
| Poisson ratio           | 20 deg C                      | 0.210       |                    |
| thermal conductivity    | 0 deg C                       | 0.856       | W/Km               |
|                         | 100 deg C                     | 0.971       | W/Km               |
| medium specific heat    | 16-100 deg C                  | 0.667       | kJ/kg K            |

## 6.3 Bibliographical Data

R.K. Wagener, L Görlich, P.W. Jungblut: Dihydrotestosterone Receptor in Human Mammary Cancer; Acta endocr.(kbh)., Suppl. 173 (1973) 65.

E. Wagner, 5. Frosch, O. Kempf: Endogenous Rhythmicity and Energy Transduction; Plant Science Letters, 3 (1974) 43-48

G.F.Deitzer, O. Kempf, S. Fischer, E. Wagner: Endogenous Rhythmicity and Energy Transduction; Planta (Berl.) 117 (29-41), 1974

E.O.R.T.C Breast Cancer Cooperative Group: Standards of the Assessment of Estrogen Receptors in Human Breast Cancer; Europ. J. Cancer Vol. pp. 379-381 (1973)

S. Berndt, U. Schwabe: Effect of Psychotropic Drugs on Phosphodiesterase and Cyclic AMP Level in Rat Brain; Brain Research 63 (1973) 303-312

E.O.R.T.C, M. Schmitt, K Mengele, E. Schueren et al: European Organisation for Research and Treatment of Cancer; European Journal of Cancer 43 (2007), 835-844

**6.4 Safety Instructions Homogenizers**

(Safety Instructions see appendix)

**6.5 EG Statement of Conformity**

(EG Statement of Conformity see appendix)

**6.6 Declaration of Decontamination**

(Declaration of Decontamination, see Appendix)

**6.7 Information and Instructions on Disposal and Repairs**

(Information and Instructions on Disposal and Repairs see appendix)

**6.8 Dimension sheet Mikro-Dismembrator S**

(Dimension sheet see appendix)

**6.9 Prospectus Mikro-Dismembrator**

(See attached information about the accessories)

Sartorius Stedim Biotech GmbH  
August-Spindler-Straße 11  
37079 Göttingen, Germany

Telephone +49.551.308.0  
Telefax +49.551.308.3289  
www.sartorius-stedim.com

Copyright by  
Sartorius Stedim Biotech GmbH,  
Göttingen, Germany.  
All rights reserved. No part  
of this publication may be  
reprinted or translated in any  
form or by any means without  
the prior written permission  
of Sartorius Stedim Biotech GmbH.

The status of the information,  
specifications and illustrations in  
this manual is indicated by the  
date given below.  
Sartorius Stedim Biotech GmbH  
reserves the right to make  
changes to the technology,  
features, specifications, and  
design of the equipment  
without notice.

Status:  
February 2008,  
Sartorius Stedim Biotech GmbH,  
Göttingen, Germany

Printed in Germany on paper that has  
been bleached without any use of chlorine  
W4A000 · KT  
Publication No.: SB6046-t08023  
Order No.: 85030-522-30

## **Safety Instructions**

**Homogenizers**

## **Content**

|          |   |          |
|----------|---|----------|
| <b>1</b> | <b>Introduction</b>                       | <b>3</b> |
| <b>2</b> | <b>Transport and Installation</b>         | <b>3</b> |
| 2.1      | Transporting Equipment                    | 3        |
| 2.2      | Installation Location Requirements        | 3        |
| <b>3</b> | <b>Operating Instructions</b>             | <b>4</b> |
| 3.1      | Working with Media                        | 4        |
| 3.2      | Initial Startup and Normal Operation      | 4        |
| 3.2.1    | Device-specific Information               | 4        |
| <b>4</b> | <b>Cleaning, Maintenance, and Service</b> | <b>5</b> |
| 4.1      | Cleaning                                  | 5        |
| 4.2      | Maintenance and Service                   | 5        |
| 4.3      | Returning Defective Units                 | 5        |
| 4.4      | Equipment Relocation Measures             | 5        |

# 1 Introduction

These Safety Instructions inform you about precautions and handling procedures for the mechanical homogenizers

- Potter S
- MSK (cell homogenizer)
- Mikro-Dismembrator U and S (laboratory scale grinding mills)

belonging to the product program of Sartorius Stedim Biotech GmbH, especially for

- transport and installation
- equipping and initial startup
- operation
- cleaning, maintenance, and service

These safety instructions are nontransferable. They are valid exclusively for the delivered equipment. Additional, specific safety regulations (i.e., legally or otherwise mandatory rules) are not dealt with by these instructions. Please consult your local authorities for further information.

Possible risks and dangers are marked with this symbol and highlighted like this paragraph. Ignoring these warnings may cause damage to the instrument or lead to other material or personal damage.

Steps to be carried out with special care or special aspects or issues to be considered are marked with safety instructions formatted as this paragraph.

[→ ..] Identifies references to contents of these instructions or other documents. Labels of illustrations, sections, or documents are in parentheses.

# 2 Transport and Installation

## 2.1 Transporting Equipment

Some of the homogenizers are quite heavy. Use only suitable transportation tools when moving homogenizers or setting up units.

Always verify that the clearance is sufficient along the transport route.

Make sure third parties are protected from injury during the transport.

## 2.2 Requirements at the site of installation

1. The units should be placed on horizontal and non-slip laboratory benches or tables. The benches or tables must be stable enough to carry the full weight of the completely equipped units.

The units may transmit vibrations to the installation surface. For this reason, install the units so that surrounding equipment is not affected.

2. The power supply must match the rating indicated on the type plate and must be grounded correctly.

### 3 Operating Instructions

1. Make sure, that only authorized people have access to the working area.
2. Personnel has to be instructed about the applicable safety instructions and other legally or otherwise compulsory regulations.

#### 3.1 Working with Media

 The utilized media and substances can be hazardous and may pose dangers that are specific to the substances and the processes and thus cannot be described here in detail. Your company should issue corresponding safety regulations and carefully instruct the operating personnel.

 Recommended basic protection equipment includes such items as suitable working clothes, gloves, goggles, and respiratory equipment, if necessary.

#### 3.2 Initial Startup and Normal Operation

1. Use only units, equipment, accessories, and spare parts released by the company Sartorius Stedim Biotech GmbH, for use with the homogenizer.
2. Check all components for damage, especially glass components such as jars. Any damaged components must not be used.
3. Switch off the unit before installing the holders and sample tubes. Make sure all parts are attached firmly.
4. Before starting and during operation, make sure the sample tubes are positioned firmly in their holders. Switch off the unit immediately if the sample container becomes loose during operation.
5. After stopping the unit, wait until the sample container or the motor has come to a complete stop before removing the tube.

##### 3.2.1 Device-specific Information

###### 1. Mikro-Dismembrator U and S

- When closing the cover, a safety switch releases the unit for operation. Never bypass this switch to start the unit while the cover is open.

###### 2. Potter S

- Reduce the plunger pressure or the speed if the overload indicator is illuminated for an extended period.
- Do not touch the plunger rod during operation. Wear adequate working clothes and head protection. Prevent hair or clothing from being caught by the rotating shaft. If necessary, switch off the unit immediately.

### 3. Cell Homogenizer MSK

 Before starting, make sure the protective cover for the pulleys at the rear of the unit is closed and secured with the wing nut. Close the carrier tube carefully. Close the protective cover at the front of the unit and fasten it with the clamping lever.

 Do not touch the belt, pulley, or the carrier tube when the unit is active. Wear suitable working clothes and head protection if necessary. Prevent hair or clothing from being caught by the belt or the pulley. If necessary, switch off the unit immediately.

 Disconnect the unit from the mains supply before changing the belt position (e.g., to change the speed).

#### When operating with CO<sub>2</sub> coolant:

 Avoid unprotected contact with cold surfaces or liquids. Risk of injury!

 Provide sufficient ventilation to prevent an accumulation of CO<sub>2</sub>. Asphyxiation hazard!

## 4 Cleaning, Maintenance, and Service

Regular cleaning and maintenance ensures the equipment's functionality and safety.

Comply with the legal safety regulations applicable to the field of application concerning the handling of microbial material and contaminated equipment.

Prior to any cleaning or maintenance work, disinfect or sterilize any biologically contaminated components according to the safety regulation.

Prior to any cleaning or maintenance work, switch unit off and pull the power cable to disconnect the unit from the mains power supply. Make sure that the unit cannot be restarted or activated accidentally.

Servicing and repairing drive systems, replacing drive belts, or working on electrical/electronic equipment should only be performed by qualified service personnel.

### 4.1 Cleaning

- Clean the surfaces of the units with standard laboratory cleaners. Do not use corrosive agents (e.g., chlorous agents).

### 4.2 Maintenance and Service

Maintenance and troubleshooting by the user should be restricted to the following tasks:

- Checking the power supply in the lab and the equipment connections.
- Checking and replacing defective fuses.

If problems or malfunctions cannot be solved or remedied, please contact your Sartorius AG representative or contact us directly:

- Sartorius AG  
Servicezentrum Nord  
Weender Landstraße 94-108  
37075 Göttingen, Germany  
Telephone +49.551.308.3729 / 3740 / 3741  
Fax +49.551.308.3730

A written cleaning and decontamination certificate has to be provided when requesting service technicians of the company Sartorius AG to perform maintenance tasks.

### 4.3 Returning Defective Units

Defective units or components can be returned to the responsible Sartorius AG representative or directly to Sartorius AG.

1. The equipment has to be clean, hygienically flawless, and carefully packaged.
2. Contaminated parts have to be cleaned, disinfected, or sterilized according to the respective and applicable safety guidelines. Compliance with the guidelines has to be documented with a cleaning and decontamination certificate to be enclosed with the shipment.
3. Negotiate the return shipment with the Service Dept. before shipping any equipment

### 4.4 Equipment Relocation Measures

1. Comply with the described cleaning and maintenance measures described to avoid biological or chemical hazards.
2. Disconnect the equipment from the mains power supply. Remove any connected accessories and transport these separately.
3. Read and comply with the additional notes under "Transporting Equipment" (2.1).

Sartorius Stedim Biotech GmbH  
August-Spindler-Straße 11  
37079 Göttingen, Germany

Telephone +49.551.308.0  
Fax +49.551.308.3289  
www.sartorius-stedim.com

Copyright by  
Sartorius Stedim Biotech GmbH,  
Göttingen, Germany.  
All rights reserved. No part  
of this publication may be  
reprinted or translated in any  
form or by any means without  
the prior written permission  
of Sartorius Stedim Biotech GmbH.

The status of the information,  
specifications and illustrations in  
this manual is indicated by the  
date given below. Sartorius  
Stedim Biotech GmbH  
reserves the right to make  
changes to the technology,  
features, specifications, and  
design of the equipment  
without notice.


Status:  
February 2008,  
Sartorius Stedim Biotech GmbH,  
Göttingen, Germany

Printed in Germany.  
Printed on paper bleached without chlorine  
W4A000 · KT  
Publication No.: SB6042az08022  
Order No.: 85030-520-66

**Declaration of Conformity**  
according to EU Machinery Directive 89/392/EU,  
Appendix II A

|   |   |   |
|---|---|---|
| Company   | <b>Sartorius Stedim Biotech GmbH</b>  |   |
| Address   | August-Spindler-Strasse 11, 37079 Goettingen; Germany<br>Phone +49.551.308.0, Fax +49.551.308.32 89<br>www.sartorius-stedim.com   |   |
|   | <p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized.</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p> |   |
| Designation of the device                               | Mikro-Dismembrator S  |   |
| Model, version  | Laboratory Grinding Mill  |   |
| Cat.-No.  | 8531609   |   |
| Relevant directives of the EC                           | EU Machinery Directive (89/392/EU) in version 93/44/EU<br>EU Directive on Electromagnetic Compatibility (89/336/EU) in version 92/31/EU<br>EU Low-Voltage Equipment Directive (73/23/EU) revised by Directive 93/68/EEC   |   |
| Applied harmonized standards                            | EN 292-1, EN 292-2<br>EN 50081-1, EN 50082-1, EN 55011<br>EN 61010-1  |   |
| Applied national standards and technical specifications |   |   |
| Date and Signature                                      |   |   |
| <u>11.01.08</u>   |    |   |
| Function of the Signatory                               | V. Niebel<br>Managing Director  | Dr. Susanne Gerighausen<br>Quality Management<br>Representative |

## Declaration about decontamination and cleaning of equipment and components

When returning equipment or components, please describe on page 2 of this form the problem(s) or fault(s) you have found. Please also indicate the remedial actions you require.

To protect our personnel, we require all equipment or components be free of biological, chemical, or radioisotopic contaminants. We will only accept such equipment or components when:

- the equipment or components have been adequately cleaned and decontaminated.
- this declaring document has been completed, signed and returned by an authorized person.

Please help us in assuring a safe, hazard-free work environment.

### A. Description of the Equipment or Component(s)

|                              |  |
|------------------------------|--|
| Description / Cat. No.       |  |
| Serial no.                   |  |
| No. of invoice/delivery note |  |
| Date of delivery             |  |

### B. Contamination / Cleaning

|  |  |
|--|--|
| Attention: Please specify exactly the biological, chemical, or radioisotopic contaminant | Attention: Please describe the cleaning and decontamination procedure/method |
| The equipment was contaminated with  | and it has been cleaned and decontaminated by                                |
|  |  |

|                              |   |                             |  |
|------------------------------|---|-----------------------------|--|
| F-13-01<br><br>Version<br>01 | Page<br>2 / 4<br><br>Effective Date<br>01.01.2008 | <h2>Return of material</h2> |  |
|------------------------------|---|-----------------------------|--|

### C. Legally binding declaration

I (we) certify that all information given in this form is correct and complete.  
 The equipment and components have been adequately decontaminated and cleaned according to the legal requirements.  
 No chemical or biological or radioisotopic risks remain that can endanger exposed persons' safety or health.

|                               |  |
|-------------------------------|--|
| Company / Institute           |  |
| Address / Country             |  |
| Tel. / Fax (with area code)   |  |
| Name of the authorized person |  |
| Position                      |  |
| Signature / Date              |  |

### D. Reason for return

wrong delivery   
  exchange   
  repair   
  modification   
  disposal   
  other

### E. Please describe the problem(s) or fault(s) you have found (for repair) and/or indicate the remedial actions you require

|               |                              |                           |  |
|---------------|------------------------------|---------------------------|--|
| F-13-01       | Page<br>3 / 4                | <b>Return of material</b> |  |
| Version<br>01 | Effective Date<br>01.01.2008 |                           |  |

## F. Reserved for Sartorius-Service Center

Notes

Please pack the unit properly and send it freight paid to your local service supplier or directly to Sartorius AG Servicezentrum Nord, Germany.

Sartorius AG  
Servicezentrum Nord  
Weender Landstraße 94-108  
37075 Göttingen, Germany  
Telephone +49.551.308.3729 / 3740 / 3741  
Fax +49.551.308.3730

|               |                              |                           |  |
|---------------|------------------------------|---------------------------|--|
| F-13-01       | Page<br>4 / 4                | <b>Return of material</b> |  |
| Version<br>01 | Effective Date<br>01.01.2008 |                           |  |

Sartorius Stedim Biotech GmbH  
August-Spindler-Straße 11  
37079 Göttingen, Germany

Telephone +49.551.308.0  
Telefax +49.551.308.3289  
www.sartorius-stedim.com

Copyright by  
Sartorius Stedim Biotech GmbH,  
Göttingen, Germany.  
All rights reserved. No part  
of this publication may be  
reprinted or translated in any  
form or by any means without  
the prior written permission  
of Sartorius Stedim Biotech  
GmbH.

The status of the information,  
specifications and illustrations in  
this manual is indicated by the  
date given below. Sartorius Stedim  
Biotech GmbH reserves the right  
to make changes to the technology,  
features, specifications, and  
design of the equipment  
without notice.


Status:  
February 2008,  
Sartorius Stedim Biotech GmbH,  
Göttingen, Germany

## Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:

Sartorius AG  
Service Center  
Weender Landstrasse 94-108  
37075 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.


Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet/manual or visit our Internet website ([www.sartorius.com](http://www.sartorius.com)) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

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**Mikro-Dismembrators**  
Laboratory ball mills  
for fine grinding

## Mikro-Dismembrators Laboratory ball mills for fine grinding

The Mikro-Dismembrator U and Mikro-Dismembrator S from Sartorius have proven their worth in many laboratories world-wide. Even the most difficult materials can be reliably disintegrated. Whereas the Mikro-Dismembrator U with its shaking frequency of  $2.000 \text{ min}^{-1}$  is suitable for most applications, the Mikro-Dismembrator S with a maximum frequency of  $3.000 \text{ min}^{-1}$  is especially well suited for fast homogenization of samples frozen in liquid nitrogen, for instance of biopsych materials in cancer diagnosis.

Thanks to the high efficiency of shaking with grinding balls, samples are rapidly disintegrated – resulting in a much faster isolation of DNA, RNA and proteins than with protocols involving manual grinding processes, thereby reducing the risk of sample degradation.

In addition, the use of enzymes, which are often employed for cell lysis, can be replaced. Thus the introduction of foreign proteins, which may interfere with subsequent analysis, is avoided. This is particularly advantageous in proteome research.

Reproducibility of the disintegration process is guaranteed by the digitally controlled shaking frequency and shaking time. This ensures comparable results even within large series of experiments – another big advantage over any manual processes.

In botany or environmental analysis, the Mikro-Dismembrators serve to grind needles, seeds and other hard plant materials. In food analysis, hair samples are pulverized to check whether animals have been treated with illegal feed additives. In medical research, pathology or toxicology, the units are used to disrupt frozen tissue samples, hair, bone or cartilage samples or teeth. Microbiologists can disintegrate frozen cell sediments. In chemistry, pigments are ground to a uniform particle size.

In combination with a large range of various shaking flasks, grinding balls and glass beads, the Mikro-Dismembrator units can be adapted to many different and individual applications.

If you have other homogenization tasks in your laboratory, Potter S homogenizer, LABSONIC® M and P ultrasonic disintegrators and Dounce hand homogenizers may be an option. Please refer to our Laboratory Product Catalogue or to our website, [www.sartorius.com](http://www.sartorius.com)

### Benefits

Rapid and highly efficient homogenization of brittle, frozen or hard materials

Fine grinding of pigments

Reproducible action due to digitally controlled shaking frequency


Homogenization time can be preset

Safety lid protects the user

Shaking flasks made from stainless steel or PTFE



Large selection of grinding balls

Pig teeth, approx. 2fold

Pig teeth, homogenized, approx. 100fold

Geranium spec, approx. 2fold

Geranium spec, homogenized, approx. 100fold

## Accessories

| Cat.-No. | Description  |
|----------|--|
| 853 1803 | <b>Shaking flasks made of stainless steel 1.4301</b><br>Shaking flask, stainless steel 1.4301, volume approx. 3 ml, with PTFE gasket and cap |
| 8531811  | Shaking flask, stainless steel 1.4301, volume approx. 5 ml, with PTFE gasket and cap   |
| 853 1820 | Shaking flask, stainless steel 1.4301, volume approx. 7 ml, with PTFE gasket and cap   |

|          |   |
|----------|---|
| 853 1838 | <b>Shaking flasks made of PTFE</b><br>Shaking flask, PTFE, volume approx. 3 ml, with cap                        |
| 853 1846 | Shaking flask, PTFE, volume approx. 5 ml, with cap  |
| 853 1854 | Shaking flask, PTFE, volume approx. 7 ml, with cap  |
| 853 1862 | Shaking flask, PTFE, volume approx. 20 ml, with cap<br>Holder 8531897 is required for using this shaking flask! |

|          |   |
|----------|---|
| 853 1943 | <b>Shaking flasks made of PTFE</b><br>Shaking flask, PTFE, volume approx. 3 ml, with screw cap                        |
| 853 1935 | Shaking flask, PTFE, volume approx. 5 ml, with screw cap  |
| 853 1927 | Shaking flask, PTFE, volume approx. 7 ml, with screw cap  |
| 853 1951 | Shaking flask, PTFE, volume approx. 20 ml, with screw cap<br>Holder 8531897 is required for using this shaking flask! |

|          |   |
|----------|---|
| 853 1889 | <b>Containers for disposable tubes, holder</b><br>Container for 3 disposable test tubes 2.2 ml $\varnothing$ 10.8×37 mm, for instance Sarstedt no. 72.608<br>Holder 8531897 is required for using this shaking flask! |
| 853 1960 | Container for 4 cryotubes<br>Holder 8531897 is required for using this shaking flask!   |
| 853 1897 | Holder for shaking flask 20 ml (8531951) and for the containers for disposable tubes (8531889 and 8531960)  |

|          |   |
|----------|---|
| 853 1900 | <b>Special accessories</b><br>Tray for microwell plates, capacity of 2 microwell plates with 96 borings, for using the Mikro-Dismembrator as a small shaker |
|----------|---|

| Cat.-No.   | Description  |
|--|--|
| <b>Grinding balls</b>                                      |  |
| 854 7505   | Made of Brazilian agate<br>Grinding ball made of Brazilian agate, Ø 10 mm, weight 1.4 g, package with 10 pieces          |
| 854 7602   | Made of PTFE, with steel core<br>Grinding ball made of PTFE with steel core, Ø 10 mm, weight 2 g, package of 2 pieces    |
| 854 6606   | Made of chromium steel<br>Grinding ball made of chromium steel, specific weight 7.85 g/ml, Ø 3 mm, package of 100 pieces |
| 854 6703   | Grinding ball made of chromium steel, specific weight 7.85 g/ml, Ø 5 mm, package of 100 pieces                           |
| 854 6800   | Grinding ball made of chromium steel, specific weight 7.85 g/ml, Ø 10 mm, package of 10 pieces                           |
| <b>Made of Tungsten carbide</b>                            |  |
| 854 7009   | Grinding ball made of Tungsten carbide, specific weight 14.5 g/ml, Ø 1mm, 1 piece  |
| 854 7106   | Grinding ball made of Tungsten carbide, specific weight 14.5 g/ml, Ø 3 mm, 1 piece                                       |
| 854 7203   | Grinding ball made of Tungsten carbide, specific weight 14.5 g/ml, Ø 5 mm, 1 piece                                       |
| 854 7408   | Grinding ball made of Tungsten carbide, specific weight 14.5 g/ml, Ø 7 mm, 1 piece                                       |
| 854 7300   | Grinding ball made of Tungsten carbide, specific weight 14.5 g/ml, Ø 10 mm, 1 piece                                      |
| <b>Glass beads</b>   |  |
| 854 1400   | Glass beads Ø 0.10–0.11 mm, bottle, approx. 570 ml   |
| 854 1507   | Glass beads Ø 0.17–0.18 mm, bottle, approx. 570 ml   |
| 854 1604   | Glass beads Ø 0.25–0.30 mm, bottle, approx. 570 ml   |
| 854 1701   | Glass beads Ø 0.40–0.60 mm, bottle, approx. 570 ml   |
| 854 1809   | Glass beads Ø ca. 1 mm, bottle, approx. 570 ml   |
| <b>Adapter set for the old model Mikro-Dismembrator II</b> |  |
| 853 1986   | Adapter set for using Mikro-Dismembrator U/S flasks with the Mikro-Dismembrator II, for shaking flasks 3 ml, 5 ml, 7 ml  |

## Technical Data

### Ordering information

|         |                      |                                  |
|---------|----------------------|----------------------------------|
| 8531722 | Mikro-Dismembrator U | 230 V/50–60 Hz                   |
| 8531730 | Mikro-Dismembrator U | 115 V/50–60 Hz                   |
| 8531609 | Mikro-Dismembrator S | 115 V/230 V/50–60 Hz convertible |

### Accessories

Shaking flasks and grinding balls or glass beads are required for the grinding process. We offer flasks with sizes of 3 ml, 5 ml, 7 ml and 20 ml. Choose from grinding balls made of agate, steel or Tungsten carbide. Glass beads can only be used for wet or liquid samples. Grinding balls have to be used for frozen or brittle material. The choice of suitable flasks and grinding balls depends on the intended application but the following table provides a general guide.

| Shaking flasks        | Stainless steel | PTFE |
|-----------------------|-----------------|------|
| <b>Grinding balls</b> |                 |      |
| PTFE   steel core     | ■               | □    |
| Agate                 |                 | ■    |
| Stainless steel       | ■               | □    |
| Tungsten carbide      | ■               | □    |

■ = very suitable, □ = suitable

### Technical data

|  |                                  |
|--|----------------------------------|
| Shaking amplitude:                         | approx. 16 mm                    |
| Shaking frequency<br>Mikro-Dismembrator U: | 50 to 2.000 min <sup>-1</sup>    |
| Shaking frequency<br>Mikro-Dismembrator S: | 50 to 3.000 min <sup>-1</sup>    |
| Accuracy of frequency:                     | +10%, -2%                        |
| Drive:                                     | brushless drive motor            |
| Power consumption:                         | 120 W                            |
| Setting the timer:                         | membrane key pad                 |
| Setting the frequency:                     | membrane key pad                 |
| Programmable timer:                        | To 98.59 min, continuous         |
| Display:                                   | LCD, alphanumeric                |
| Admissible temperature range:              | +10°C to +60°C                   |
| Admissible humidity:                       | 10 to 90% non condensing         |
| Power supply:                              | 230 V/50–60 Hz<br>115 V/50–60 Hz |
| Dimensions (W×H×D):                        | 297×259×205 mm                   |
| Weight (net):                              | 19 kg                            |

All units are delivered without shaking flasks and other accessories.

Sartorius AG  
Weender Landstrasse 94-108  
37075 Goettingen, Germany

Phone +49.551.308.0  
Fax +49.551.308.3289

[www.sartorius.com](http://www.sartorius.com)

Sartorius BBI Systems GmbH  
Schwarzenberger Weg 73-79  
34212 Melsungen, Germany

Phone +49.5661.71.3400  
Fax +49.5661.71.3702

[www.sartorius-bbi-systems.com](http://www.sartorius-bbi-systems.com)

Sartorius North America Inc.  
131 Heartland Blvd.  
Edgewood, New York 11717, USA

Phone +1.631.254.4249  
Toll-free +1.800.3687178  
Fax +1.631.254.4253

Sartorius BBI Systems, Inc.  
2800 Baglyos Circle  
Bethlehem, PA 18020, USA

Phone +1.610.866.4800  
Fax +1.610.866.4890

Sartorius Ltd.  
Longmead Business Park  
Blenheim Road, Epsom  
Surrey, KT19 9 QQ, U.K.

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Phone +81.3.3740.5407  
Fax +81.3.3740.5406

Specifications subject to change  
without notice.

Printed in Germany on paper that has been  
bleached without any use of chlorine.

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The reproduction, distribution and utilization of this document as well as the communication of its contents to others without explicit authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.
Copyright by Sartorius Stedim
Sartorius Stedim Singapore Pte. Ltd.
Sartorius Stedim Australia Pty. Ltd.