TOPLOADER BALANCES











# The World's Most User-Friendly Balances...

Our revolutionary Super Hybrid Sensor (SHS) technology improves response speed, accuracy and minimizes maintenance costs. SHS combines the best of magnetic force restoration and single point parallelogram load sensing technologies. Its unique design provides the fastest response speed in its class-just 1 second! And if the time ever comes when the SHS needs repair, you just replace one of three parts-there is no need to replace the whole sensor! This translates to Minimum Maintenance Costs.

(An additional benefit of our new SHS is that the GX automatically performs self calibration using the internal mass when the balance detects ambient temperature changes.)

# Super Hybrid Sensor





# **Motor-driven Internal Calibration**

# **Internal Calibration**





The GX internal calibration with an internal mass saves time and money. External calibration requires a great deal of skill, time and an expensive external mass. With our internal calibration with internal mass, there is no downtime while the technician calibrates the halance

- Automatic Self Calibration Self-calibrates automatically when the balance detects ambient temperature changes.
- One Touch Automatic Calibration Lets you calibrate the balance on demand with just the push of one button.
- Calibration is essential for accuracy in weighing. For example, in high resolution balances, changes in sensitivity drift can be caused by changes in the ambient temperature. Take the GX 600 for example, the weight might change by ±0.010g with a change of 10°C. (Sample's weight: 500g)
- The GX's internal calibration mass is adjusted to an "OIML F1 class level." by the standard balance calibrated with an "OIML E2 class level mass". (The GX's internal calibration mass is not "the standard mass for legal metrology" nor "OIML's level of calibration mass." Instead it is the mass for adjustment.)

# Just 1 Second — Ultra Fast **Response Speed**

- Ultra Fast Response Speed

   Just 1 Second.\* Ideal for installation on filling machines.
- Ultra High Performance with minimum maintenance costs.

#### Standard RS-232C



Incorporates bi-directional communication with a PC, printer or other peripheral device.

# Standard Windows Communication WinCT **Tools Software WinCT**



A free copy of WinCT-A&D's powerful and flexible data collection software tool for virtually instant connection to a PC and a network of other devices.





# **TOPLOADABLE BALANCES**

### **Splash Proof Keyboard and Display**


Full Digital Calibration

Balances protected from dust and spills up to IP54 level

# **GLP/GMP/GCP/ISO Compliance**

Full Digital Calibration

Allows GLP or LIMS balance management by outputting the Balance ID number and data used to calibrate the balance. The data can be output to A&D's AD–8121 printer or a computer, indicating date, time, Balance ID number, serial number and calibration data.



#### Standard RS-232C

FDCW

Bi-directional communication with a PC or connects directly to a printer

#### **Commands**

**EDC**II

A PC can control the balance by sending commands to the balance.

- Commands to query weighing data
- Commands to control the balance
- Commands to control the memory function
- Commands to control the comparator function

#### **ID Number**

FDCU

The balance ID number can be set. It is used to identify the balance when GLP is used. The ID number is memorized and maintained once it is fixed, even if the AC adapter is removed.

#### **Display Auto Power Off Function**

**EDCIII** 

Display turns off after 10 minutes of inactivity (can be disabled when undesired)

#### **Data Memory Function** (Patent Pending)



Format the Data Memory Function for:

- **GX** 200 sets of weighing data
  - 100 sets of weighing data with Time & Date
  - 50 sets of GLP Data-Time & Date, Calibration Data, ID number and Serial number
  - 20 sets of Unit Weight for Counting Mode
- GF 40 sets of weighing data
  - 20 sets of Unit Weight for Counting Mode

#### Time & Date (GX only)

FDCW

Standard Time & Date Function complies with GLP and Interval Weighing in the Data Memory Function

## **Auto Self Checking**

FDCW

Automatically checks itself when setting the Automatic Adjustable Environment

#### **Standard Underhook**

FDCW

Ideal for density determination and weighing magnetic substances.

#### **Density Measurement Function**

Picture shown with GX-13

Density Measurement Kit



# **Interval Time Setting**

FOCU

Weighing intervals of 2, 5, 10, 30 seconds and 1, 2, 5, 10 minutes in the Data Memory Function

#### **Auto Power On Function**

FOCU

Plug it in and it turns itself on without pressing ON/OFF key

#### **Auto Re-Zero Function**

FDCW

Re-Zeroes automatically after data output

#### **Large Square Weighing Pan**

FDCW

Standard large square weighing pan on all models 128 x 128mm for 0.001g resolution models 165 x 165 mm for 0.01g and 0.1g resolution models

# Multiple Weighing Units & Programmable Unit

テクこと

Uses Programmable or standard units of measure (g, ct, dwt etc.)

#### **Front Air Bubble Leveling**

アクこと

Easily check the balance's level from the front

#### **Quick Reference Card**

**50**(1)

A fast, convenient operating guide installed at the bottom of the balance





# **GX/GF Series Specifications**

|  |                          |  |  | GX-400 | GX-600  | GX-800         | GX-1000         |                               | GX-2000         |               | GX-4000 | GX-6100          | GX-6000                                  | GX-8000 |  |
|--|--------------------------|--|--|--------|---------|----------------|-----------------|-------------------------------|-----------------|---------------|---------|------------------|--|---------|--|
|  |                          | GF-200*  | GF-300*  | GF-400 | GF-600  | GF-800         | GF-1000         | GF-1200*                      | GF-2000*        | GF-3000       | GF-4000 | GF-6100          | GF-6000*                                 | GF-8000 |  |
| Weighing capacity                          |                          | 210g   | 310g   | 410g   | 610g    | 810g           | 1010g           | 1210g                         | 2100g           | 3100g         | 4100g   | 6100g            | 6100g                                    | 8100g   |  |
| Minimum weighing value                     |                          | 0.001g   |  |        |         |                | 0.01g           |                               |                 |               | 0.1g    |                  |  |         |  |
| Repeatability (Standard Deviation)         |                          | 0.001g   |  |        |         |                |                 | 0.01g                         |                 |               |         |                  | 0.1g                                     |         |  |
| Other units of measure                     |                          | Decimal ounce, Troy ounce, Pennyweight, Carat, Momme, Grain unit, Pound, Pound/Ounce |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Linearity                                  |                          | ±0.002g  |  |        | ±0.003g | 610g           | 610g            | ±0.02g ±0.03g                 |                 |               | ±0.03g  | ±0.1g            |  |         |  |
| Stabilization time (typical at FAST)       |                          | Approx. 1 second  Approx. 1.5 seconds  |  |        |         |                |                 |                               |                 |               |         | Approx. 1 second |  |         |  |
| Sensitivity drift (10°C to 30°C/50°F-80°F) |                          | ±2ppm/oC (When automatic self calibration is not used)                               |  |        |         |                |                 |                               |                 |               |         |                  | ±5ppm/°C (auto-self califbration<br>OFF) |         |  |
| Operating temperature                      |                          | 5°C to 40°C (41°F to 104°F), 85%RH or less (No condensation)                         |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Sensing method                             |                          | Super Hybrid Sensor (SHS) *SHS: Patent issued and pending                            |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Display type                               |                          | Vacuum Fluorescent Display (VFD)   |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Display refresh rate                       |                          | 5 times/second or 10 times/second  |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Counting Mini                              | imum unit mass           | 0.001g   |  |        |         |                | 0.01g           |                               |                 |               | 0.      | 1g               |  |         |  |
| mode Num                                   | nber of samples          | 10, 25, 50 or 100 pieces   |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Percent Mini                               | imum 100% reference mass | 0.100g 1.00g   |  |        |         |                |                 | 10.0g                         |                 |               |         |                  |  |         |  |
| mode Mini                                  | imum 100% display        |  | 0.01%, 0.1%, 1% (Depends on reference mass stored) |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Standard serial I/F                        |                          | RS-232C Interface with Windows Communication Tools (WinCT: included as CD-ROM)       |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |
| Weighing pan                               |                          | 128 x 128mm (5.04 x 5.04inch)  |  |        |         |                |                 | 165 x 165mm (6.50 x 6.50inch) |                 |               |         |                  |  |         |  |
| Calibration                                |                          |  |  |        |         |                | Motor-driven I  | nternal Calibr                | ation (GX only  | )             |         |                  |  |         |  |
| External dimensions                        |                          |  |  |        |         | 210(W          | /) x 317(D) x 8 | 6(H)mm (8.2                   | 7 x 12.48 x 3.3 | 9inch)        |         |                  |  |         |  |
| AC adapter                                 |                          |  |  |        | Confir  | m that the ada | apter type is c | orrect for the                | local voltage a | nd power rece | eptacle |                  |  |         |  |
| Power consumption                          |                          |  |  |        |         |                | Approx. 11V     | A (supplied by                | AC adapter)     |               |         |                  |  |         |  |
| Net weight                                 |                          | Approx. 4.6kg(GX)/3.8kg(GF) Approx. 5.1kg(GX)/4.3kg(GF)                              |  |        |         |                |                 |                               |                 |               |         |                  |  |         |  |

<sup>\*</sup>NTEP Class II versions available on request

Notes: The internal mass is equal to (approximately) 500grams. Standard breeze break included for GX-200/400/600 and GF-200/300/400/600

## **Options**

| GX-04* | Comparator Output with a Buzzer/RS-232C/ |
|--------|--|
|        | Current Loop Output                      |

GX-06\* Analog Output/Current Loop Output
GX-10 Breeze Break for GX-200/400/600 and
GF-200/300/400/600

GX-11 Glass Breeze Break for GX-2000/4000/6100/ 6000/8000 and GF1200/2000/3000/4000/ 6100/6000/8000

GX-12 Animal weighing pan for GX-400/600/2000/4000/ 6100/6000/8000 and GF-300/400/600/1200/2000/ 3000/4000/6100/6000/8000

**GX-13** Density Determination Kit for GX-200/400/600 and GF 200/300/400/600

#### Accessories

AD-8121A Dot Matrix Compact Printer

AD-8118B Universal Printer
AD-8920 Remote Display

AD-1682 Rechargeable Battery Unit

(Maximum operating time: 8hours)

**AX-K01710-200** RS-232C Cable (25P-9P) **AX-SW-128** Foot Switch

AX-073003691-S Dust Cover for GX-200/400/600 and


GF-200/300/400/600

**AX-073003692-S** Dust Cover for GX-2000/4000/6100/

6000/8000 and GF1200/2000/3000/

4000/6100/6000/8000

## **Physical Dimensions (mm/inches)**




#### Pan Size (mm/inches)

| Model  | Α        | В        |
|--|----------|----------|
| GX-200/400/600<br>GF-200/300/400/600                                 | 128/5.03 | 128/5.03 |
| GX-2000/4000/6100/6000/8000<br>GF-1200/2000/3000/4000/6100/6000/8000 | 165/6.49 | 165/6.49 |

Approx. 1 second

# Dimensions With Optional Glass Breeze Break (mm/inches)



Specifications subject to change without notice.



1555 McCandless Drive Milpitas, CA 95035 (800) 726-3364 or (408) 263-5333 Email: scales@andweighing.com www.andweighing.com



<sup>\* 04&</sup>amp; 06 cannot co-exist \* Factory installed option