

# ÄKTA™ pure

ÄKTA pure is a flexible and intuitive chromatography system (Fig 1) for fast purification of proteins, peptides, and nucleic acids from microgram to gram levels of target product. ÄKTA pure is a reliable system where hardware and UNICORN™ software are designed to work together with columns and media to meet any purification challenge. The system supports a wide range of chromatography techniques and meets the automation requirements needed to deliver the highest purity. The system is configurable and can be upgraded at any time with a wide range of options to further increase its capabilities depending on your purification needs. ÄKTA pure is the product of over fifty years expertise in protein research and three decades of experience in the development of ÄKTA purification systems.

## ÄKTA pure offers the following benefits for you:

- Modular system design with a large range of options to allow flexibility in purification of proteins and peptides
- Customizable system that is easy to upgrade as your research needs develop
- Reliable system with components and integrated features based on the proven design of ÄKTA avant
- UNICORN 6 software provides simple, intuitive, and flexible preprogrammed method templates and total system control to simplify your job
- Predefined method settings for all GE Healthcare lab-scale chromatography columns

## System overview

ÄKTA pure chromatography system is a highly versatile, modular system with a number of design features to facilitate reliable purification.

The system consists of the ÄKTA pure instrument and UNICORN 6 control software. The system is modular in design with all valves, monitors, and columns mounted on the forward facing wet side of the system, which faces front to allow easy interaction with the instrument modules (Fig 2). Additional components such as valves, monitors, and sensors



**Fig 1.** ÄKTA pure is a flexible chromatography system for the reliable purification of proteins, peptides, and nucleic acids at laboratory scale.

from the wide range of optional modules can be easily added to the available positions. Multiple rails for attachment of column holders and equipment are located at the front and on the side of the instrument. A buffer tray on the top of the instrument provides a large storage area for vessels and bottles. The instrument control panel shows the system state and allows the possibility to interact with the run (pause/continue) by the touch of a button.

The system weighs only 48 kg in basic configuration and 53 kg when fully equipped with options. The relatively low weight enables easier placement in the lab and the system dimensions allow it to fit conveniently into a standard cold cabinet for work with labile samples.



(A)

Rails for attachment of column holders and accessories (also on the left-hand side)



Verify gradients with Conductivity monitor C9n for verifying gradients

Round fraction collector F9-R for collection in various tube types



Column valve V9-Cs with built-in bypass and reverse-flow functionality


UV monitor U9-L for detection at 280 nm

Tubing organizer



Inlet valve V9-IAB allows selection between two A and two B positions in a single valve


High-performance system pumps



(B)

Sample pump S9 and Sample inlet valve V9-IS for automated sample loading of multiple samples

Flexible fraction collector F9-C for collection in deep-well plates and/or tubes




Column valve V9-C enables media screening and measures the pressure drop over the column

UV monitor U9-M for triple wavelength detection

pH valve V9-pH enables in-line pH monitoring, calibration, and storage.

Inlet automation valves V9-IA and V9-IB gives 14 inlet positions and enables condition screening. Air sensors are integrated to protect the columns

Position for additional modules



**Fig 2.** Two examples of system configurations for ÄKTA pure showing positions of modules on the front panel and flow paths for each. (A) A basic system configuration for convenient protein purification. (B) System configured for high level of automation.

Regardless of configuration, ÄKTA pure always comes with two high-performance system pumps, system pressure monitor for column protection, mixer, injection valve, and UV monitor. ÄKTA pure has a wide range of optional modules to allow a large number of possibilities. The system flow path is designed to minimize band-broadening effects, and all wetted materials used in the flow path are biocompatible and resistant to commonly used solvents. The instrument front is designed with empty module positions where optional valves and monitors can be mounted to enable a flexible configuration of the flow path. Examples of two ÄKTA pure system configurations are shown in Figure 2.

UNICORN 6 control software allows a fast and easy start to creating methods and starting runs. UNICORN 6 eliminates the need for programming skills as creation of chromatography methods is done by simple drag-and-drop operations. In addition, the software is modular allowing the addition of features such as Column Logbook and Design of Experiments (DoE) functionality for method development. Licensing options for remote access to the system and/or for creating methods or evaluating results gives even greater convenience.

ÄKTA pure system components and available options are described in the following sections in more detail.

## ÄKTA pure standard components

### System pump

The two system pumps are based on the modern technology developed for ÄKTA avant systems. The robust construction delivers reproducible flow rates at both low and high back pressures, allowing short separation times.

Each pump consists of one pair of pump heads which deliver low-pulsation flow to the mixer. The continuous and accurate flow rates generated enable reproducible isocratic or gradient elution. The system pumps provide a flow rate range of up to 25 ml/min (up to 50 ml/min when packing columns) at operating pressures of up to 20 MPa. A system pressure monitor is connected to the pumps that continuously measures system pressure and enables the flow rate to be automatically adjusted to avoid reaching any defined pressure limit.

### Mixer

The mixer ensures homogeneous buffer composition during gradient runs. The choice of mixer chamber size depends on the flow rate and buffers used. A 1.4-ml chamber size is provided as standard. Chamber sizes of 0.6 and 5 are available as options. A larger mixer volume chamber is required for higher flow rates or difficult-to-mix buffers. An inline filter is mounted on the mixer. The filter is easy to change, and the mixer is also easily changed by snapping it in or out of the mixer holder.

### Injection valve

The injection valve allows for a variety of sample application techniques using capillary loops or Superloop™. The novel valve design eliminates the need for replumbing when changing between various sample application techniques. Capillary loops can be filled manually, via a syringe, or with a sample pump; the same sample application options apply to use of Superloop. Capillary loops can also be filled using the system pump.

Sample can also be applied to the column directly using an optional sample pump or the system pump.

### UV monitoring

ÄKTA pure is equipped with either a fixed wavelength UV monitor or a variable multiwavelength UV and visible spectrum monitor.


The fixed-wavelength (280 nm) UV monitor (U9-L) incorporates LED technology, which is durable and reliable and is ready-to-use at start-up. Moreover, UV monitor U9-L does not heat the sample. The monitor is available with a 2 mm flow cell as standard and an optional 5 mm flow cell when higher sensitivity measurements are required.

To determine protein separation at different wavelengths, UV monitor U9-M is designed for multiwavelength detection in the UV and visible spectrum from 190 to 700 nm.

UV monitor U9-M allows monitoring of up to three wavelengths simultaneously (Fig 3 and 6). For optimum performance when purifying samples with different protein concentrations, there are three flow cell path lengths (0.5, 2 [default], and 10 mm) available. The optimal flow cell design, together with fiber optic technology, provides a high signal-to-noise ratio without causing any local heating of the UV flow cell. The monitor contains a high-intensity xenon lamp with a long lifetime that requires minimal start-up time. Every time the instrument is switched on, the monitor is automatically calibrated.

Monitoring with multiple wavelengths can be used to detect contaminants, specifically labeled proteins, or target molecules that do not absorb light at 280 nm. To demonstrate this, molecular weight standards were monitored at 214, 280, and 340 nm wavelengths. Detection at 214 nm reveals peptide bonds of all proteins and may be useful if the concentration and extinction coefficient at 280 nm is low for the target protein. Ferritin, a multimeric iron-storage protein, showed stronger absorbance at 340 nm than the other proteins due to the high number of ferric ions in the center of the molecule (Fig 3).

**Column:** Superdex 200 10/300 GL  
**Sample:** Molecular weight standards for gel filtration  
**Sample volume:** 100  $\mu$ l  
**Eluent:** PBS (10 mM sodium phosphate, 140 mM NaCl, 2.7 mM KCl, pH 7.4)  
**Flow rate:** 0.5 ml/min  
**System:** ÄKTA pure



**Fig 3.** Gel filtration (GF) with multiwavelength detection (214, 280, and 340 nm) of proteins using ÄKTA pure with UV monitor U9-M. The column used was Superdex™ 200 10/300 GL. The peaks observed on the chromatogram are 1) ferritin (M<sub>r</sub> 440 000), 2) aldolase (158 000), 3) conalbumin (75 000), 4) ovalbumin (44 000), 5) carbonic anhydrase (29 000), 6) ribonuclease A (13 700), and 7) aprotinin (6 500).

Both UV monitor U9-L and UV monitor U9-M can be combined with a second UV monitor U9-L to give increased application capabilities such as multistep applications or for using small and large flow-cells to simultaneously detect both low and high protein concentrations.

### Conductivity monitor

The conductivity monitor measures conductivity of buffer and samples for on-line monitoring of the true gradient. An integrated temperature sensor corrects for variations in conductivity due to the temperature. The conductivity monitor has a broad reading range and is therefore able to monitor conductivity in all different chromatographic techniques.

## ÄKTA pure optional modules for enhanced automation

### Sample application options

The optional sample pump\* (Fig 4) is designed to allow automatic sample application directly to a column or indirectly via a capillary loop or Superloop. Using the sample pump saves time by eliminating laborious sample application steps and is especially useful when handling large sample volumes. The pump consists of two pump heads and is based on the same pump principle as the system pumps. Pump purging and air removal can easily be performed automatically. The sample pump is equipped with a pressure sensor for control of the sample flow rate to protect the column while preventing pressure stops and minimizing the time for sample loading. Using the sample pump, samples can be loaded at flow rates of up to 50 ml/min.

\* To be released.

The optional sample inlet valve\* is intended to be used with the sample pump. It allows fast, automatic loading of up to 7 different samples. The integrated air sensor ensures safe and complete sample application without the need to preprogram the sample volume. The valve has seven sample inlet positions plus a dedicated buffer inlet for filling the sample pump with solution before the sample is introduced and for washing out the valve and pump between runs. During sample application, the air sensor detects when sample has been completely loaded so that the method can continue to the next step without air being introduced into the flow path or column.



**Fig 4.** ÄKTA pure sample pump.

### Buffer selection

ÄKTA pure can be equipped with two different types of inlet valves that allow selection of buffers and wash solutions. Valves with multiple inlets ensure that cleaning reagents can be permanently online, which means that columns and system can be cleaned conveniently at regular intervals.

Inlet selection valve V9-IAB comprises four inlet positions in a single valve giving a convenient solution for automation of buffers and post-run cleaning of columns and system when performing basic chromatography.

The inlet automation valves (V9-IA, V9-IB) provide up to up to 2 × 7 inlets; multiple inlets enable automatic screening of buffer and reagent conditions. Each of the inlet automation valves is equipped with an integrated air sensor, which ensures safe exclusion of air from the system. If air is detected, the system can be paused so that the air can be purged before it enters the flow path.


### Column control

A column valve can be connected to the system and used to control the flow to the column. ÄKTA pure can be equipped with one of two different column valves.

Column control valve V9-Cs allows connection of one column and has an integrated bypass function, which enables washing of the system without the need to remove the column. The column control valve also allows reverse flow through the column for fast and effective elution of strongly bound proteins, sharper bands, and a concentrated target molecule eluent.

Column selection valve V9-C also has the integrated bypass and reverse-flow functions. Connection of up to five columns for automatic column switching is possible using this valve. Connection of multiple columns minimizes manual intervention and reduces further the risk of introducing air into the column.

The column selection valve has two integrated pressure sensors; the first sensor measures pressure before the column, enabling protection of the column hardware while the second measures the pressure after the column. The pressure drop over the column ( $\Delta p$ ) is calculated by measuring the difference between the two pressure readings and can be used to protect the packed media bed (Fig 5).




**Fig 5.** For increased operational safety, the column selection valve enables continuous measurement of pre-column (Pre-CP) and post-column pressure (Post-CP) during runs. The pressure difference over the packed media bed ( $\Delta p$ ) is calculated from the two signals.

The flexibility of Column valve V9-C for connection of up to five columns was demonstrated in a column scouting study using columns for hydrophobic interaction chromatography (HIC). Five columns from HiTrap™ HIC Selection Kit were connected to ÄKTA pure and used for column scouting for optimization of purification conditions of S-aminotransaminase in clarified *E. coli* extract. UV monitor U9-M was used for detection of the protein at two wavelengths. Chromatograms of the five separate HIC runs are shown in Figure 6. Eluted fractions were analyzed using GF and SDS-PAGE (data not shown).

The  $A_{420}$  signal specifically monitors the target protein. The columns giving the sharpest and most symmetrical peaks at  $A_{420}$ , as well as the highest possible purity, were selected for subsequent optimization and scale-up experiments. HiTrap Phenyl FF (high sub) 1 ml and HiTrap Butyl FF 1 ml gave the most promising results under the conditions used, and HiTrap Phenyl FF (high sub) 1 ml was selected for further optimization in this case.

**Columns:** Five columns from HiTrap HIC Selection Kit  
**Sample:** Supernatant after precipitation with 2 M ammonium sulfate (AS) at room temperature of extract of *E. coli* expressing S-aminotransaminase (adjusted to 1.5 M AS)  
**Sample volume:** 2 ml  
**Buffer A:** 1.5 M ammonium sulfate, 50 mM sodium phosphate, pH 7.0  
**Buffer B:** 50 mM sodium phosphate, pH 7.0  
**Flow rate:** 1 ml/min  
**UV cell:** 10 mm  
**System:** ÄKTA pure equipped with column valve V9-C and loop valve V9-L



**Fig 6.** Column scouting for purification of S-aminotransaminase expressed in *E. coli*. Column selection valve V9-C allowed the connection of five HiTrap HIC columns to ÄKTA pure for this evaluation. UV monitor U9-M was used for multiwavelength detection. From this scouting, HiTrap Phenyl FF (high sub) 1 ml was selected for use in further scale-up studies.

## pH monitoring

An optional pH valve includes an integrated pH electrode that enables in-line pH monitoring during the run. The pH monitor is easily calibrated by injection of calibration buffer directly into the valve with the pH electrode mounted. A flow restrictor is connected to the pH valve and can be automatically included in the flow path to generate a back pressure that prevents the formation of air bubbles in the UV flow cell. The pH valve is used to direct the flow to the pH electrode and flow restrictor, or alternatively, to bypass one or both. Bypassing the pH electrode means that it can be stored and kept in place on the valve at all times.

## Outlet valves

Two different valve options are available to direct the flow to the fraction collector, waste, or other outlet ports. Outlet control valve V9-Os allows connection of a fraction collector and one outlet for collection of, for example, flowthrough. Outlet fractionation valve V9-O enables connection of more than one fraction collector and the 10 available outlets allow collection of large fractions.

## Fraction collection

ÄKTA pure can be equipped with the round Fraction Collector F9-R (Fig 7) or with the flexible Fraction Collector F9-C\* (Fig 8). Both fraction collectors are controlled through UNICORN control software. Fraction collection can be based on time, volume, or automatic peak recognition. Automatic peak recognition minimizes cross-contamination and unwanted eluent can be diverted to the waste. For increased capacity, two units of Fraction Collector F9-R or one F9-R and one F9-C can be connected together.

Fraction Collector F9-R provides a basic option with high capacity. A variety of racks is available to allow the use of 3, 8, 15, and 50 ml tubes. To minimize spillage, the Drop sync function can be used for flow rates up to 2 ml/min. Drop sync minimizes spillage by timing fraction changes between drops.



Fig 7. Fraction Collector F9-R allows collection in 3, 8, 15, or 50 ml tubes.

Fraction Collector F9-C\* provides flexibility, high capacity, and security. The fraction collector is equipped with a variety of cassettes that can hold tubes (3, 8, 15, and 50 ml) as well as deep well plates (24-, 48-, and 96-well), which means that samples can be collected in any format needed. Six cassettes can be loaded into the fraction collector in any combination that fits the user's needs (Fig 8). As an alternative to using 6 cassettes, loading capacity can be maximized by using a large tube rack for 50 ml tubes or a bottle rack for 250 ml bottles. Upon loading, the type of cassette is automatically detected by a sensor and the tube/bottle configuration is confirmed, eliminating mistakes in sample handling. Cassettes designed for tubes are equipped with a function that locks tubes into place when discarding liquid waste. Later, the tubes can be easily unlocked and discarded.

\* To be released.

The cassettes can also be used for convenient storage of fractions or holders for sample tubes and are easy to handle and clean. The fraction collector is covered, protecting samples from dust contamination. The top of the fraction collector can be used for placement of accessories and equipment.



Fig 8. Fraction Collector F9-C holds cassettes for a variety of tubes from 3 to 50 ml as well as 24-, 48-, and 96- deep-well plates.

Fraction collector F9-C has two beneficial features that minimize cross-contamination and spillage during fraction collection. Drop sync can be used for flow rates up to 2 ml/min and minimizes spillage by timing fraction changes to occur between drops. At higher flow rates, the Accumulator function temporarily holds the liquid flow during the time it takes to move to the next tube or well. The system can automatically change between the two modes for optimal performance.

## Additional module options

ÄKTA pure is a fully modular system that can be further expanded to increase system capability and productivity. Due to the accessibility and design of the modules, they are easily changed, which allows quick and efficient customization.

Versatile valve V9-V is a general four-position valve that can be used to tailor the flow system to specific tasks, for example, for multistep purification schemes.

Mixer bypass valve V9-M is used for bypassing the mixer if samples are loaded through the system pump.

Loop valve V9-L allows the use of up to five loops and can be used for collection of intermediate fractions when performing multistep purification or for automated purification of up to five different samples. The loop valve can, for example, also be used for holding reagents.

Up to two extra 8-position inlet valves can be deployed to expand on buffer and sample inlet capacity. Up to four additional air sensors can be placed in the flow path to enhance security, for example, before the inlet valves or before the column.

I/O-box E9 provides a means of connecting external interfacing equipment such as detectors. I/O-box E9 receives analog or digital signals from, or transfers analog or digital signals to external equipment that needs to be incorporated in the system. Two I/O-box E9 units can be connected to ÄKTA pure. A list of available additional valves and other options is found in Ordering information.

## UNICORN 6 control software

UNICORN control software is used on all ÄKTA systems and can be used from lab scale up to full-scale production. UNICORN contains the tools needed for beginners and experienced users to perform all types of chromatography, from setting up and running a method to evaluating the data.

Operation of ÄKTA pure is controlled by UNICORN 6 system control software, which provides built-in knowledge for planning, controlling, and analyzing the system and chromatography results. For greater convenience, licenses enabling remote control of ÄKTA pure are also available.

### Key features of UNICORN 6 are:

- Intuitive software design: simple, intuitive, and flexible method creation using predefined phases (steps)
- Database storage: robust data storage allows easy access to data, data security, and data integrity
- Interactive Process picture: shows the current open flow path and offers an intuitive way to control the run manually
- Column Logbook (optional): valuable tool to keep track of individual column and run data for traceability and operational security
- Design of Experiments (DoE, optional): integrated tool for experimental design provides more precise information in fewer experiments for cost- and time-efficient method development

### Method Editor

The Method Editor contains all the instructions used for controlling the chromatographic run. In UNICORN 6, the Method Editor has a user friendly graphical interface for easy viewing and editing of the method and run properties. Figure 9 shows a screenshot of the Method Editor with customizable panes that provide a comprehensive overview of the run.

In UNICORN 6, methods are built by using phases. Each phase reflects a step in the chromatography run, such as equilibration, sample application, or wash. The Method Editor contains predefined methods for different chromatography techniques and maintenance procedures, as well as a library of predefined phases for creating or editing your own methods.




Fig 9. Method Editor has customizable panes that give a comprehensive overview of the method.

A method is created or edited simply by dragging-and-dropping phases from the Phase Library into the Method Outline. In the Phase Properties pane, specific run parameters are set, and these settings are automatically programmed in the Text Instructions pane. When selecting a column type, parameter settings for that column type (e.g., default flow rates and pressure limits) are automatically set.

The Phase Properties pane allows for easy editing of the method, while more advanced users may program the method directly in the Text Instructions pane. For quick editing, the toolbar contains useful icons such as Undo, Redo, and Save.

### System Control

The System Control module is used to start, view, and control a method run. The module consists of three panes that provide an overview of the status of the run. The Run Data pane presents current data in numerical values, while the Chromatogram pane illustrates data as curves during the entire method run. The Process picture pane displays the current flow path during the run and can be used to control the run (Fig 10). Color indication incorporated in the Process picture shows the current open flow path with flow, closed flow path, or open flow path without flow. Real-time data from monitors are also displayed in the Process picture.




Fig 10. The Process picture in UNICORN 6 shows currently active and inactive flow paths, and provides a fast and easy way to control the system.


## Column Logbook

To increase operational safety, an optional feature of the software is the Column Logbook. This practical tool keeps track of important run data related to individual columns to provide traceability and operational security. Many prepacked columns from GE Healthcare are barcode labeled, and individual columns are identified using a 2-D barcode scanner, or the information may be entered manually into UNICORN 6. UniTag sticker labels, with preprinted 2-D matrix codes (barcodes), are available for labeling columns that are not already barcode labeled (e.g., empty columns).

By tracking individual columns, information is recorded for each run regarding run data such as total number of runs and maximum pressures. Notification limits can be set, for instance, to define the number of times the column may be run between cleanings, and the user is notified when it is time for column maintenance. In the Column Logbook, clicking on the Column History icon for a particular column provides a list of all runs that have been performed with that column.

## Design of Experiments (DoE) functionality

Design of experiments (DoE) functionality can be added as an option. DoE is a powerful tool for method optimization. DoE provides an efficient and structured approach where selected parameters are varied simultaneously such that a large data set can be obtained from few experiments (Fig 11). Since DoE is integrated seamlessly in UNICORN 6, scouting methods are automatically generated from DoE schemes, allowing for fast and efficient method optimization.



**Fig 11.** The Design of Experiments (DoE) module, which is optional in UNICORN 6, provides an efficient approach to optimization, capturing more information in fewer experiments.

## Accessories

ÄKTA pure accessories include column holders and clamps for attaching columns, flasks, and tubing to the system (Fig 12). Tubing kits in 0.25-, 0.5-, and 0.75-mm inner diameter dimensions allow optimization of the flow-path for various objectives and connection of any laboratory-scale column from GE Healthcare.



**Fig 12.** ÄKTA pure accessories include holders and clamps for attaching columns, flasks, and tubing to the system.

## Prepacked columns complete the package

GE Healthcare offers an extensive range of columns for purification from microgram to hundreds of milligram levels of target protein and for almost every chromatography technique (Fig 13). This range includes HiTrap, HiPrep™, HiScreen™, and HiLoad™, columns for preparative chromatography. Tricorn™ columns are also available for high-resolution semipreparative purifications at microgram scale as well as for protein characterization. In addition to prepacked columns, empty Tricorn, XK, and HiScale™ columns for packing with chromatography media of your choice are available.

### Columns for microgram-scale purification

Tricorn GL and PE columns are high-performance columns prepacked with media for a variety of chromatography techniques (Fig 13). The column design allows even distribution of liquid eluent over the entire column cross-section, which enables high-resolution purification at micro- and milligram scale. Tricorn GL columns are manufactured in glass to facilitate visual inspection of the media bed while the PEEK tube and filter of PE columns are designed to withstand greater pressure.

### Columns for milligram-scale purification


HiTrap 1 and 5 ml columns are prepacked with a wide range of media for purification using various chromatography techniques (Fig 13). The columns can be connected in series for greater capacity. Further scale-up can be achieved with HiPrep 20 ml columns.



RESOURCE™ columns are designed for high-resolution purification of proteins at high flow rates. The columns are prepacked with SOURCE™ media that have high particle size uniformity and stability to allow high flow rates at low back pressure.

HiScreen columns are prepacked with a wide range of robust BioProcess™ media to allow repeated use with highly reproducible results. Designed for scalable method optimization, HiScreen columns have a 10 cm bed height and can easily be connected in series to achieve a 20 cm bed height.

HiLoad columns are prepacked glass columns with Superdex prep grade media designed for high-resolution GF applications.



Tricorn	RESOURCE	HiTrap	HiScreen	HiPrep	HiLoad
High quality and high resolution (GF, IEX)	Fast purifications with good resolution (IEX, HIC, RPC)	Easy to use for a broad range of applications (AC, DS, IEX, HIC)	Optimized for method and process development (AC, IEX, HIC)	Fast and easy scale-up (AC, GF, IEX, HIC)	Preparative gel filtration with high resolution (GF)

**Fig 13.** Columns for use with ÄKTA pure system for different scales of purification. AC = affinity chromatography, DS = desalting, GF = gel filtration, IEX = ion-exchange chromatography, HIC = hydrophobic interaction chromatography, RPC = reversed-phase chromatography.

## System specifications

Control system	UNICORN 6 (version 6.3 or later)
Dimensions (W × H × D)	535 × 630 × 470 mm
Weight (excluding computer)	Up to 53 kg
Power supply	100–240 V, ~50–60 Hz
Power consumption	300 VA (typical), 25 VA (power-save)
Enclosure protective class	IP 21

## System pump

Pump type	Piston pump, metering type
Flow rate setting	0.001 to 25 ml/min (normal range) 0.01 to 50 ml/min (column packing flow)
Flow rate accuracy	± 1.2% (conditions: 0.25 to 25 ml/min, < 3 MPa, 0.8 to 2 cP)
Pressure range	0 to 20 MPa (2900 psi)
Viscosity range	0.35 to 10 cP (5 cP above 12.5 ml/min)

## Sample pump\*

Pump type	Piston pump, metering type
Flow rate setting	0.01 to 50 ml/min
Flow rate accuracy	± 2%
Pressure range	0 to 10 MPa (2900 psi)
Viscosity range	0.7 to 10 cP

## Mixer

Mixing principle	Chamber with a magnetic stirrer
Mixer volume	0.6, 1.4 (included with system), or 5 ml
Gradient flow rate range	0.1 to 25 ml/min
Gradient composition accuracy	± 0.6% (conditions: 5 to 95% B, 0.6 to 25 ml/min, 0.2 to 2 MPa, 0.8 to 2 cP)

## Valves

Type	Rotary valves
Number of valves	Up to 12
Functions	Standard: Injection valve Options: Inlet selection, mixer bypass, loop selection, column selection, pH, outlet, versatile

## Pressure sensors

Placement of sensors	Standard: After system pump Options: After sample pump, pre-column, post-column
Range	0 to 20 MPa (2900 psi)
Accuracy	± 0.02 MPa or ± 2%, whichever is greater

## Module options

### Inlet valves

Inlet A	1, 2, or 7 inlets
Inlet B	1, 2, or 7 inlets
Sample inlet	0, 1, or 7* inlets
Additional inlets	Up to 16

\* To be released

### UV monitors

Number of monitors	Up to 2 (two U9-L or one U9-L and one U9-M)
Wavelength range	UV monitor U9-L: 280 nm UV monitor U9-M: 190 to 700 nm in steps of 1 nm, up to 3 simultaneous wavelengths
Absorbance range	-6 to +6 AU
Linearity	Within $\pm 2\%$ at 0 to 2 AU
Operating pressure	0 to 2 MPa (290 psi)
Flow cells: U9-L	Standard: Optical path length 2 mm Cell volume 2 $\mu$ l Options: Optical path length 5 mm Cell volume 6 $\mu$ l
Flow cells: U9-M	Standard: Optical path length 2 mm Cell volume 2 $\mu$ l Options: Optical path length 10 mm Cell volume 8 $\mu$ l Optical path length 0.5 mm Cell volume 1 $\mu$ l

### Conductivity monitor

Conductivity reading range	0.01 mS/cm to 999.99 mS/cm
Accuracy	$\pm 0.01$ mS/cm or $\pm 2\%$ , whichever is greater (within 0.3 to 300 mS/cm)
Operating pressure	0 to 5 MPa (725 psi)
Flow cell volume	22 $\mu$ l
Temperature monitor range	0°C to 99°C
Temperature monitor accuracy	$\pm 1.5$ °C within 4°C to 45°C

### pH monitor

pH reading range	0 to 14
Accuracy	$\pm 0.1$ pH unit within pH 2 to 12,
Operating pressure	0 to 0.5 MPa (72 psi)
Flow cell volume	76 $\mu$ l

### Round fraction collector, F9-R

Number of fraction collectors	Up to 2 (Two Round fraction collector F9-R or one F9-R and one Flexible fraction collector, F9-C*)
Number of fractions	Up to 350
Vessel types	3, 8, 15, or 50 ml tubes
Fraction volumes	0.1 to 50 ml
Spillage-free mode	Drop sync
Flammable liquids	Yes
Delay volume (UV – dispenser head)	205 $\mu$ l (86 $\mu$ l with optional tubing kit)

### Flexible fraction collector, F9-C\*

Number of fraction collectors	Up to 2 (One Flexible fraction collector, F9-C and one Round fraction collector, F9-R)
Number of fractions	Up to 576
Vessel types	3, 8, 15, or 50 ml tubes; 24-, 48-, or 96- well deep-well plates; 250 ml bottles
Fraction volumes	0.1 to 250 ml
Spillage-free mode	Drop sync, Accumulator, or Automatic
Flammable liquids	Cleaning only (not fractionation)

### Outlet valves

Number of outlets	Valve V9-Os: 3 (waste, frac, 1 outlet position) Valve V9-O: 12 (waste, frac, 10 outlet positions)
Fraction volumes	0.01 to 100 000 ml
Delay volume (UV – outlet valve)	125 $\mu$ l (66 $\mu$ l with optional tubing kit)

### Air sensors

Number of sensors	Up to 7
Placement of in built-in sensors	Inlet valve V9-IA, Inlet valve V9-IB, Sample inlet valve V9-IS*
Sensing principle	Ultrasonic

### I/O-box E9

Number of ports	2 analog in, 2 analog out 4 digital in, 4 digital out
Analog range	In $\pm 2$ V Out $\pm 1$ V

\* To be released

## Ordering information

Product	Code number
ÄKTA pure L	29-0182-24
ÄKTA pure M	29-0182-26
ÄKTA pure L1	29-0182-25
ÄKTA pure M1	29-0182-27
ÄKTA pure M2	29-0182-28
ÄKTA pure System Handbook, printed copy (digital included)	29-2827-26
UNICORN 6.3 Workstation license with DVD	29-0209-27
UNICORN 6.3 remote license without DVD	29-0203-51
UNICORN 6.3 dry license without DVD	29-0203-55
UNICORN 6.3 DoE concurrent license	29-0203-66
UNICORN 6.3 Column logbook, network license	29-0203-84

### Standard system modules and accessories

Mixer	
Mixer chamber 0.6 ml	28-9561-86
Mixer chamber 1.4 ml (incl. with system)	28-9562-25
Mixer chamber 5 ml	28-9562-46
Online filter kit	18-1027-11

### Optional system modules and accessories

#### Valves/valve kits containing tubing

Inlet valve kit V9-IA	29-0122-63
Inlet valve kit V9-IB	29-0123-70
Inlet valve kit V9-IAB	29-0113-57
Inlet valve V9-X1	28-9572-27
Inlet valve V9-X2	28-9572-34
Mixer valve kit V9-M	29-0113-54
Loop valve V9-L	29-0113-58
Column valve V9-Cs	29-0113-55
Column valve kit V9-C	29-0113-67
pH valve kit V9-pH	29-0113-59
Versatile valve V9-V	29-0113-53
Outlet valve kit V9-O	29-0122-61
Outlet valve kit V9-Os (1 outlet)	29-0113-56

#### UV monitor

UV monitor U9-L	29-0113-60
UV flow cell U9-0.5, 0.5 mm for U9-M	28-9793-86
UV flow cell U9-2, 2 mm for U9-M (incl. in system with U9-M)	28-9793-80
UV flow cell U9-10, 10 mm for U9-M	28-9563-78
UV flow cell 2 mm for U9-L (incl. with U9-L)	29-0113-25
UV flow cell 5 mm for U9-L	18-1128-24

#### Fraction collector

Fraction collector F9-R	29-0113-62
Tube Rack with 175 positions for 12 mm vials, bowl, tube support, holder and guide	19-8684-03
Tube Rack with, 95 positions for 10-18 mm vials	18-3050-03
Tube Rack with 40 positions for 30 mm vials, bowl, tube support, holder and guide	18-1124-67

#### Additional air sensors

Air sensor L9-1.2 mm	28-9565-02
Air sensor L9-1.5 mm	28-9565-00
Adapter for air sensor	28-9563-42

#### I/O-box E9

I/O-box E9	29-0013-61
------------	------------

#### Tubing kits

Tubing Kit 0.5 standard	29-0113-27
Tubing Kit 0.25	29-0113-28
Tubing Kit 0.75	29-0113-29

#### Barcode labels and scanner

UniTag (1 sheet with 108 labels)	28-9564-91
Barcode scanner 2-D with USB	28-9564-52

#### Holders

Column holder rod	28-9562-70
Tubing holder spool	28-9562-74
Column holder	28-9562-82
Tubing holder comb	28-9562-86
Flexible column holder	28-9562-95
Column clamp o.d. 10-21 mm	28-9563-19
Bottle holder	28-9563-27
Tube holder (5-pack)	28-9543-29
Multi-purpose holder	29-0113-49
Rail extension	29-0113-52
Loop holder	29-0113-50

### Related literature

	Code number
UNICORN 6 control software, Data file	28-9573-46
Validation Support File UNICORN software, Data file	28-9626-50
ÄKTA pure: Protein purification your way, Brochure	29-0211-97
Purification of a miniature recombinant spidroin protein expressed in <i>E. coli</i> using ÄKTA pure system, Application note	29-0211-98
Purification and immobilization of a transaminase for the preparation of an enzyme bioreactor, Application note	29-0211-99
Prepacked chromatography columns for ÄKTA systems, Selection guide	28-9317-78

For local office contact information, visit  
[www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

[www.gelifesciences.com/aktapure](http://www.gelifesciences.com/aktapure)

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden



GE, imagination at work, and GE monogram are trademarks of General Electric Company.

ÅKTA, BioProcess, HiLoad, HiPrep, HiScale, HiTrap, HiTrap, RESOURCE, SOURCE, Superdex, Superloop, Tricorn, and UNICORN are trademarks of GE Healthcare companies.

Any use of UNICORN software is subject to GE Healthcare Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

© 2012 General Electric Company—All rights reserved.  
First published Aug. 2012.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Limited  
Amersham Place  
Little Chalfont  
Buckinghamshire, HP7 9NA  
UK

GE Healthcare Europe, GmbH  
Munzinger Strasse 5  
D-79111 Freiburg  
Germany

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327  
Piscataway, NJ 08855-1327  
USA

GE Healthcare Japan Corporation  
Sanken Bldg., 3-25-1, Hyakunincho  
Shinjuku-ku, Tokyo 169-0073  
Japan