CERTOMAT® IS

Inkubations-Schüttler

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Operating Manual

Mode d'emploi

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Operating Manual

Introductory Note

This Operating Manual describes the CERTOMAT® IS at the time this document was printed (see rev.-no. in the bottom line). Sartorius Stedim Biotech GmbH reserves the right to modify the equipment of its products as well as the order information without notice.

The information given in this document has been carefully estimated. However, this document is not necessarily complete. Due to the continuous development of the unit it may be possible that you find information which does not apply to your device. If you find missing or wrong information or if you need more detailed information concerning any of the features of the unit please send us a note.

The CERTOMAT® IS is just an example of laboratory equipment products from the Sartorius Stedim Biotech GmbH.

Our range of equipment includes shakers and incubator shakers, centrifuges, ultrapure water systems, air samplers and fermenters/bioreactors.

You will find further information on our homepage www.sartorius-stedim.com

If you are interested in more information about the complete product program or about specific equipment, please contact us.

Sartorius Stedim Biotech GmbH
Weender Landstraße 94 - 108
37075 Goettingen, Germany
Telephone +49.551.308.3118
Telefax +49.551.308.3918
info@sartorius-stedim.com
www.sartorius-stedim.com

EG-Declaration of Conformity

With the declaration of conformity, which is attached or which accompanies the equipment, Sartorius Stedim Biotech GmbH confirms the conformance of the CERTOMAT® IS with the specified guidelines.

General Safety Notes

- Laboratory seniors or the persons responsible for the use of the equipment must inform the personnel about the dangers at the operation of the CERTOMAT® IS. Furthermore they have to guarantee that the relevant safety notes are observed. The personnel must have the necessary qualification for the operation of the equipment and its handling.
- Please take special care regarding the operating information. Improper use of the unit is thus avoided and you have full claim for guarantee. Please contact us if you have questions.

Important safety notes in this manual are marked with this symbol. Such instructions apply additionally to the other safety instructions for the field of application and at the working place.

- Do not use devices with the wrong potential version.
- Disregard of the safety instructions may for example result in damages at the equipment or can cause other material and eventually also person damages.

Contents

Introductory Note

EG Declaration of Conformity

General Safety Notes

1	Product Information	4
1.1	Applications	4
1.2	Design Features	4
1.3	Available Versions and Accessories	5
1.3.1	Trays	5
1.3.2	Illumination Unit	5
1.3.3	Alarm contact	5
1.3.4	Analog exit	5
2	Delivery and Installation	6
2.1	Unpacking and Checking of the Delivered Parts	6
2.2	Installation of the CERTOMAT® IS	6
2.3	Transport Protection	6
3	Operating Information	7
3.1	Safety Notes	7
3.2	Preparing the Tray and Mounting the Assemblies	8
3.2.1	Preparing the Tray	8
3.2.2	Mounting the Required Tray	8
3.3	Commissioning	9
3.3.1	Mounting and set-ups	9
3.3.2	Function and Menu Overview	11
3.3.3	Manual Menu "MAN"	11
3.3.4	Program menu "PROG"	15
3.3.5	Setup menu "SETUP1"	25
3.3.6	Set-up menu "SETUP2"	31
3.4	Connection of External Devices	33
3.4.1	Analog Exits (ANALOG OUT)	33
3.4.2	RS 232 Interface	33
3.4.3	Collective Alarm (optionally)	33
3.5	Maintenance and Cleaning	34
3.5.1	Cleaning	34
3.5.2	Malfunction Correction	34
3.6	Warranty Regulations and Service	35
4	Technical Data, Ordering Information	36
4.1	Technical Data	36
4.1.1	Mechanical Design	36
4.1.2	Mains Supply	36
4.1.3	Operating Data	36
4.1.4	Analog Output	36
4.1.5	Collective alarm	36
4.2	Ordering Information	37
4.2.1	Configuration of Shakers	37
4.2.2	Accessories	37
4.3	Safety Instructions Laboratory Shakers	37
4.4	EG Declaration of Conformity	37
4.5	Declaration of Decontamination	37
4.6	Information and Instructions on Disposal and Repairs	37
4.7	Dimension sheet CERTOMAT® IS	37
4.8	Prospectus CERTOMAT® IS	37

1 Product Information

1.1 Applications

The CERTOMAT® IS is an easy to operate, efficient and durable Incubation Shaker for the universal use in the various tasks of biological and chemical laboratories.

The basic equipment of the CERTOMAT® IS includes shaking drive and integrated heating. The shaking drive includes a brushless external rotor motor and the noiseless, robust belt-drive. Shaking drives with different lifting capacity and for different shaking speeds are available, which can be equipped with trays of different size as well as with mounting systems and the accessories of the CERTOMAT® shaker program.

A plexiglas cap which can be opened upwards and is led over a gas pressure spring makes the assembly and the disassembly of shaking boxes or the changing of the tray easier and makes the visual control possible. For the operation from 5 °C over room temperature and below it a refrigerating set must be integrated. An illumination unit is optionally available. For more information see chapter 4, „Technical Data and Ordering Information“. The customer can thus use the cabinet according to his specific requirements.

The bottom is designed as a spillage tray made of stainless steel. If containers get damaged during shaking operation this prevents spillage from escaping unchecked. The spillage tray also facilitates the cleaning of the shaker.

1.2 Design Features

- Counter weight for mass balancing of the load applied; vibrations, which can be transmitted to the shaker unit and working place, i.e. due to resonance effects, will be minimized.
- Shaking amplitudes of 25 mm and 50 mm; delivery setting of shaking amplitude can be changed to the alternate setting in-place at customer's site by authorized service, if necessary.
- Especially designed mounting of the tray support for tension-free operation; providing precisely reproducible rotation at minimum vibrations even at heavy loads.
- New, simple and reliable tray locking.
- Motor of external rotor type with brushless drive and robust, noiseless Poly-V belt drive; providing noiseless and long-term, maintenance free operation.
- Alphanumeric display with menu-controlled programming.
- Visualization in different languages.
- Microprocessor-controlled motor and temperature regulation.
- MEMORY-function at power failure.
- Memorization of alarms which occurred, according to kind and duration.
- Analog exit 0 ... 10 V for the documentation of temperature and speed development.
- Serial interface RS 232 for documentation of temperature and shaking speed development.

1.3 Available Versions and Accessories

1.3.1 Trays

The CERTOMAT® IS can be used with trays of the following size:

- Trays: type E/EU (420 x 420 mm)¹⁾
- Trays can be delivered being finally equipped with sets of clamps for Erlenmeyer flasks of different size, universal mounting sets or slew able test tube racks. The universal mounting sets allow the assembly of various vessels, such as bottles, beakers or also separation funnels. Further information you will find in the separate documentation “trays, mounting systems and accessories for the CERTOMAT® -shaker.

1.3.2 Illumination Unit

- As further option a lighting device with 2 fluorescent lamps of 8 W each is available.

1.3.3 Alarm contact

- As further option a contact without potential is available for the external alarm message.

1.3.4 Analog exit

Optionally the analog exit can be provided also as current exit 0 ... 20mA / 4 ... 20 mA.

¹⁾ Trays of modified design for new locking; rebuilding of existing trays for this locking on request. Please contact the Sartorius AG customer service.

2 Delivery and Installation

2.1 Unpacking and Checking of the Delivered Parts

1. For the transport to the working place please note the weight of the shaker and the accessories included in the delivery. Use only suitable transport means.
2. During unpacking carefully check whether the delivery is complete, according to your order. Please check carefully all parts for damages, that may be caused by transport.
3. The mains supply at the working place must correspond to the mains specifications of the shaker's type plate and must have an earthed conductor. Admissible mains supplies are 230 V, 50 Hz, 1.2 kW for instance. Check on the type plate whether the shaker has the correct power ratings.
 - Modifications of the CERTOMAT® IS for different mains supplies or changes of the mains cable according to local standards will be done and checked in the factory or by the Sartorius Stedim Biotech GmbH customer service in place. The unit will be labeled correspondingly.
4. If the CERTOMAT® IS delivered to you has the wrong mains specifications, if any parts are damaged or are missing or in case of malfunctioning, you should return your claim note immediately to your Sartorius Stedim Biotech GmbH representative or directly to Sartorius Stedim Biotech GmbH.

2.2 Installation of the CERTOMAT® IS

1. You can place the CERTOMAT® IS on any sufficiently big and stable laboratory working place, e. g. common laboratory tables. Please observe the dimensions, the empty weight and the intended load. There must be sufficient room for a convenient operation.

Dimensions: W x H x D = 540 x 564 x 680 mm

Weight: 75 kg without loading
(heating version)

Load: depending on type of tray,
 mounting system and shaking
 vessels

2. Make sure, that the base is stable enough to carry the weight of the completely equipped shaker and can resist vibrations caused by its operation.

The tray action causes vibrations at the shaker, which cannot always be compensated by the internal balance weight as regards loads and shaking speeds. Especially at heavy loads (> 10 kg) and high shaking speeds the vibrations can be transmitted to the working place.

3. The working place should have an antiskid surface. For stable standing of the shaker during operation you may clean the table and the stands with a wet cloth before starting operation.
4. After placement of the shaker at the working place carefully adjust the unit horizontally. An unbalanced placement can result into faster wear and tear of the internal drive, especially if the unit is often applied with heavy loads and/or at high shaking speeds.
5. Place equipment around in such a way that it cannot be interfered by the shaker's action.

2.3 Transport Protection

The CERTOMAT® IS will be delivered without special transport protection. The device is ready for use after unpacking, connecting to the power supply and equipping with the intended trays and mountings.

3 Operating Information

3.1 Safety Notes

Media used with incubation shakers like the CERTOMAT® IS can be hazardous. This is especially true for burning, inflammable or infectious media, for instance. Since the risks and the correlated safety requirements depend on the media, they cannot be discussed herein.

our company must issue applicable safety guidelines, distribute them to the laboratory staff and ensure that they are strictly followed.

The following safety instructions have to be observed:

Please compare the existing supply voltage with the type plate on the device before you start it. If the device does not correspond with your supply voltage you cannot start it.

Make sure that the mounting place (laboratory table etc.) can bear the weight of the shaker under load and can compensate the vibrations which are produced during the operation.

Do only use clamps and racks which hold safely the respective shaking vessels at all possible operation modes of the shaker. Check the position of the vessels before you start the shaking operation. The shaking vessels must not hit against each other during operation since they might be damaged.

Please always adjust the shaker to the STOP position before you put vessels on the shaker table or in the vessel racks or take them out.

If vessels break during the operation or if media are spilled, please stop the shaking operation immediately, remove broken vessel parts and clean the unit.

Please care for a safe standing during the operation. Stop the shaker if it moves over the mounting place or if vibrations influence the mounting place in such a way that other devices are disturbed. In such case you should reduce the shaking speed.

If vibrations cannot be compensated and there is the danger that the shaker moves over the mounting place you have to check regularly if it is still in a safe operation state.

Only personnel which is qualified and authorized to use the shaker has access to the working place and is allowed to operate the unit.

It is generally recommended to use suited work wear and personal safety equipment, such as gloves, safety glasses and eventually a gas mask.

After unpacking or a change of the mounting place (e. g. cooling room <-> laboratory) you should wait until the unit has warmed up to room temperature before you start it. Otherwise humidity may condense in the unit and eventually cause function errors.

Only operate the shaker at the following ambience conditions:

Ambient temperature: + 10 ... + 60 °C

Ambient temperature at

integrated cooling unit: + 10 ... + 30 °C

Relative air humidity: + 10 ... + 90 %, not condensing.

The CERTOMAT® IS is preadjusted for a load of 10 kg if it is delivered. The internally mounted balance weight corresponds with the respective shaking amplitude of the delivery state.

If you intend to use the shaker at the alternative shaking amplitude please contact our service concerning the reconstruction of the drive unit and the balance weight.

3.2 Preparing the Tray and Mounting the Assemblies

3.2.1 Preparing the Tray

1. Trays can be delivered finally equipped with sets of clamps for Erlenmeyer flasks of different size or as universal tray without assemblies. If you are going to use the existing equipment of the tray continue with step 3.
2. In order to change the equipment of the tray you have to release the parts you do not need and install the accessories necessary for your purposes. Take care at the arrangement of the mountings that the vessels can be distributed evenly if you load the tray.
- Clamps for Erlenmeyer flasks of different size can be individually combined, depending on the type of tray and the prepared borings. For other vessels (such as beakers, bottles, separation funnels, for instance) special mounting sets are available. For details see chapter 4 "Ordering information".
3. Place the shaking vessels into their corresponding holders. We recommend to place single vessels, etc., in the centre of the tray. Place other vessels evenly around from the centre outwards to avoid unbalances.
4. Since the mountings must properly locate the intended vessels at any time during operation under any of the intended operating conditions, they should be well suited to the vessels.

 The shaking containers must not hit against each other during the shaking operation. They must not get loose unintentionally.

3.2.2 Mounting the Required Tray


Fig. 1: Mounting a universal tray to the tray support (the tray is shown without mounting systems and containers)

1. Place the tray onto the shaking table. Push it with the indents (1a/1b) against the rear knurled screws (1c/1d) and the clamp (2a) of the shaking table.

 Please observe the weight of the tray if you have already loaded it with shaking vessels.

2. Push down (2b) the front edge of the tray so that it properly fits into the front clamp (2c).
3. The borings (3a/3b) in the tray will precisely fit onto the locating pins of the tray support (3c/3d) and locate the tray.

 Please take care that the locating pins (3c/3d) are placed evenly in the tray. Only then the tray is situated correctly on the shaker table.

 If necessary you can adjust the locating pins as follows:

1. Release the screws of the locating pins (3c/3d) with the hexagon socket screw key.
2. Eventually agitate the tray slightly. The locating pins (3c/3d) locate themselves in the centre of the borings (3a/3b) of the tray.
3. Screw down the screws of the locating pins (3c/3d) carefully again after the tolerance balancing.
4. In order to improve the safe position of the tray – especially at high loads in connection with high speeds – the provided flat headed screws M5 x 8 must be screwed in the prepared borings.

3.3 Commissioning

3.3.1 Mounting and set-ups

1. If you have changed the mounting place (e. g. cooling room <-> laboratory) please wait until the device has warmed up to room temperature before you start it. Otherwise the air humidity can condensate in the device and can lead to functional disturbances.

Following the transport of a CERTOMAT® IS with integrated cooling unit, the device may only be put into operation after a waiting period of an hour. This avoids danger to the cooling unit.

2. Please connect the CERTOMAT® IS with the mains cable to a suited socket.
3. By pressing the mains switch (1) the electronics are supplied with line voltage.

After switching on the device initializes. After approx. 5 sec. the background illumination of the display as well as the green operating LED light. The manual menu „MAN“ is active. The device is now ready for operation.

In pressing the ESC-key you change from the manual menu to the main menu.


Fig. 2: Frontpanel with key pad

1. mains switch
2. alphanumeric display
3. function field [F1-F5]
4. [+/-]-key for the set-up of process values (temperature, speed)
5. [ESC] -key for menu jump back
[ENTER] -key for confirmation of inputs / changes
6. Operation-LED (green)
7. [START] -key
8. ALARM-LED (red)
9. [STOP] -key
10. display field MODE

Checking the door contact switch

 Whenever you start to operate the CERTOMAT® IS, you have to check the proper function of the door contact switch!

 To check the function of the switch, proceed as follows:

 After turning on the CERTOMAT® IS and initialisation, the unit is in the Manual-Menu „MAN“.

- In the manual mode „MAN“ enter these values via the keypad: (see also chapter 3.3.3)

F1 speed [l/min.]: 50

F3 time [h : min]: 00:10

- Press the START-key to start the shaker.

 Open the hood of the CERTOMAT® IS. The CERTOMAT® IS must stop.

 Once the preset time has expired (see chapter 3.3.5), the alarm LED on the front panel lights up. If activated, an alarm signal is generated in addition.

- Press ENTER to acknowledge the alarm.
- Close the hood of the CERTOMAT® IS. The shaker re-starts.

 If the shaker continues to run inspite of the openend hood, the unit is defective. In this case the CERTOMAT® IS must not be started!

Repairs can be carried out by the Sartorius AG service representative or by authorized personnel of a customer-owned workshop. Please inform your local Sartorius AG sales representative or contact:

Sartorius AG
Service-Zentrum Nord
Weender Landstraße 94-108
37075 Goettingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Telefax +49.551.308.3730

3.3.2 Function and Menu Overview

In the menu overview you will find the menu points which are available in the main menu. In the chapters described in the following you will find the set-ups of the respectively selected menu item.


Fig. 3: Main menu (mene overview)

3.3.3 Manual Menu "MAN"

- In activating the key [F4] "MAN" in the main menu the "manual menu" is activated.
- Back jump from the "manual menu" to the main menu with [ESC].


Fig. 4: main menu

Abb. 5: MAN: manual menu

The following set-ups in the "manual menu" (fig. 5) can be carried out:

F1 **SPEED:** **Set-up of speed**

- After pressing of [F1] the input field "SPEED" in speed/min. flashes.
- In pressing the [+/-] -keys the speed is adjusted.
- Confirm input with [ENTER].

F2 **TEMP:** **Set-up of the temperature**

- After pressing of [F2] the input field "TEMPERATURE" flashes in °C.
- In pressing the [+/-] -keys the temperature is adjusted.
- Confirm input with [ENTER]

F3**TIME:****Set-up of the TIMER**

- After pressing of [F3] the “TIMER“ [TIME] can be adjusted.
- In pressing the keys [+/-] the desired operating time in [h : min] is adjusted.
- Confirm input with [ENTER].
- The operating time can be set between 00h:01min and 99h:59min as well as “CONT.“ (continuous operation).

F4**COOLING:****Switching the cooling on and off (only at “UHK version “)**

- After pressing of [F4] the cooling is switched on / off.

Switching the cooling on is only possible at temperature set values < 38.5°C!

MODE**MODE:****Display field**

- The respective “active“ menu is indicated, here: „MAN“ for manual menu.

- With the [F5] key you can mutually change to another display field in order to indicate all values to be set (only at option “illumination”).


Fig. 6: adjustable parameters

< F5 >


Fig. 7: further set-up possibilities

F2

LIGHT:

**Switching the illumination on and off
(only at option “illumination”)**

- After pressing of [F2] the illumination is switched on / off.

MODE

MODE:

Display field

- The respective “active” menu is indicated, here “MAN“ for manual menu.
- In pressing the key [START] you change to the operating mode “RUN“.
- The previously adjusted values for:

Speed:	SPEED
Temperature:	TEMP
Lighting:	LIGHT
Cooling:	COOLING
Timer:	TIME

are activated now.

- With the key [F5] you can alternately change to another display field in order to indicate all actual values.


Fig. 8: Operating mode RUN

< F5 >


Fig. 9: Operating mode RUN

- The "TIMER" runs. The respective rest operating time is indicated.

In the operating mode the ACTUAL values are indicated.
A change of the values is not possible here!

- In pressing the function keys [SPEED], [TEMP], [COOLING] and [LIGHT] a little bit longer the parameter flashes and the SET-value is indicated. After the release of the function keys the ACTUAL-values are indicated again.
- The ACTUAL value display of the temperature refers to the air temperature in the incubation room.
- In pressing the key [STOP] you jump back to the manual menu "MAN".
- Heating, cooling, light and shaker are switched off.
- A new set-up can be effected.
- After the adjusted operating time is finished, you jump automatically back to the manual menu "MAN".
- After pressing of the [STOP]-key the cooling is still running for 3 minutes. Switching the cooling on again is only possible after 3 minutes.

3.3.4 Program menu "PROG"

- In activating the key [F5] "PROG" the "program menu" is activated. "PROG1" flashes!
- Jump back to the main menu with [ESC].


Fig. 10: main menu


Fig. 11: PROG1-4: program menu

- With [F1] [F5] the respective program is selected. After the selection the selected program flashes in the display.
- In the "MODE"-field the selected program appears automatically, here "PROG1".
- The programming according to an example is explained in point 3.3.4.2 "Example for programming".

3.3.4.1 Program definition

- After pressing of the [ENTER] key (alternately press the flashing function key for approx. 2 sec.) you change to the program definitions “PRSTEP . .STEP1 . .STEP4”
- In this menu the different “STEPS“ can be selected within the program “PROG1.....PROG4“.
- At programs protected by passwords the respective password must be input for the editing/modification of present programs.


Fig. 12: PROG1-4: program menu


Fig. 13: program definition

- In fig. 13 (program definition) is indicated, whether the PRSTEP is “ON” or “OFF”. Furthermore the step duration is indicated.

Selection of the different “STEPS“ (fig. 13):

F1 for: **PRSTEP**

- After pressing of [F1] the set-ups for “PRSTEP“ can be effected.

F2 for: **STEP1**

- After pressing of [F2] the set-ups for “STEP1“ can be effected.

F3 for: **STEP2**

- After pressing of [F3] the set-ups for “STEP2“ can be effected.

F4 for: **STEP3**

- After pressing of [F4] the set-ups for “STEP3“ can be effected.

F5 for: **STEP4**

- After pressing of [F5] the set-ups for “STEP4“ can be effected.

MODE **MODE:** **Display field**

- The respective “active” PROG is indicated, here “PROG1”.

With the keys [F1 to F5] the adjustable parameters of the steps can be selected.


Fig. 14: program definition



Abb. 15: adjustable parameters

For each “STEP“ the following parameters can be adjusted (fig. 15):

F1

SPEED:

Set-up of the speed

- After pressing of [F1] the input field “SPEED“ in speed/min. flashes.
- In pressing the keys [+/-] the speed is adjusted.
- Confirm input with [ENTER].

F2

TEMP:

Set-up of the temperature

- After pressing of [F2] the input field “TEMPERATURE“ in °C flashes.
- After pressing of the keys [+/-] the temperature is adjusted.
- Confirm input with [ENTER].

F3

TIME:

Set-up of the STEP-duration [h : min]

- After pressing of [F3] the input field “TIME“ flashes.
- In pressing the keys [+/-] the time intervall (duration) is adjusted.
- Confirm input with [ENTER].
- At the STEP duration values between 00h:01min (immediate change to the next step) and 99h:59min as well as “CONT.“ (continuous operation until the manual break with STOP) can be adjusted.

You cannot define a firm operating time for the PRSTEP. It can only be switched “ON“ or “OFF“!

- PRSTEP in operating mode: the conditions selected in the PRSTEP from the pressing of the [START] key until the actual program start [STEP1] are activated.
- PRSTEP not operating: the shaker, tempering, cooling and light are deactivated until the program start [STEP1].

F4

COOLING:

**Switching the cooling on and off
(only at "UHK version")**

- After pressing of [F4] the cooling is switched on / off.

☞ Switching the cooling on is only possible at temperature set values < 38.5°C!

MODE

MODE:

Display field

- The respective "active" STEP is indicated.
- With the [F5] key you can alternately change to another display field in order to set up further parameters (only at option "illumination").


Fig. 16: adjustable parameters

Fig. 17: further set-up possibilities

F2

LIGHT:

**Switching the illumination on and off
(only at "option illumination")**

- After pressing of [F2] the illumination is switched on / off.

MODE

MODE:

Display field

- The respective "active" STEP is indicated, here "STEP1"

- In pressing the key [ESC] you jump back to the previous menus (until the main menu)


Fig. 18: set-up possibilities within a "STEP"


Fig. 19: program definition


Fig. 20: program menu


Fig. 21: main menu

3.3.4.2 Example for programming

- Altogether five different programs can be memorized in the device. Each program consists of one "PRSTEP" and four further steps which can be processed one after the other. For each step the values for cooling, light, speed and temperature can be set independently from one another. Each step has a defined duration time.
- If only one or two steps are required the remaining steps can be defined with a duration time 00:00 so that the program is directly stopped afterwards. A repeat function for endless loops is available.
- The "PRSTEP" can only be switched on or off (no operating time). It defines the conditions (cooling, light, speed and temperature) from the pressing of the key [START] until the program start (STEP1).
- At "PRSTEP" switched off, the cooling, light, speed and temperature are switched off until the program start.

- The main purpose of use for the "PRSTEP" is the date/time exactness for the provision of cultures, e. g. after the weekend or the next morning.

Advantages:

- You have to define a program only once with operating times and you can start it afterwards at any time. The time of the program end always remains the same.
- Comfortable start of night or weekend programs during the day or before the weekend.
- Conditions are reproducible before the start of the program because these are held by the device according to the settings in the "PRSTEP".
- Automatic program start without the personnel being present.

Example:

On Monday morning at 9:00 o'clock any cultures have to be ready. The program (STEP1...STEP4) for the production of these cultures needs altogether 38 h. So the program would have to be started manually without use of the "PRSTEP" on Saturday evening at 7:00 o'clock..

The "PRSTEP" makes it possible to choose any time for the start and to keep the material until the program start in the desired conditions.

PROG1: PRSTEP: ON at 20°C, speed 0, light off, cooling on

Step 1: 7 h at 30° C, speed 200, light on, cooling off

Step 2: 10 h at 40° C, speed 100, light off, cooling off

Step 3: 20 h at 10° C, speed 60, light on, cooling on

Step 4: 1 h at 10° C, speed 200, light on, cooling on

Function	PRSTEP	STEP1	STEP2	STEP3	STEP4
Cooling	on	off	off	on	on
Light	off	on	off	on	on
Speed [rpm]	0	200	100	60	200
Temperature [°C]	20	30	40	10	10
TIME (operating time)	ON	7 h	10 h	20 h	1h
Day	Saturday	Sunday	Sunday	Monday	Monday
Time	19:00	02:00	12:00	08:00	09:00
	START PRSTEP	Program START			Program end

Fig. 22: Example for programming

3.3.4.3 Password protection for programs

- The “password protection“ (code) for programs serves for the protection of programs which have been once established for this purpose against unauthorized access. The content of the program cannot be modified without knowing the password. Optionally also the “starting“ and „stopping“ of the programs can be protected by this password.
- After pressing of the [ENTER] key you change to the menu where you input the “password“ (code) for the respective active program.
- For safety reasons the “password“ must be input twice if the password is determined.


Fig. 23: program definition of PROG1


Fig. 24: determination of password for PROG1

- In pressing the function keys [F1] ... [F5], [+], [-], [START] or [STOP] as admissible signs the “password“ (code) (visible “★“ on the display) is given in.
- The “password“ (code) must have at least 4 and at most 7 places.
- After the input of the password the “password“ (code) must be confirmed with the [ENTER] key. Then you change automatically to the confirmation of the password which has been input.
- After the repeated input of the password the “password“ (Code) must be confirmed once again with the [ENTER] key in order to secure the “password input“. At the same time the back jump to the “program definition menu“ is effected.
- In pressing the [ESC] key the “password determination“ can be stopped. All “password inputs“ which have been effected until this time become invalid.

3.3.4.4 Cancelling of a program password

- Select in the program definition menu the program in which you want to cancel the password (password knowledge is necessary).
- In pressing the [ENTER] key you change from the program definition menu to the menu for the password determination.


Fig. 25: program definition of PROG1


Fig. 26: determination of password for PROG1

- Deleting / deactivating of the “password” (code) is carried out with the [ENTER] key as first password sign. You change then to the confirmation of the deletion / deactivation.
- The repeated input of [ENTER] causes the deletion of the existing “password” (code) and the automatic back jump to the program definition menu is effected.
- In pressing the [ESC] key the deletion / deactivation of the “password“ (code) can be stopped.

3.3.4.5 Starting of a program

- In pressing the [START] key afterwards you change to the setup mode of the program selected before, here program 1 "PROG1".


Fig. 27: PROG1-4: program menu


Fig. 28: start menu

The following parameters in the start menu (fig. 28) can be adjusted:

TIME:

Set-up of the time (start time of the program)

- After pressing of [F1] the input field "TIME" flashes.
- In pressing the [+/-]-key the time is adjusted.
- Confirm input with [ENTER].

DAY:

Set-up of the date (start date of the program)


- After pressing of [F2] the input field "DAY" flashes.
- After pressing of the [+/-] keys the date is adjusted.
- The current date and the current time are respectively preadjusted.
- Confirm input with [ENTER].

REPEAT:

Set-up of repeat

- After pressing of [F3] the automatic program repeat is activated / deactivated.
- At activated program repeat the "STEP 1" to "STEP 5" of the program is processed in an endless loop.
- The "PRSTEP" is cycled only at the first time.

- The respective “active” program is indicated, here “PROG1“.
- In pressing the [START] key afterwards you change to the operating mode “RUN“. The program activated before is processed.
- With the [F5] key you can alternately change to another display field in order to indicate all actual values!


Fig. 29: RUN: operating mode





Fig. 30: RUN: operating mode

- At the TIME-display the rest operating time of a “STEP” is indicated. During the “PRSTEP” processing no time is indicated (display “PRSTEP“).
- In the operating mode the ACTUAL VALUES are indicated and the started program runs automatically.
- In pressing the respective function keys the SET VALUE is indicated.
- After having released the function key the ACTUAL VALUES are indicated again.
- In pressing the [STOP] key you jump back to the program menu “PROG“. Heating, cooling, light and shaking drive are switched off. The program processing is stopped.

During the program processing it is indicated which program and which „STEP“ are being processed at the moment, e.g.:

- For the programs “PROG2-PROG5“ the set-ups are used analogously to the description for program 1 “PROG1“.

3.3.5 Setup menu "SETUP1"

- In pressing the [F3] key you change to the "setup 1 menu".
- All set-ups and modifications are memorized internally and remain active also after voltage loss or switching off and on!


Fig. 31: main menu


Fig. 32: setup1 menu

Display contrast

- Can be modified in pressing the keys [+] and [-].

Attention:

In pressing the key [-] the contrast is so much reduced that no information is visible at the display!

The following parameters in the "setup1 menu" (fig. 32) can be adjusted:

F1 LANG.: Set-up of the menu language

- After pressing of the [F1] key you come to the menu in which you adjust the menu language (German, English, French).
- An arrow marks the respectively active language.
- In pressing the keys [+] / [-] the arrow can be moved.
- In pressing the [ENTER] key the respectively adjusted language is confirmed.
- In pressing the [ESC] key you change again to the "setup1-Menu". The language set-up before is now shown in the menus.

F2**DIAG:****Diagnosis of the alarm messages which occurred**

- The alarms which occurred are indicated on a list. The following information is separately disposed for each error type:

Error type	e.g. speed or temperature alarm
Error counter	Number of the alarms which occurred since the last acknowledgment
Date / time / alarm start	Date and time of the alarm start (only of the alarms which occurred last)
Date / time / alarm end	Date and time of the alarm end (only of the alarms which occurred last)

- An arrow marks the respectively active list entry.
- In pressing the keys [+/-] the arrow can be moved.
- In pressing the [ENTER] key the respectively active list entry is cancelled. Afterwards the alarm information for this kind of error is not available any more.
- If no alarm occurred or if all list entries were cancelled [“....<ESC>....“] appears on the display. The diagnosis menu can only be quit with [ESC].

F3**ALERTS:****Alarm definition**

In a list the alarm parameters for temperature, speed, door contact and voltage loss can be adjusted.

- An arrow marks the respectively active parameter.
- In pressing the keys [+/-] the arrow can be moved.
- In pressing the [ENTER] key the YES / NO definitions are directly changed or the numeric values can be released for modification (value flashes).
- In pressing the keys [+/-] the respectively flashing value can be set.
- In pressing the [ENTER] key the adjusted value is confirmed.
- In pressing the [ESC]-key you jump back to the SETUP1-menu.

The following alarm parameters can be set up:

Temperature 1	ON	The temperature alarm is active
	OFF	The temperature alarm is inactive
Signal	ON	Acoustic signal active in case of an alarm
	OFF	Acoustic signal inactive in case of an alarm
Time OK	40 [s]	Duration of time, during which the set temperature must have been at least reached after a temperature error, before the next error can be registered (here 40 seconds).
Time ERR	40 [s]	Duration of time, during which the temperature error must at least occur, before the error is indicated (here 40 seconds).
min.	3,0 [°C]	Permitted temperature deviation downwards, where no error is indicated (here 3,0°C)
max.	3,0 [°C]	Permitted temperature deviation upwards, where no error is indicated (here 3,0°C)
Door contact	ON	The door contact alarm is active
	OFF	The door contact alarm is inactive
Signal	ON	Acoustic signal is active if the alarm occurs
	OFF	Acoustic signal is inactive if the alarm occurs
Time OK	5 [s]	Duration of time, during which the plexiglass cover must be at least closed, before the error is registered again (here 5 seconds)
Time ERR	2 [s]	Duration of time, during which the plexiglass cover must be at least opened, before the error is indicated (here 2 seconds)
Speed	ON	The speed alarm is active
	OFF	The speed alarm is inactive
Signal	ON	Acoustic signal active at alarm message
	OFF	Acoustic signal inactive at alarm message
Time OK	20 [s]	Duration of time, during which the set speed after a speed error must have been at least reached again, before the next error can be registered (here 20 seconds)
Time ERR	20 [s]	Duration of time, during which the speed error must at least occur, before the error is indicated (here 20 seconds)
min.	10 [U/min]	Permitted speed deviation downwards, where no error is indicated (here 10 U/min.)
max.	10 [U/min]	Permitted speed deviation upwards, where no error is indicated (here 10 U/min.)
Power failure	ON	The power failure alarm is active
	OFF	The power failure alarm is inactive
Signal	ON	Acoustic signal is active if the alarm occurs
	OFF	Acoustic signal is inactive if the alarm occurs

- In the case of an alarm there is an optical display via alarm-LED or an acoustic message via buzzer (if active). At the same time the alarm is taken over in the diagnosis and can be read-out after the program end in the point SETUP / DIAG.
- A message to outside can be sent over a contact without potential via the SUB-D socket (see point 3.4.3).
- In pressing the [ENTER] key the optical and acoustic messages are acknowledged as well as the contact without potential is set back.
- In opening the cap the shaking drive is stopped and the heating is switched off, independent of the fact, if the door contact alarm is ON or OFF. If the cap is closed again, the shaking drive and the heating are operating automatically again.
- In the case of voltage loss during the program processing or in the operation mode „RUN“ the time which has already been processed is safed. If there is voltage again, the device starts automatically. The remaining time of the program is processed. This happens independent of the fact if the voltage loss alarm is active or not.
- The alarm definition for the voltage loss serves only for the activation of the optical/acoustic alarm message or the entry in the diagnosis after the voltage has returned.

F4

TIME:

Set-up of time and date

- In pressing the key [F4] you can change to the display field “set-up of time / date”.


Fig. 33: Setup-Menü



Fig. 34: Zeit- Datumseinstellung

The following parameters for the set-up of “time and date” (fig. 34) can be adjusted:

F1

DAY:

Set-up of the week day

- After pressing of [F1] the input field “DAY“ flashes.
- In pressing the keys [+/-] the desired week day can be set up.
- Confirm input with [ENTER].

F2

MONTH:

Set-up of the month

- After pressing of [F2] the input field “MONTH“ flashes.
- In pressing the keys [+/-] the desired month can be set up.
- Confirm input with [ENTER].

F3

YEAR:

Set-up of the year

- After pressing of [F3] the input field “YEAR“ flashes.
- In pressing the keys [+/-] the desired year can be set up.
- Confirm input with [ENTER].

F4

HOUR:

Set-up of hours

- After pressing of [F4] the input field "HOUR" flashes.
- In pressing the keys [+/-] the desired hour can be set up.
- Confirm input with [ENTER].

F5

MINUTE:

Set- up of the minutes

- After pressing of [F5] the input field "MINUTE" flashes.
- In pressing the keys [+/-] the desired minute can be set up.
- Confirm input with [ENTER].

MODE

MODE:

Display field

- The respective "active" menu is indicated, here "SETUP1"

3.3.6 Set-up menu "SETUP2"

- In pressing the [F5] key you change to the "SETUP2-menu" for the input of the device parameters.
- All set-ups and modifications are memorized internally and remain active also after voltage loss or switching off and on!


Fig. 35: Setup1 menu


Fig. 36: Setup2 menu

The following "device parameters" (fig. 36) can be adjusted:

F1

DEVICE:

Activation of different device options

Device options can be activated or deactivated in a list as well as device information can be called.

- An arrow marks the respectively activated parameter.
- In pressing the keys [+/-] the arrow can be moved.
- In pressing the [ENTER] key the YES / NO-definitions are directly modified or numeric values are released for modification (value flashes).
- In pressing the keys [+/-] the respectively flashing value can be adjusted.
- In pressing the key [ENTER] the adjusted value is confirmed.
- In pressing the [ESC] key you jump back to the SETUP2 menu.

Temp.Offset	0,0	Actual value correction for the 1.PT100 (for temperature regulation). Calibrating values can be adjusted. The digit sign can be modified by pressing the [F1 ... F5] keys.
Light	YES	Illumination active
	NO	Illumination inactive
August 2002	111	Software version (here version 1.11)
Ser. No.	39	Hardware serial number (here no. 39)

- In pressing the [F2 und F3] keys you change to menus where you input internal parameters.


Fig. 37: Setup2 menu


Abb. 38: Password input menu

- These menus are reserved for service purposes and are protected by passwords.
- Modifications of these parameters can only be carried out by authorized personnel.
- Improper modifications may destroy the unit.

F4	RESET
-----------	--------------

- In pressing the [F4] key a little bit longer the device is reset.

- All set-ups, programs, passwords etc. will definitely get lost!

3.4 Connection of External Devices

On the backside of the device the CERTOMAT® IS has two Sub-D sockets for analog signal exits with 9 poles, RS 232 interface as well as a collective alarm contact.


Fig. 39: Rear of the shaker
 1. Supply entry socket
 2. fuse T8A for the 230 V - version
 3. RS232 interface
 4. 9-poles socket „Analog Out“ for recorder connection, etc.,
 Pin 1 = Signal + /Temperature
 Pin 2 = Gnd / Temperature
 Pin 6 = Signal + / rpm
 Pin 7 = Gnd / rpm
 Pin 4 = collective alarm
 Pin 9 = collective alarm

3.4.1 Analog Exits (ANALOG OUT)

- For the external registration of speeds and temperatures, e. g. with a recorder, the CERTOMAT® IS has analog signal exits (0 ... 10 V).

The conditions are as follows:

- Speed: 0 V corresponds with 0 rpm.
10 V corresponds with 400 rpm.
- Temperature: 0 V corresponds with 0°C
10 V corresponds with 100 °C

☞ The moving load at the voltage entries must not be lower than 10 kΩ.

- Current signals (0 ... 20 mA) or (4 ... 20 mA) can optionally be made for disposal.

☞ The moving load at the current entries must not be higher than 500 kΩ.

3.4.2 RS 232 Interface

- The Sub-D socket with 9 poles "RS 232" is currently a Sartorius AG service interface.

3.4.3 Collective Alarm (optionally)

- A contact without potential is made available on PIN 4/9 via a Sub-D socket "ANALOG OUT" with 9 poles. By help of this contact, alarm situations can be reported externally.
- Contact loading capacity: 230 V AC, 0.5 A (Ohm load).

3.5 Maintenance and Cleaning

Except for necessary cleaning and the replacement of a defective drive belt or ball bearing, the CERTOMAT® IS is maintenance-free. The motor is a service-free brushless motor of external rotor type. All bearings are lubricated for their lifetime.

Disconnect the mains cable from laboratory grounded socket before cleaning and/or maintenance and service. Unintended starting of the shaking drive during cleaning and maintenance can cause personnel hazards.

3.5.1 Cleaning

1. The interior surface and the components built-in consist of stainless steel or powder coated steel sheet. For cleaning, we recommend usual household cleaning agents or alcohol. Do not use aggressive cleaning agents, e. g. cleaning agents with chlorine.
2. After incubation of living microorganisms and cells the CERTOMAT® IS may need to be disinfected. For disinfection we recommend disinfectants based on alcohols, such as meliseptol, see ordering information in the supplement. After cleaning and disinfection with such agents you should thoroughly dry and ventilate the incubation room.

You should remove broken glass, spilled liquids and other impurities immediately after they occur. Please take care that impurities do not dry.

3. If the shaking incubator is operating in a very damp or dusty environment, the drive belt should be checked regularly (about twice a year) for deposits or for wear and tear.
4. Take care that no foreign bodies can enter the air-circulation system. Remove dust from the outer covering grids of the built-in ventilator regularly.

3.5.2 Malfunction Correction

In case the device should not start or if it runs irregularly:

1. Check at the mains supply whether the cable is properly connected and laboratory mains supply is working. If other devices in the laboratory are malfunctioning at the same time, a distorted mains supply may be the cause. Check and stop such a kind of interference.
2. Check the fuses on the rear of the shaker, see position (3) and position (4), figure 37. A melted or burned-off fuse link can be detected visually.

Please note that the regular disinfection of the plexiglass canopy with alcoholic solutions may lead to the brittleness of the canopy. Such a treatment should only be carried out if it is necessary for the operation of the unit.

For replacement note the correct voltage/current of the fuses on the corresponding table.

3. If the malfunctioning continues please contact your Sartorius AG service representative.

At uneven and/or noisy operation or if the operation stops unintentionally, the drive belt can be defective. To check the condition of the drive belt, you have to remove the shaker table as well as the spillage tray:

1. Unscrew the 4 screws of the tray support located in the centre and remove the table carefully.

The tray support is heavy. Handle with care.

2. You can see the built-in shaker. If the drive belt is o.k., reassemble the shaking table and fasten the screws carefully.
3. If the drive belt is defective or shows wear it must be replaced. Replacement of the drive belt, other service measures and repairs at components inside the device, especially at the motor and the electronic control, must be done by authorized service personnel only.
4. If such a defect occurs, please contact your Sartorius AG service representative in charge or call directly:

Sartorius AG
Service-Zentrum Nord
Weender Landstraße 94-108
37075 Goettingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Telefax +49.551.308.3730

3.5.2.1 Maintenance of the Illumination Unit

Maintenance of the illumination unit is restricted to checking and replacing the fluorescent tubes. If other malfunctions occur, please contact your service representative of Sartorius AG.

3.5.2.2 Malfunctions of the Shaking Drive

- If the shaking drive runs noisy or irregularly the belt may be dirty. If the tray support does not turn although the shaking drive is working (you can hear the motor run), the drive belt may be damaged. If the motor is inoperative although the power supply and all necessary set-ups are available, the motor may be damaged.
1. For checking the belt drive the tray and tray support must be disassembled. Unscrew the 4 mounting screws in the centre of the tray support. Afterwards you remove the stainless steel spillage tray carefully. Then you can see the built-in shaker.
 2. At the assembly of the shaking table, please fasten the screws carefully crosswise.
- If the drive belt is damaged (torn) or the motor is suspected to be defective, please contact your service representative of Sartorius AG.

3.5.2.3 Malfunctions of the Optional Refrigerator

- The compressor of the refrigerating unit has an excess temperature protection, which opens if a certain ambient temperature is exceeded. It can take up to one hour until the excess temperature protection closes again and the device continues to operate.
- If the excess temperature protection cuts-off the refrigerator often or repeatedly, you should contact your service representative of Sartorius AG.

3.6 Warranty Regulations and Service

- If not agreed otherwise, Sartorius Stedim Biotech GmbH offers a warranty according to the General Terms and Conditions of Business, valid at the time of delivery, which you can obtain together with the order confirmation or the device or on request.
- The warranty applies to malfunction of the equipment, defective material and workmanship and resulting malfunctions. It includes the repair or the replacement of damaged parts. Not included in the warranty are materials of consumption and parts, which have become defective due to normal wear.
- If you use the device in a corrosion environment or if you use aggressive media you have to check the lastingness of the device. Malfunctions due to improper handling as well as defects and malfunctions due to aggressive media are not subject to the warranty.
- Claims for warranty have to be proved by a respective receipt of delivery. The warranty will lapse if any modifications or unauthorized repairs are performed on the unit by the purchaser or unauthorized third parties.
- You may return defective equipment to the Sartorius AG. Here repairs are carried out according to the valid maintenance and repair conditions.

Any units to be sent to us for repair have to be in perfect hygienic condition, clean and well packed. If parts were polluted by media and media components, they have to be cleaned, decontaminated, disinfected or sterilised, according to the valid safety rules of this particular application, e.g. for chemical or biological safety.

The sender has to prove that all safety guidelines have been observed. For this purpose you have to fill in the enclosed decontamination declaration and attach it to the apparatus. No apparatus will be repaired without declaration of decontamination or description of the measures taken.

Transport damages or subsequent cleaning or disinfection, if required, will be borne by the sender.

4 Technical Data, Ordering Information

4.1 Technical Data

4.1.1 Mechanical Design

Cabinet, size	W x H x D = 540 x 564 x 680 mm
Weight (without tray)	75 kg, heating version
Material	sheet steel case with scratch resistant varnish
Trays, type/size	type E/EU (420 x 420 mm) ¹⁾
max. load	depending on tray and mounting system applied

4.1.2 Mains Supply

Mains supply	230 V / 50 Hz via shock-proof plug
Power consumption	1200 W included cooling
Illumination	2 x 8 W
Fuses	fuse links 2 x T8A/250 V (230 V-version)
Radio interfer. suppress	according to EN 55011 and EN 55014

4.1.3 Operating Data

Orbital motion	Ø 25 mm or 50 mm, depending on the delivered device version; can be changed to alternate amplitude in-place at customers site by authorized service representatives of Sartorius AG	
Shaking speed	40 ... 400 rpm, according to the kind of loading	
Rotation speed accuracy	max. deviation ±1% from max value	
Incubation temperature	heating version	RT +8°C ... +60°C
	heating and cooling version	RT -10°C .. +60°C
Ambient temperature	+10 ... +60°C (+10 ... +30°C at heating and cooling version)	
Ambient humidity	10 ... 90 %, not condensing	

4.1.4 Analog Output

Socket, 9-poles Sub-D	0 ... 10 V (min. moving load = 10 kΩ) (Pin configuration see point 3.4, fig. 39) 0 ... 20 mA or 4 ... 20 mA (optionally) (max. moving load = 500 Ω)
Exactness	+/-0,3 V

4.1.5 Collective alarm

Collective alarm (optionally)	contact without potential (closer), max. 230 VAC (0,5 A Ohm load) via SUB-D-socket „analog out“ Pin 4/9
-------------------------------	--

(We reserve the right to carry out dimensional and construction modifications).

¹ Trays of modified design for new locking; rebuilding of existing trays for this locking on request. Please contact the Sartorius AG customer service.

4.2 Ordering Information

4.2.1 Configuration of Shakers

Cat.-no.	Features and Specifications
	Incubation shaker CERTOMAT® IS / 25 mm Version 230 V / 50 Hz
8864829	Heating version (UH)
8864845	Heating and cooling version (UHK)
	Incubation shaker CERTOMAT® IS / 50 mm Version 230 V / 50 Hz
8864926	Heating version (UH)
8864942	Heating and cooling version (UHK)

4.2.2 Accessories

Cat.-no.	Features and Specifications
	Finally equipped trays (with clamps for Erlenmeyer flasks)
	Tray E 420 x 420 mm
8853533	Tray incl. 39 clamps for Erlenmeyer flasks 100 ml
8853568	Tray incl. 20 clamps for Erlenmeyer flasks 250 ml
8853584	Tray incl. 14 clamps for Erlenmeyer flasks 500 ml
8853606	Tray incl. 9 clamps for Erlenmeyer flasks 1000 ml
	Universal - Trays
8853002	type EU (420 x 420 mm)
on request	Rebuilding of trays type E/EU of existing CERTOMAT® -shakers to new locking of the CERTOMAT® S-II

4.3 Safety Instructions Laboratory Shakers

(Safety Instructions Laboratory Shakers see appendix)

4.4 EG Declaration of Conformity

(EG Declaration of Conformity see appendix)

4.5 Declaration of Decontamination

(Declaration of Decontamination, see Appendix)

4.6 Information and Instructions on Disposal and Repairs

(Information and Instructions on Disposal and Repairs see appendix)

4.7 Dimension sheet CERTOMAT® IS

(Dimension sheet see appendix)

4.8 Prospectus CERTOMAT® IS

(Prospectus CERTOMAT® IS see appendix)

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Goettingen, Germany

Telephone +49.551.308.0
Telefax +49.551.308.3289
www.sartorius.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission
of Sartorius Stedim Biotech GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications, and
design of the equipment
without notice.

Status:
February 2008,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany

Printed in Germany on paper that has
been bleached without any use of chlorine
W4A000 · KT
Publication No.: SB6029-t08024
Order No.: 85030-520-53

Safety Instructions

Laboratory Shakers and Incubation Shaking Cabinets

Content

1	Introduction	3
2	Transport and Installation	3
2.1	Transporting Equipment	3
2.2	Requirements at the site of installation	3
3	Operating Instructions	4
3.1	Working with Media	4
3.2	Initial Startup and Normal Operation	4
4	Cleaning, Maintenance, and Service	5
4.1	Cleaning	5
4.2	Maintenance and Service	5
4.3	Returning Defective Units	5
4.4	Relocating equipment	5

1 Introduction

These Safety Instructions contain our recommendations for precautions in handling CERTOMAT® laboratory shakers and incubation shaking cabinets from the product program of the company Sartorius Stedim Biotech GmbH, especially for

- transport and installation
- equipping and initial startup
- operation
- cleaning, maintenance, and service

These safety instructions are non-transferable and meant exclusively for the equipment delivered. They are not valid for other equipment. Additional, specific safety regulations (i.e., legally or otherwise mandatory rules) are not dealt with by these instructions. Please consult your local authorities for further information.

Possible risks and dangers are marked with this symbol and highlighted like this paragraph. Ignoring these warnings may cause damage to the instrument or lead to other material or personal damage.

Steps to be carried out with special care or special aspects or issues to be considered are marked with safety instructions formatted as this paragraph.

[→ ..] Identifies references to contents of these instructions or other documents. Labels of illustrations, sections, or documents are in parentheses.

2 Transport and Installation

2.1 Transporting Equipment

The shakers and especially the incubation shaking cabinets are heavy. Use only suitable transportation and lifting tools when moving or setting up units.

Make sure clearances are sufficient along the transport route and that the combined weight of lifting and transportation tools does not exceed the maximum load permitted for the floor or ground.

Use only qualified personnel to transport the equipment. Make sure personnel is protected from injury during transport.

Do not remove any transport protection before the equipment has been transported to the site of installation.

2.2 Requirements at the site of installation

1. The floor or the laboratory bench where the equipment is to be installed must be able to support the weight of the fully equipped units and systems.
2. The installation place should be sufficiently dimensioned and feature a non-slip surface.
3. The equipment has to be easily accessible for operation, maintenance, and service tasks.
4. The shakers and the incubation shaking cabinets may transmit vibrations to the installation surface especially when carrying a heavy load. For this reason, install the shakers so that surrounding equipment is not affected.
5. The power supply must match the rating indicated on the type plate and must be grounded correctly.

Use a level to align the equipment after setting it up!

Before stacking equipment (max. 3 units), make sure that the bottom most device is properly aligned horizontally!

3 Operating Instructions

1. Make sure, that only authorized people have access to the working area.
2. Personnel has to be instructed about applicable safety instructions and other legally or otherwise compulsory regulations.

3.1 Working with Media

The utilized media and substances can be hazardous and may pose dangers that are specific to the substances and the processes and thus cannot be described here in detail. Your company should issue corresponding safety regulations and carefully instruct the operating personnel.

Recommended basic protection equipment includes such items as suitable working clothes, gloves, goggles, and respiratory equipment, if necessary.

3.2 Initial Startup and Normal Operation

1. Use only units, equipment, accessories, and spare parts released by the company Sartorius Stedim Biotech GmbH for use with the shaker.
2. Check all components for damage, especially glass components such as jars. Any damaged components must not be used.
3. Never use the shaker without the tray. Switch the unit off before installing/uninstalling trays, mounting sets, or additional accessories.
4. When mounting the trays, leave enough space between the unit and other objects, especially in case of trays extending over or beyond the shaker table or bench. Make sure that trays are properly installed and fastened.
5. Fasten the tray mountings and accessories properly. Distribute loads evenly on the tray when loading the shaker with the shaking vessels (Erlenmeyer flasks, bottles, etc.).
6. In case of incubation shaking cabinets, make sure the tray is firmly positioned; close the doors before starting the shaking process.
7. Before starting the shaker and during normal operation, check that the vessels are firmly positioned in their holders. The vessels must not become loose, rattle, or touch one another or the walls of the shaker. Switch the unit off immediately if vessels become loose in their holders.
8. Operate the equipment only within the permitted range of operating parameters. Comply with the equipment instructions and the technical data.

Head the lateral movement of the trays and the resulting impact or crushing hazards, especially with large trays, in case of low clearances to other equipment or objects or large amplitudes and heavy loads.

Never put fingers between tray support and frame when the unit is active. The space between these components is narrow due to constructional constraints. Risk of injury!

9. After switching off the shaker, wait until the shaker support has come to a complete standstill before removing or adding vessels.

4 Cleaning, Maintenance, and Service

Regular cleaning and maintenance ensures the equipment's functionality and safety.

 Comply with the legal safety regulations applicable to the field of application concerning the handling of microbial material and contaminated equipment.

 Prior to any cleaning or maintenance work, disinfect or sterilize any biologically contaminated components according to the safety regulations.

 Prior to any cleaning or maintenance work, switch the unit off and pull the power cable to disconnect the unit from the mains power supply. Make sure that the unit cannot be restarted or activated accidentally.

Servicing and repairing drive systems, replacing drive belts, or working on electrical/electronic equipment should only be performed by qualified service personnel.

4.1 Cleaning

1. Clean the surfaces of the units with standard laboratory cleaners. Do not use corrosive agents (e.g., chlorous agents).
2. Broken glass or spilled liquids should be removed as quickly as possible.
3. For incubator shakers special hints are given in the operating manual.

4.2 Maintenance and Service

Maintenance and troubleshooting by the user should be restricted to the following tasks:

- Checking the voltage supply in the lab and the equipment connections.
- Checking and replacing defective fuses.
- For incubation shaking cabinets: Replacing defective light fittings.

If problems or malfunctions cannot be solved or remedied, please contact your Sartorius AG representative or contact us directly:

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Goettingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Fax +49.551.308.3730

A written cleaning and decontamination certificate has to be provided when requesting service technicians of the company Sartorius AG to perform maintenance tasks.

4.3 Returning Defective Units

Defective units or components can be returned to the responsible Sartorius AG representative or directly to Sartorius AG.

1. The equipment has to be clean, hygienically flawless, and carefully packaged.
2. Contaminated parts have to be cleaned, disinfected, or sterilized according to the respective and applicable safety guidelines. Compliance with the guidelines has to be documented with a cleaning and decontamination certificate to be enclosed with the shipment.
3. Negotiate the return shipment with the Service Dept. before shipping any equipment.

4.4 Relocating equipment

1. Comply with the described cleaning and maintenance measures described to avoid biological or chemical hazards.
2. Disconnect the equipment from the mains power supply. Remove vessels and additional equipment or components inside of the unit or connected to the unit.
3. If applicable, reinsert the transport protection. Read and comply with the additional notes under "Transporting Equipment" (2.1).

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Goettingen, Germany

Telephone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission of
Sartorius Stedim Biotech GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications, and
design of the equipment
without notice.



Status:
February 2008,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany

Printed in Germany on paper that has
been bleached without the use of chlorine
W4A000 · KT
Publication No.: SB6043az08022
Order No.: 85030-520-67

Declaration of Conformity
according to EU Machinery Directive 89/392/EU,
Appendix II A

Company	Sartorius Stedim Biotech GmbH	
Address	August-Spindler-Strasse 11, 37079 Goettingen; Germany Phone +49.551.308.0, Fax +49.551.308.32 89 www.sartorius-stedim.com	
	<p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized.</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p>	
Designation of the device	CERTOMAT® IS	
Model, version	Incubation Shaking Cabinet	
Cat.-No.	8864829, 8864837, 8864926, 8864934, 8864845, 8864853, 8864942, 8864953	
Relevant directives of the EC	EU Machinery Directive (89/392/EU) in version 93/44/EU EU Directive on Electromagnetic Compatibility (89/336/EU) in version 92/31/EU EU Low-Voltage Equipment Directive (73/23/EU)	
Applied harmonized standards	EN 292-1, EN 292-2 EN 50081-1, EN 50082-1, EN 55011 EN 61010-1	
Applied national standards and technical specifications		
Date and Signature		
Function of the Signatory	V. Niebel Managing Director	 Dr. Susanne Gerighausen Quality Management Representative

Declaration about decontamination and cleaning of equipment and components

When returning equipment or components, please describe on page 2 of this form the problem(s) or fault(s) you have found. Please also indicate the remedial actions you require.

To protect our personnel, we require all equipment or components be free of biological, chemical, or radioisotopic contaminants. We will only accept such equipment or components when:

- the equipment or components have been adequately cleaned and decontaminated.
- this declaring document has been completed, signed and returned by an authorized person.

Please help us in assuring a safe, hazard-free work environment.

A. Description of the Equipment or Component(s)

Description / Cat. No.	
Serial no.	
No. of invoice/delivery note	
Date of delivery	

B. Contamination / Cleaning

Attention: Please specify exactly the biological, chemical, or radioisotopic contaminant	Attention: Please describe the cleaning and decontamination procedure/method
The equipment was contaminated with	and it has been cleaned and decontaminated by

F-13-01 Version 01	Page 2 / 4 Effective Date 01.01.2008	<h2>Return of material</h2>	
------------------------------	---	-----------------------------	--

C. Legally binding declaration

I (we) certify that all information given in this form is correct and complete.
 The equipment and components have been adequately decontaminated and cleaned according to the legal requirements.
 No chemical or biological or radioisotopic risks remain that can endanger exposed persons' safety or health.

Company / Institute	
Address / Country	
Tel. / Fax (with area code)	
Name of the authorized person	
Position	
Signature / Date	

D. Reason for return

wrong delivery
 exchange
 repair
 modification
 disposal
 other

E. Please describe the problem(s) or fault(s) you have found (for repair) and/or indicate the remedial actions you require

F-13-01	Page 3 / 4	Return of material	
Version 01	Effective Date 01.01.2008		

F. Reserved for Sartorius-Service Center

Notes

Please pack the unit properly and send it freight paid to your local service supplier or directly to Sartorius AG Servicezentrum Nord, Germany.

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Goettingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Fax +49.551.308.3730

F-13-01 Version 01	Page 4 / 4 Effective Date 01.01.2008	Return of material	 sartorius stedim biotech
--------------------------	---	---------------------------	---

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephone +49.551.308.0
Telefax +49.551.308.3289
www.sartorius-stedim.com

Copyright by Sartorius Stedim
Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission
of Sartorius Stedim Biotech
GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below. Sartorius Stedim
Biotech GmbH reserves the right
to make changes to the technology,
features, specifications, and
design of the equipment
without notice.


Status:
February 2008, Sartorius Stedim
Biotech GmbH,
Goettingen, Germany

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.


For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:

Sartorius AG
Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.


Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet/manual or visit our Internet website (www.sartorius.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

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
CERTOMAT[®] IS
For controlled benchtop
incubation

CERTOMAT® IS For controlled bench top incubation

The CERTOMAT® IS from Sartorius has proven its value in many laboratories world-wide. In microbiology, this benchtop shaking incubator allows bacteria or yeast cells to grow in a temperature-controlled environment from 8°C above ambient temperature up to +60°C.

It is thus ideally suited for the temperature range most commonly used in biology laboratories working with microorganisms. With an optional integrated cooling, the temperature range can be lowered to 10°C below ambient temperature.

Due to its compact design, the CERTOMAT® IS can also be used by molecular biologists having to expand their cultures beyond the petri dish size. As an additional benefit this unit can be programmed to operate for specified periods using defined ranges of speed and temperature.

For example, a pre-culturing stage and 4 successive incubation stages can easily be programmed.

Incubation of cultures on the bench has never been easier than with the CERTOMAT® IS.

Its small footprint and superior design make it the ideal incubator shaker for every laboratory.

For applications requiring larger incubation capacities you can expand using our CERTOMAT® BS-1 and CERTOMAT® BS-T incubator shakers.

Benefits

Compact design

Fully programmable with up to 5 programmes, password-protected

Menu language selectable by the user

Integrated spill tray

Contact switch stops the shaker when lid is opened

Integrated active cooling system (UHK versions only)

Alarm points may be set by the user

RS 323 Interface


Accepts 3 liter flasks

Accessories

Reference	Description
	Tray type E (420×420 mm) equipped with stainless steel clamps for Erlenmeyer flasks
885 3533	39 clamps for flasks 100 ml
885 3568	20 clamps for flasks 250 ml
885 3584	14 clamps for flasks 500 ml
885 3606	9 clamps for flasks 1,000 ml
	Tray type E (420×420 mm) equipped with plastic clamps for Erlenmeyer flasks
885 3688	39 clamps for flasks 100 ml
885 3666	19 clamps for flasks 250 ml
885 3677	14 clamps for flasks 500 ml
	Universal tray to be completed with clamps, racks or mounting system Type EU (420×420 mm)
885 3002	
	Stainless steel clamps for tray type EU
885 4505	for Erlenmeyer flasks 25 ml (max. 49)
885 4513	for Erlenmeyer flasks 50 ml (max. 48)
885 4521	for Erlenmeyer flasks 100 ml (max. 24)
885 4556	for Erlenmeyer flasks 250 ml (max. 17)
885 4572	for Erlenmeyer flasks 500 ml (max. 12)
885 4599	for Erlenmeyer flasks 1,000 ml (max. 8)
885 4610	for Erlenmeyer flasks 2,000 ml (max. 4)
885 4629	for Erlenmeyer flasks 3,000 ml (max. 4)
885 4564	for Fernbach flasks 450 ml (max. 8)
885 4600	for Fernbach flasks 1,800 ml (max. 4)
885 4640	for Fernbach flasks 2,800 ml (max. 4)
	Plastic clamps reinforced with glass fibre
885 4700	for Erlenmeyer flasks 100 ml (max. 24)
885 4711	for Erlenmeyer flasks 250 ml (max. 17)
885 4722	for Erlenmeyer flasks 500 ml (max. 12)
885 4733	for Erlenmeyer flasks 1,000 ml (max. 8)
	Hinged racks for test tubes (4 racks max. on tray EU)
885 3134	for 64 tubes Ø 14 mm
885 3142	for 42 tubes Ø 16 mm
885 3150	for 36 tubes Ø 18 mm
885 3169	for 33 tubes Ø 20 mm
885 3185	for 18 tubes Ø 25 mm
885 3177	for 16 tubes Ø 30 mm
	Hinged racks for centrifugation tubes (4 racks max. on tray EU)
885 3088	for 42 tubes Ø 16 mm
885 3096	for 36 tubes Ø 18 mm
885 3193	for 33 tubes Ø 20 mm
885 3240	for 16 tubes Ø 30 mm

**Reference****Description**

885 0321


Holders for microtiter plates, stainless steel
for 1 standard 96-well plate or deepwell plate
standard plates: max. 12 holders on EU tray
deepwell plates: max. 9 holders on EU tray

886 4497
886 0416
886 4470

Sticky tape for universal trays
Standard, roll of 50 m, 30×1 mm
Premium, roll of 10 m, 30×1 mm, repeated use
Anti-skid layer, 380×450 mm, for individual cut

885 4238
885 4254

Universal mounting system
Basic element B-2 for tray EU
Universal clamping rod for basic elements B-2

885 4262
885 4270

**Mounting set for separation funnels
in combination with universal mounting system (not shown)**
for 5 separating funnels 50 or 100 ml
for 3 separating funnels 250, 500 or 1,000 ml

886 1005
886 1013
886 1021
886 1022

**Shaking flasks, DURAN, Erlenmeyer type,
3 baffles at 120°, straight rim**
Erlenmeyer flasks 300 ml, pack of 10
Erlenmeyer flasks 500 ml, pack of 10
Erlenmeyer flasks 1,000 ml, pack of 10
Erlenmeyer flasks 2,000 ml, pack of 10

886 1064
886 1072
886 1080

**Shaking flasks, DURAN, Erlenmeyer type,
3 baffles at 120°, straight rim, connector GL 14**
Erlenmeyer flasks 300 ml, pack of 10
Erlenmeyer flasks 500 ml, pack of 10
Erlenmeyer flasks 1,000 ml, pack of 10

886 0998

**Shaking flasks, DURAN, Erlenmeyer type,
3 baffles at 120°, narrow neck for plug**
Erlenmeyer flasks 500 ml, pack of 10

Technical Data

Ordering information

CERTOMAT® IS version with circulation/heating (UH)		
230 V/50 Hz	886 4829	CERTOMAT® IS/25 mm
230 V/50 Hz	886 4926	CERTOMAT® IS/50 mm
115 V/60 Hz	886 4837	CERTOMAT® IS/25 mm
115 V/60 Hz	886 4934	CERTOMAT® IS/50 mm
CERTOMAT® IS version with circulation/heating/cooling (UHK)		
230 V/50 Hz	886 4845	CERTOMAT® IS/25 mm
230 V/50 Hz	886 4942	CERTOMAT® IS/50 mm
115 V/60 Hz	886 4853	CERTOMAT® IS/25 mm
115 V/60 Hz	886 4953	CERTOMAT® IS/50 mm

Technical data

Mechanical	Shaking speed: Amplitude: Accuracy: Motor:	40 to 400 rpm 25 mm and 50 mm, convertible ±1% of max. value brushless motor, rotor external
Temperature	Temperature range: Version with cooling: Air circulation:	ambient +8 °C to +60 °C ambient -10 °C to +60 °C 1 fan, 180 m ³ /h
Trays	Suitable trays: Adjustment of speed and temperature: Display:	type E and EU (420×420 mm) foil key pad LCD alphanumeric, 4 lines, 20 characters each
Safety	Memory: Alarm: Protection:	automatic restart after power failure acoustic and on LCD IP 21
Operating conditions	Temperature: Humidity:	from +10 °C to +45 °C from 10% to 90%, not condensing
Documentation	Analog output:	for speed and temperature
Electrical	Heating capacity: Cooling capacity:	650 W 300 W, cooled versions only
Dimensions and weight	External (L × H × D): Internal (L × H × D): Weight:	540×565×680 mm 505×370×510 mm 62 kg

All instruments are delivered without trays and other accessories.

For incubating cells or mixing liquids a tray is needed together with additional accessories to hold shaking flasks, separation funnels or tubes.

Sartorius AG
Weender Landstrasse 94-108
37075 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius.com

Sartorius BBI Systems GmbH
Schwarzenberger Weg 73-79
34212 Melsungen, Germany

Phone +49.5661.71.3400
Fax +49.5661.71.3702

www.sartorius-bbi-systems.com

Sartorius North America Inc.
131 Heartland Blvd.
Edgewood, New York 11717, USA

Phone +1.631.254.4249
Toll-free +1.800.3687178
Fax +1.631.254.4253

Sartorius BBI Systems, Inc.
2800 Baglyos Circle
Bethlehem, PA 18020, USA

Phone +1.610.866.4800
Fax +1.610.866.4890

Sartorius Ltd.
Longmead Business Park
Blenheim Road, Epsom
Surrey, KT19 9 QQ, U.K.

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Phone +81.3.3740.5407
Fax +81.3.3740.5406

Specifications subject to change
without notice.

Printed in Germany on paper that has been
bleached without any use of chlorine.

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The reproduction, distribution and utilization of this document as well as the communication of its contents to others without explicit authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.
Sartorius Stedim Singapore Pte. Ltd.
Sartorius Stedim Australia Pty. Ltd.