


Parallel GC for Fast Refinery Gas Analysis




Background information

- ◆ **Refinery gas is a collection of various gas streams produced in refinery processes**
- ◆ **Refinery gas composition is very complex, so determining individual components of each gas stream is a challenge**
- ◆ **Analysis of refinery gas is essential for process and product quality control**
- ◆ **Sample complexity requires increased separation power**
- ◆ **Labs must increase their efficiency or productivity**
- ◆ **The Agilent 7890A GC supports a third parallel detector (TCD) for simultaneous detection across three channels, which provides a complete analysis of refinery gas in 6 minutes**



Fast RGA Configuration




Column 1 HayeSep Q 80/100, 0.5m
 Column 2 HayeSep Q 80/100, 6ft
 Column 3 Molsieve 5A 60/80, 6ft
 Column 4 HayeSep Q 80/100, 3ft


Column 5 Molsieve 5A 60/80, 8ft
 Column 6 DB-1, 2m x 0.32mm x 5um
 Column 7 HP-AL/S, 25m x 0.32mm x 8um

Seven Columns in GC Oven

Five valves, seven columns, and still a lot of room for people to work within the oven.




Enhance Gas Analysis with Union Connector



... Peak shape enhancement

Column Union Connector



Fast analysis GC conditions

Oven program: 60 °C (1min), 20 °C/min, 100 °C, 30C/min 190C

Split temperature: 150 °C, Valve box: 120 °C

FID for Hydrocarbons

Columns:

(1) DB-1, 2 m x 0.32 mm x 5 um

(2) HP-AL/S 25 m x 0.32 mm x 8 um

Carrier Gas: He @ 3.3 mL/min, constant flow mode (12.7 psi at 60 °C)

Split Ratio: 1:60

Loop: 25 uL

Back TCD for Permanent Gases

Columns:

(1) HayeSep Q 80/100, 0.5m

(2) HayeSep Q, 80/100, 6ft

(3) Molsieve 5A, 60/80, 6ft

Carrier Gas: He @ 25 mL/min (36 psi at 60 °C), constant flow mode

Pre-column Flow: 22 mL/min at 60 °C (7 psi), constant pressure mode

Loop: 500 uL

2nd TCD for Hydrogen

Columns:

(1) HayeSep Q 80/100, 3ft


(2) Molsieve 5A, 60/80, 8ft

Carrier Gas: N₂ 24 mL/min, (26 psi at 60 °C), constant flow mode




Pre-column Flow: 24 mL/min at 60 °C (7 psi), constant pressure mode

Loop: 500 uL



Application and Result- RGA standards






Application and Result- NGA standards






Application and Result- Customer's sample




Flexibility for Hydrocarbon Analysis



H2S and COS Analysis




H₂ calibration (4-point)



Excellent Repeatability for Both RT and Area

Compounds	RT			Area		
	Average	Std. Dev.	n=6 RSD%	Average	Std. Dev.	n=6 RSD%
C6+	0.99648	0.00031	0.03	59.01	1.10	1.86
methane	1.50780	0.00046	0.03	490.02	1.45	0.30
ethane	1.70788	0.00052	0.03	807.40	2.35	0.29
ethylene	1.95732	0.00071	0.04	472.31	1.31	0.28
propane	2.41706	0.00075	0.03	1950.35	5.96	0.31
cyclopropane	3.18506	0.00075	0.02	145.62	0.45	0.31
propylene	3.26195	0.00072	0.02	732.90	2.01	0.27
i-butane	3.64883	0.00055	0.02	885.04	3.15	0.36
n-butane	3.79161	0.00070	0.02	682.13	2.59	0.38
propadiene	3.86098	0.00095	0.02	109.08	0.65	0.60
acetylene	3.96990	0.00120	0.03	348.17	2.39	0.69
t-2-butene	4.47301	0.00106	0.02	507.88	2.59	0.51
1-butene	4.57118	0.00110	0.02	332.39	2.03	0.61
i-butylene	4.67529	0.00121	0.03	260.95	1.95	0.75
c-2-butene	4.76367	0.00112	0.02	403.80	3.47	0.86
i-pentane	5.03923	0.00090	0.02	45.03	0.05	0.11
n-pentane	5.14583	0.00099	0.02	69.23	0.40	0.58
1,3-butadiene	5.25906	0.00122	0.02	485.49	3.66	0.75
propyne	5.36385	0.00155	0.03	101.08	0.41	0.40
t-2-pentene	5.58664	0.00121	0.02	82.85	0.66	0.79
2-methyl-2-butene	5.68220	0.00117	0.02	62.54	0.61	0.98
1-pentene	5.75553	0.00126	0.02	39.57	0.38	0.96
c-2-pentene	5.83970	0.00131	0.02	59.08	0.50	0.85
CO2	2.18561	0.00221	0.10	2040.33	2.37	0.12
O2	2.72634	0.00060	0.02	930.68	6.53	0.70
N2	3.25170	0.00044	0.01	22500.18	68.87	0.31
CO	4.61692	0.00083	0.02	903.09	2.77	0.31
H2	0.9869	0.00099	0.10	16097.38	106.53	0.66

Reporting



The Macro program will give a report in:

- ◆ Mole%, Weight% or Volume%, or any combination of the three
- ◆ Heat values for the gas analyzed and other standard calculations
- ◆ Calculations are made using ASTM, GPA or ISO formulas

A library contains all physical constants for the most common compounds found in refinery gas

Reporting

Reports can be calculated using formulas given in the ASTM/GPA or ISO standards

Calculation Procedure

Choose the calculation procedure you wish to use

ASTM/GPA

ISO

OK Cancel

ASTM/GPA Calculation Procedures at 60°F

Choose the Report Style you wish to use

Main Report Style

Mole%

Weight%

Volume%

Choose the Base Pressure for Calculation

Base Pressure

Choose the Base Pressure you wish to use

Units must be PSI absolute

14.69

Single Line Totals for a real gas

Gross Heat Value per CF

Gross Heat Value per Pound

Gross Heat Value per Gallon

Gallons per 1000 CF of gas

Mean Molecular Weight

Net Heat Value per CF

Net Heat Value per Pound

Net Heat Value per Gallon

Compressibility Factor

Specific Gravity

ISO Calculation Procedures

Choose the Temperatures for Calculation

Calculation temperature @ 101.325 kPa

Melting Temp

0°C 15°C 20°C

Combustion Temp

0°C 15°C 20°C 25°C

Choose the Report Style you wish to use

Main Report Style

Mole%

Weight%

Volume%

Single Line Totals for a real gas

Superior Heat Value (Molar)

Superior Heat Value (Mass)

Superior Heat Value (Volume)

Mean Molecular Weight

Relative Density

Wobbe Index (Superior)

Inferior Heat Value (Molar)

Inferior Heat Value (Mass)

Inferior Heat Value (Volume)

Compressibility Factor

Density

Wobbe Index (Inferior)

Reporting

- The reports use some of the features of the ChemStation report layout screen, which allows you to customize the reports by adding or removing text
- The main body of the report where the compounds are listed will always give the compound name, RT, area, and then Mole%, Weight% and Volume%, depending on what you have chosen.

Example of an ASTM/GPA report-Main report

Compound	RT	Area	Mole%	Wt%	Vol%
C6+	0.675	224.3	0.41	1.22	0.38
Hydrogen	0.986	16189.3	15.81	1.11	15.94
Methane	1.508	493.7	6.36	3.55	6.39
Ethane	1.710	812.6	5.39	5.65	5.39
Ethylene	1.962	474.8	3.18	3.11	3.18
Carbon Dioxide	2.186	2041.2	3.16	4.85	3.17
propane	2.422	1960.2	8.54	13.11	8.45
Oxygen	2.726	1006.5	2.30	2.56	2.31
cyclopropane	3.192	144.9	0.53	0.77	0.52
Nitrogen	3.249	22697.7	28.68	27.99	28.89
propylene	3.269	734.1	3.22	4.73	3.20
i-butane	3.655	887.1	2.87	5.82	2.81
n-butane	3.797	678.8	2.23	4.51	2.17
propadiene	3.868	111.1	0.97	1.35	0.96
methane	3.948	2811.5	6.35	3.55	6.39
acetylene	3.978	343.6	1.82	1.65	1.82
trans-2-butene	4.480	501.1	1.63	3.18	1.58
1-butene	4.578	328.9	1.05	2.06	1.03
Carbon Monoxide	4.616	911.8	1.59	1.56	1.61
2-methylpropene	4.683	258.0	0.85	1.67	0.83
cis-2-butene	4.771	395.3	1.29	2.52	1.26
i-pentane	5.045	44.5	0.10	0.26	0.10
n-pentane	5.151	57.6	0.13	0.32	0.12
1,3-butadiene	5.266	476.0	1.54	2.90	1.51

Real Gas Values Dry

Gross Heat Value per CF 1073.97 BTU

Net Heat Value per CF 989.17 BTU

Mean Molecular Weight 28.70

Specific Gravity 0.9949



Agilent New fast RGA Highlights



One 7890A GC configured with three parallel channels with simultaneous detection provides a comprehensive, fast, and high-resolution analysis of refinery gas with one injection in six minutes




Use of optimized columns allows faster hydrocarbon and permanent gas analysis with the same oven temperature program



A third TCD can be used for improving hydrogen detection and linearity by using nitrogen as a carrier



New, easy-to-use tubing connectors based on Capillary Flow Technology are used to connect valves and capillary columns to improve chromatographic performance including peak shape



Macro program provides automated gas properties calculation, extensive reporting options

Advantages of the Agilent RGA solution

▣ Consistent results / highest precision / highest efficiency

- The FID's wide dynamic response range improves the analysis of samples containing high- and low-concentration components
- The Agilent TCD uses a single filament flow switching design that eliminates the need for a reference column, improving stability, and the small dead volume improves sensitivity
- 0.001 psi EPC GC for consistent valve switch timing setup
- New union connector based on Capillary Flow Technology is used for the connection of capillary columns to valves, enhancing the peak shape
- Valve in separate heated zone to ensure no compromising or cold spots
- 3rd TCD can be operated simultaneously with the other two detectors, giving three parallel channels for gas analysis

▣ Proven solutions backed by world-class experts

- Our experts have invested a combined total of over 200 membership years in ASTM- the world's most trusted source for standards development
- Our 40 years of trusted industry partnership experience helping make labs around the world more productive