



850-DS
Dissolution Sampling Station

Operator's Manual

Notices

Manual Part Number

70-9074 Rev D
March 2018

Copyright

© Agilent Technologies, Inc. 2018

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Agilent Technologies, Inc.
3501 Stevens Creek Blvd.
Santa Clara, CA 95052 USA

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Content

1 Safety 7

Electrical Hazards 8

Warning 9

Caution 9

Note 9

Information Symbols 10

2 Introduction 11

Serial Number Format 13

3 Setting Up the 850-DS 15

Initial Setup 16

Equipment Ratings / Environmental Requirements for Installation 17
Specification 18

Equipment, Parts, and Accessories As Shipped 19

Parts and Accessories 19

Connections 20

Tubing 21

Power Cord 22

Apparatus Control Cable 22

Synchronized Start 23

Initial Power Up 24

Sample Tray Installation 24

Updating 850-DS Language Screens 25

Optional Filter Module	27
Changing Trays / Needle Replacement	28
4 Operating the 850-DS	29
Startup Screen	30
Main Screen	32
Menu Screens	36
Methods	36
Method Properties (Apparatus 1/2/5/6)	38
Sample Collection Parameters (Apparatus 1/2/5/6)	40
Method Properties (Apparatus 3/7)	43
Timepoint Number Selection	46
Timepoint Parameter Entry (Apparatus 1/2/5/6)	47
Timepoint Parameter Entry (Apparatus 3/7)	48
Run Method Screen 1	49
Run Method Screen 2	50
Run Method Screen 3	52
Run Method Screen 4	53
Run Method Screen 5	53
850-DS Method Transfer (to another 850-DS)	54
Alarms	56
Calibration	57
Cleaning	61
Diagnostics	64
Diagnostics Screen 2	66
Diagnostics Screen 3	67
Diagnostic Screen 4 (Apparatus 1/2/5/6)	69
Diagnostic Screen 4 (Apparatus 3/7)	71
Instrument Settings	72
Instrument Settings Screen 2	74
Manufacturing Configuration	78
Reports	79
User Access	80

Content

Dialog Screens (Blue)	83
Warning Screens (Yellow)	85
Error Screens (Red)	88

5 Maintenance and Troubleshooting 91

Preventive Maintenance 92

Clean System	92
Periodic Maintenance	92

Printer Maintenance 93

Removing a Thermal Paper Roll	93
Inserting a Thermal Paper Roll	95
Thermal Printer Test	98

Obtaining Warranty and Other Services 99

This page was intentionally left blank, except for this message.



1

Safety

Electrical Hazards 8

The Agilent 850-DS Dissolution Sampling Station has been designed and tested so that when used properly you have an accurate, fast, flexible, and safe instrument.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

Operation of a Agilent 850-DS involves the use of aqueous liquids and various pharmaceutical dosage forms. Unskilled, improper, or careless use of this instrument can create shock hazards, fire hazards, or other hazards which can cause death, serious injury to personnel, or severe damage to equipment and property.

Information on safety practices is provided with your instrument and operation manuals. Before using your instrument or accessories, you must thoroughly read these safety practices.

Observe all relevant safety practices at all times.

Electrical Hazards

The dissolution apparatus contains electrical circuits, devices, and components operating at dangerous voltages. Contact with these circuits, devices, and components can cause death, serious injury, or painful electric shock.

Panels or covers that are retained by fasteners which require the use of a tool for removal may be opened only by Agilent-trained, Agilent-qualified, or Agilent-authorized service engineers. Consult the manuals or product labels supplied with the dissolution apparatus to determine which parts are operator-accessible.

Application of the wrong supply voltage, connection of the instrument to an incorrectly wired supply outlet, or lack of proper electrical grounding can create a fire hazard or a potentially serious shock hazard and could seriously damage the instrument and any attached ancillary equipment.

Always use a three-wire outlet with ground connection which is adequately rated for the load. The installation must comply with local, state, and federal safety regulations.

Do not connect the instrument to the main power supply until you have made sure that the operating voltage is correctly set for the main power supply in the specific outlet in your laboratory to which the equipment will be connected.

Warning

WARNING

A 'Warning' message appears in the manual when failure to observe instructions or precautions could result in death or injury.

Read all warnings and cautions carefully and observe them at all times.

Caution

CAUTION

A 'Caution' message appears in the manual when failure to observe instructions could result in damage to equipment (Agilent supplied and / or other associated equipment).

Note

NOTE

A 'Note' appears in the manual to give advice or information.

Information Symbols



A depressed power button in the front of the unit (with green light) indicates main power on.



An extended power button in the front of the unit (with green light off) indicates main power off.



Indicates single-phase alternating current



Indicates the product complies with the requirements of one or more European Union (EU) directives.



Indicates specific equipment meets standards of safety. These products are safe for use in the workplace for North America.



Indicates that this product must not be disposed of as unsorted municipal waste.

All Agilent products that are subject to the WEEE directive shipped after August 13, 2005 are compliant with the WEEE marking requirements. Such products are marked with the "crossed out wheeled bin" WEEE symbol in accordance with European Standard EN 50419.

For more information on collection, reuse, and recycling systems, please contact your local/regional waste administration, your local distributor, or Agilent.



Indicates the product complies with regulatory compliance requirements of New Zealand and Australia.



2

Introduction

Serial Number Format 13

The 850-DS allows you to autonomously collect samples from a dissolution apparatus. Samples may be collected into various sample holders including 15 mL test tubes, 2 mL vials, and 96-well plates. The vials and well plates may be conveniently placed on a sample tray that is compatible with Agilent HPLC autosamplers. The 850-DS provides two basic modes of collecting samples from a dissolution apparatus: manual control and method control. Manual control provides full control of the instrument through the Diagnostics screens. You can independently control all of the features of the instrument.

CAUTION

In manual control as it is possible to pump volumes greater than those that the sample vessels may hold.

Method control allows you to enter all the parameters relevant to a sample collection; prime, sample, and purge volumes, sample collection times, and self cleaning. If using an Agilent dissolution apparatus, the 850-DS will fully control the dissolution apparatus through a serial interface.

WARNING

The dissolution apparatus contains electrical circuits, devices, and components operating at dangerous voltages. Contact with these circuits, devices and components can cause death, serious injury, or painful electric shock.

CAUTION

Panels or covers that are retained by fasteners which require the use of a tool for removal may be opened only by Agilent-trained, Agilent-qualified, or Agilent-authorized service engineers.

Serial Number Format

The serial number contains 10 characters and follows this syntax:

CC1234xxxx

Syntax Code	Meaning	Description
CC	Country of origin	2 alpha characters matching the required trade designation for the country of origin
12	Year of manufacture	'16' for 2016, '17' for 2017, etc.
34	Week of manufacture	'01' for week 1 to '52' for week 52

This page was intentionally left blank, except for this message.



3

Setting Up the 850-DS

Initial Setup	16
Equipment, Parts, and Accessories As Shipped	19
Connections	20
Optional Filter Module	27
Changing Trays / Needle Replacement	28

Initial Setup

Complete the following sections to initially set up the Agilent 850-DS and all other system components.

- 1 Prepare an area next to a dissolution apparatus where the Agilent 850-DS will be placed. The preferred location is on the right side of the dissolution apparatus.

Height	40 cm (60 cm with Filter Module)
Width	39.0 cm (48.0 cm with filter waste bin)
Depth	60.0 cm

- 2 Remove the Agilent 850-DS and all other system components from the packing material and place on a bench or table.
- 3 Inspect the equipment and accessories to ensure there has been no damage during shipment.
- 4 Orient the Agilent 850-DS and all other system components appropriately.

Equipment Ratings / Environmental Requirements for Installation

- Electrical Supply: 100-240 V~ $\pm 10\%$, 50/60 Hz, max. 2 A
- Temperature: 5 °C to 40 °C
- Humidity: max relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40°C
- Altitude: 2000 m
- Pollution Degree: 2
- Installation Category: II
- Indoor use only

Specification

Parameter	Specification
Operating Environment	5 to 40 °C
Humidity (non-condensing)	Not more than 80% RH
Power requirements	Voltage: 90-250 V, 50-60 Hz Current: 2.5-1.0 A
Dimensions	Height: 40.0 cm Height with filter module: 60 cm Width: 39.0 cm Width with filter waste bin: 48.0 cm Depth: 60.0 cm
Weight	27 kg
Interface	7" color touchscreen (English, Simplified Chinese, Spanish, Japanese, German, Hungarian, Portuguese, Russian) or PC controlled
Pumping system and speed	Rotary piston syringe pump, 8 syringes with adjustable 6-12 mL/min.
Sampling accuracy	2.5% *
Sample volume per vial/tube	0.1-14 mL (up to 28 mL with dual sample)
Sample frequency	Minimum of 2 minutes (under certain conditions)
Number of samples per test	Up to 36 (with manual tray exchange after every 12 samples)
Maximum run time	999 hours

* The Sampling Accuracy specification is derived from the following conditions: 10 mL volume, DI H₂O, 10 mL/min pump speed and 3-second Aspiration Dwell Time. This percentage may change based on the sampling parameters or volume, product formulation and type of dissolution media or when in line filters are used.

Equipment, Parts, and Accessories As Shipped

Parts and Accessories

Locate the following items, as applicable, for your system configuration:

Base Unit Components		
Part Number	Quantity	Description
G7930-64000	1	850-DS base unit
K1005-02066	1	Rack for 15mL test tubes
K3090-00155	7	1/4-28 Male to 1/8" barb fitting
3090-0160	8	1/4-28 Male to 1/4-28 Male tubing adapter
33-9019	3	25' of 1/8" silicone tubing
K1003-00232	1	2.5mm Hex key
K1003-00233	1	3mm Hex key
K1003-00234	1	1/4" slotted screwdriver
K1003-00235	1	#1 Phillips screwdriver
1150-7997	1	Stylus for touchscreen
17-5001	1	Box of 100 15mL test tubes
5075-0042	1	DB25 communications cable
77-0001	1	Safety Manual and Documentation CD
Power cord	1	Country-specific power cord

Optional Accessory Kits		
Part Number	Quantity	Description
K1005-05187	1	850-DS printer kit
G7931A	1	850-DS filtration module
K1001-01202	1	Kit for 2, 54-position LC trays, 2mL vials
K1001-01203	1	Kit for 100-position LC tray, 2mL vials
K1001-01204	1	Kit for well plate LC tray, 2 mL x 96 wells
K1005-05212	1	Sample tray for 2mL vials (8 x 12)

Connections

Complete the following sections to connect the necessary tubing and cables for the Agilent 850-DS.




Figure 1. 850-DS Tubing Connections

Tubing

- 1 Connect the three 1/4-28 male to 1/8" barb fittings to the three ports on the left side, nearest the rear of the instrument.
- 2 Using a length of 1/4" silicone tubing, connect the waste fitting to a suitable waste collection container. If media replacement will be used, use a length of 1/4" silicone tubing to connect to a suitable media container. The media must be kept at a suitable temperature by a device provided by the user. If the self-cleaning function will be used, use a length of 1/4" silicone tubing to connect to a suitable rinse media container.

CAUTION

Ensure that the waste connection is not capable of siphoning back into the instrument. The inlet tubing must never fall below the level of waste in the container.

The media bottles should be located at the same bench level as the 850-DS or lower.

- 3 Connect up to 8 sample lines from a dissolution apparatus to the input connections on the left side of the Agilent 850-DS.
- 4 Connect up to 8 return lines from a dissolution apparatus to the output connections on the left side of the Agilent 850-DS. If applicable, connect the male-to-male tubing adapters to the fittings separately; then, attach the return tubing lines. This will prevent possible damage caused by excessive twisting of the tubing.

NOTE

If connecting a dissolution apparatus using less than 8 positions, attach a barb fitting to the inlet and outlet ports of the unused positions. Cut appropriate lengths of silicone tubing and attach them to each barb fitting and route the tubing to an appropriate waste container.

Power Cord

CAUTION

Ensure the correct voltage has been supplied prior to connecting the power cables.

- 1 Attach the provided power cord to the back of the instrument.
- 2 Plug the cord into an electrical outlet of suitable voltage.

Apparatus Control Cable

If connecting to an Agilent, Varian or VanKel dissolution apparatus, connect the DB25 communications cable from the RS232 port on the right rear of the Agilent 850-DS to the 25-pin port on the rear of the dissolution apparatus.

For 850-DS setup with an online UV Dissolution system, refer to the Operator's Manual or the help files in the software for the applicable system.

Synchronized Start

To utilize the Synchronized Start feature of the 850-DS, connect a 9-9 pin null modem cable (5075-0449) from the AUX RS-232 port of one 850-DS to another 850-DS. The desired method may be initiated on both units by pressing **Synchronized Start** on the Start Options screen on one of the 850-DS units. The identical method to be executed must be loaded in the same method slot of each sampling station.

The screenshot shows the 'Start Options' window. It has a blue header bar with the title 'Start Options'. Below the header, there are two main panels. The left panel is titled 'Documentation' and contains four text input fields: 'Operator', 'Lot', 'Batch', and 'Note'. The right panel is titled 'Start Options' and contains three radio buttons: 'Instant' (selected), 'Vessel Temperature', and 'Bath Temperature'. Below these is a checkbox labeled 'Time Delayed' followed by two empty text input fields. At the bottom of the window, there is a button labeled 'SYNCHRONIZED START' with a blue icon, a 'CANCEL' button with a red 'X' icon, and a 'RUN' button with a blue play icon.

Figure 2. Start Options

Initial Power Up

- 1 Make sure there are no obstructions inside the sample tray area of the Agilent 850-DS
- 2 Make sure the door is fully closed.
- 3 Turn on power using the front panel power switch. The unit will execute a self test and initialization.
- 4 Verify the unit passed the initialization and is at the Main screen.

NOTE

Before operating the Agilent 850-DS, a sample collection tray with suitable test tubes, vials or well plates must be inserted. Failure to do so may cause liquids to be pumped into the base of the unit.

Sample Tray Installation

- 1 From the 850-DS main screen, press **Eject** to eject the sample tray platform. The sample tray platform moves to the front of the 850-DS. The movement forward will pause and the unit will beep to alert the user that the tray platform will extend outside the unit.

CAUTION

Do not attempt to insert or remove the sample tray until movement of the platform has completed.

- 2 Load the 12 x 8 sample tray (for 16 x 100 mm test tubes) onto the platform.
- 3 Press **Eject**. The sample tray platform retracts the sample tray to the home position in the 850-DS.

Updating 850-DS Language Screens

The 850-DS firmware can support several languages including:

- English
- Simplified Chinese
- Spanish
- Japanese
- German
- Hungarian
- Portuguese
- Russian

English is the default language setting. To update the language, the file containing the new screens should be placed on an appropriate SD card. Contact your Agilent representative to obtain the necessary files.

Copying files to the SD card

NOTE

You must use a 2GB or smaller SD card formatted as FAT16 or FAT32 to complete this procedure.

Copy the files listed below containing the desired language screens from the PC to the SD card. No additional files may be on the SD card.

- Config.ini
- Flash2.ini
- SLCD-850.BIN

Downloading screens to the 850-DS firmware

- 1 Power off the 850-DS.
- 2 Insert the SD card in the slot at the right front of the instrument and press the card in fully. You may feel a click as the card seats completely.
- 3 Power on the 850-DS and wait approximately 2 minutes for the images to automatically load. When complete, one of the Main screen images will display on the screen.
- 4 Power off the 850-DS, wait 30 seconds, and then power it on again. The update is now complete and the 850-DS should start normally.

Optional Filter Module

NOTE

The 850-DS Filter Module should only be installed by a certified Agilent representative.



Figure 3. Filter Module and Filter Plates

The Agilent 850-DS Dissolution Sampling Station can be optionally equipped with a filtering option that greatly improves dissolution sample processing through the use of specifically designed filter plates. The plates are consistent with filtration membranes and housing materials currently used for dissolution sampling. They are manufactured and sold by GE-Whatman specifically for use with the Agilent 850-DS. Each plate contains eight 25 mm disks incorporated on a flat cartridge. The compact arrangement of the filter disks on the plate also make it much easier for the automated equipment to process and filter samples over traditional automation-ready filters.

The Filter Module, when installed, may be easily disabled if not required for all dissolution methods. Media replacement may also be used with the Filter Module installed. Contact Agilent for assistance with filter validation or visit www.gelifesciences.com for complete filter plate availability and ordering.

Changing Trays / Needle Replacement

The 850-DS can support various sample tray configurations. The standard needle block can support two trays - for test tube or 2 mL vial sample collection. Other alternative sample trays including Agilent HPLC sample trays and 96-well plates are also compatible. Use of these trays requires a modification to the needle block assembly of the 850-DS. This conversion should only be performed by qualified Agilent personnel. These sample tray kits are meant to be installed for extended use and are not designed for frequent exchanges.



4

Operating the 850-DS

Startup Screen 30

Main Screen 32

Menu Screens 36

Startup Screen




Figure 4. Startup Screen

The Startup screen provides initial information after you power on the instrument. The firmware version, LCD version, serial number and release date are displayed while a self-test is executed in the background. If the self-test completes successfully, a green check mark will be displayed for a few seconds and then the instrument will progress to the Main screen. If an error is encountered during the self-test, a red X will be displayed along with a short text description of the component causing the error. After approximately 10 seconds, the Startup screen will change over to the Main screen.

The following errors may be reported: DataFlash, Real Time Clock (RTC), Needle Motor, tray Motor, Syringe Pump Motor, Filter Clamp, and Filter Motor. If you encounter any of these errors, please contact Agilent for service See “Obtaining Warranty and Other Services” on page 99..




Figure 5. Startup Screen with Error

Main Screen

The Main screen is the home screen for all instrument operations. The screen will configure itself to display the appropriate sample format based on selections from the Instrument screen. The Main screen displays the following information:

- instrument model
- menu bar
- system date
- system time
- current user mode
- screen lock state
- leak notification
- alarm notification
- service due notification
- if connected, status of the dissolution apparatus (except when the apparatus is idle)

The sample tray may be ejected from this screen by pressing **Eject**. The tray may be retracted by pressing **Eject** again.

Several configurations are illustrated on the following pages.




Figure 6. Main Screen in Idle Mode




Figure 7. Main Screen with 2 mL Vial, 54-position Blocks Installed




Figure 8. Main Screen with 96-well Plates Installed




Figure 9. Main Screen with 2 mL, 100-position HPLC Tray Installed




Figure 10. Main Screen with BIO-DIS (Apparatus 3)

CAUTION

Ensure the front door is open prior to ejecting the sample tray. Failure to do so may damage the door or cause injury.

Menu Screens

Option	Function
Methods	See "Methods" on page 36..
Alarms	See "Alarms" on page 56..
Calibration	See "Calibration" on page 57..
Cleaning	See "Cleaning" on page 61..
Diagnostics	See "Diagnostics" on page 64..
Instrument	See "Instrument Settings" on page 72..
Mfg Config	See "Manufacturing Configuration" on page 78..
Reports	See "Reports" on page 79..
User Access	See "User Access" on page 80..

Methods

Press **Methods** from the Main screen to create sample-collection methods.

The screenshot shows the 'Method Properties' screen with a grid of 10 numbered boxes. The first two boxes contain the following information:

Box Number	Method Name	Date	Time
1.	ACETAMINOPHEN	01/01/2013	12:30 PM
2.	APP 2 METHOD	01/02/2013	2:30 PM
3.		--/--/----	--:--
4.		--/--/----	--:--
5.		--/--/----	--:--
6.		--/--/----	--:--
7.		--/--/----	--:--
8.		--/--/----	--:--
9.		--/--/----	--:--
10.		--/--/----	--:--

At the bottom right of the screen are two buttons: 'NEXT' (with a right arrow icon) and 'RETURN' (with a checkmark icon).

Figure 11. Method Properties Screen

The first four screens of Method Properties show memory locations 1 through 40. Ten methods are displayed per page and the method name and last save time and date are shown for all non-blank methods. Blank memory slots are indicated by a series of dashes.

Select a storage location for the method to be entered by pressing one of the **numbered memory locations**.

NOTE

Only the memory-location numbers are touchable on screen. The software will not respond if the method details in the box next to each memory location are touched.

Method Properties (Apparatus 1/2/5/6)

NOTE

This screen is applicable to Apparatus 1/2/5/6 only. For Apparatus 3 or 7, advance to **“Method Properties (Apparatus 3/7)”** on page 43.

Method Properties

Name: ACETAMINOPHEN

Vessel Temp: 37.0 °C
Bath Temp: 37.2 °C
Spindle Speed: 50.0 RPM
Media Volume: 900.0 mL
Rotation Start Delay: 3 s

Profile Print Interval: 00:30

☐ Apparatus 1
☒ Apparatus 2
☐ Apparatus 5
☐ Apparatus 6

☐ Initial Temp
☐ Initial Sample
☐ Final Temp

☒ Enable Final Spin
RPM: 250.0
Duration: 1:00

IMPORT CLEAR CANCEL NEXT SAVE

Figure 12. Instrument Parameters

Option	Function
Name	30-character field for the name of this method.
Vessel Temp	Desired vessel temperature for this method.
Bath Temp	Desired bath temperature for this method.
Spindle Speed	Desired spindle speed for this method.
Media Volume	Volume of media initially in the vessels.
Rotation Start Delay	Time after the start of the method before the spindle rotation begins.
Profile Print Interval	Specifies the frequency that dissolution apparatus operational parameters are printed on the optional printer.
Apparatus 1/2/5/6	Select the appropriate apparatus type .
Initial Temp	Takes a temperature reading of the vessels prior to dropping a dosage at the start of a method.
Initial Sample	Takes a sample prior to dropping a dosage at the start of a method.
Final Temp	Takes a temperature reading of the vessels at the end of the method.
Enable Final Spin	After the final timepoint, this option runs the spindles at a high speed in order to try to get all remaining drug into suspension.
RPM	Rotational speed used for the Final Spin.
Duration	Duration the final spin will be executed.
Clear	Erases all values entered and returns them to default values.
Import	Allows user to import method from SD card.
Cancel	Exits the method entry screens and returns to the Main screen without saving any parameters.
Save	Saves all method parameters and returns to the Main screen.

Sample Collection Parameters (Apparatus 1/2/5/6)

Method Properties

Sample Volume mL

Prime Volume mL

Purge Volume mL

Waste Drop Volume mL

Dual Sample ☐

Estimated Min Transfer Time

☐ Media Change

☐ Full Media Change

☐ Media Addition

☒ Enable Cleaning Cycle

☒ PreFill Tubes / Vials

PreFill Volume mL

Figure 13. Sample Collection Parameters - Screen 1

Option	Function
Sample Volume	Volume to be collected into the sample test tubes / vials.
Prime Volume	Volume used to prime all of the sample collection lines prior to collecting a sample.
Purge Volume	Volume used to purge all of the sample collection lines after collecting a sample.
Waste Drop Volume	Volume used to rinse the sample needles prior to delivering each sample. This volume, in addition to the sample volume, should be accounted for if media replacement is enabled.
Dual Sample	When enabled, collects two consecutive samples for each timepoint. Uses the sample volume specified above. This option also halves the total number of timepoints available.
Estimated Min Transfer Time	This is a calculated estimate, based on the previously entered sample collection parameters, of how frequently samples may be taken.
Prefill Tubes / Vials	Use this option to fill each test tube or sample vial with the specified volume (in mL). This volume will be taken from the dissolution vessels and automatically replaced using the Media Replacement channel. The number of rows filled is based on the number of time points of the method. Ensure the media replacement line is connected and sufficient volume is available in the container to maintain initial vessel volume.
Media Change	When enabled, allows you to manually intervene for a media modification.

Option	Function
Full Media Change	Signifies that you will be removing the dosage and performing a complete replacement of media. A dialog box will display to prompt you to begin and to keep track of the total elapsed time. This option will pause the elapsed time of the dissolution method during the change.
Media Addition	Signifies that you will be adding a quantity of media to each vessel. A dialog will display to prompt you to begin as well as to keep track of the total elapsed time. This option will continue the elapsed time of the dissolution method during the change.
Enable Cleaning Cycle	When enabled, executes a cleaning cycle after the test has finished. The parameters for the cleaning routine are taken from the stored values on the Cleaning screen (" Cleaning " on page 61). A validated cleaning method should be performed after each dissolution test. This enables the 850-DS to perform at a high level for an extended period of time while reducing ongoing maintenance costs.
Clear	Erases all values entered and returns them to default values.




Figure 14. Sample Collection Parameters - Screen 2

Option	Function
Enable Filtration	If a Filter Module is installed, this enables filtering samples prior to collection. Ensure an adequate amount of 850-DS filter plates are installed prior to beginning the method.
Samples Per Filter	Specifies how many times to use filters before they are discarded. The filter plates are designed for single use only. The Samples Per Filter setting should be determined for each product and method as part of the system validation. The maximum number of uses per plate is 5 and a warning message will appear for any value greater than 2.
Type	25-character field used to describe what type of filter is being used. Typically porosity and membrane type are described.
Enable Media Replacement	When enabled, pumps a specified amount of media back to the vessels after each timepoint. Media for this step is provided by the user via an external beaker or vessel.
% of Sample	Specifies that the volume to replace will be a percentage of the sample volume.
mL	Specifies that the volume to replace will be this absolute. For complete replacement, the sample + waste drop volume should be entered here.

Method Properties (Apparatus 3/7)

Figure 15. Sample Collection Parameters (Apparatus 3/7)

Option	Function
Sample Volume	Volume to be collected into the sample test tubes / vials.
Prime Volume	Volume used to prime all of the sample collection lines prior to collecting a sample.
Purge Volume	Volume used to purge all of the sample collection lines after collecting a sample.
Waste Drop Volume	Volume used to rinse the sample needles prior to delivering each sample. This volume, in addition to the sample volume, should be accounted for if media replacement is enabled.
Dual Sample	When enabled, collects two consecutive samples for each timepoint. Uses the sample volume specified above. This option also halves the total number of timepoints available.
Enable Filtration	If a Filter Module is installed, this enables filtering samples prior to collection. Ensure an adequate amount of 850-DS filter plates are installed prior to beginning the method.
Samples Per Filter	Specifies how many times to use each filter before they are discarded. The filter plates are designed for single use only. The Samples Per Filter setting should be determined for each product and method as part of the system validation. The maximum number of uses per plate is 5 and a warning message will appear for any value > 2.
Filter Type	25-character field used to describe what type of filter is being used. Typically porosity and membrane type are described.

Option	Function
Enable Cleaning Cycle	When enabled, executes a cleaning cycle after the test has finished. The parameters for the cleaning routine are taken from the stored values on the Cleaning screen (" Cleaning " on page 61). A validated cleaning method should be performed after each dissolution test. This enables the 850-DS to perform at a high level for an extended period of time while reducing ongoing maintenance costs.
Clear	Erases all values entered and returns them to default values.
Save	Saves all method parameters and returns to the Main screen.




Figure 16. Prefill Tubes / Vials Apparatus 3 / 7

Option	Function
Prefill Tubes / Vials	Use this option to fill each test tube or sample vial with the specified volume (in mL). For Apparatus 3 / 7 units, this volume will be taken from the Row specified and automatically replaced using the Media Replacement channel. The number of sample tray rows filled is based on the number of time points of the method. Ensure the media replacement line is connected and sufficient volume is available in the container to maintain the initial vessel volume.

Timepoint Number Selection

Up to 36 timepoints may be specified per method. This screen and the next allow you to define when samples are taken. Timepoints must be consecutively numbered. Entering zero for the time will delete a timepoint.

Method Properties		Sample Time		RPM	
1.	000:15:00	50.0			
2.	000:30:00	50.0			
3.	000:45:00	50.0			
4.	001:00:00	50.0			
5.	---	---			
6.	---	---			
7.	---	---			
8.	---	---			
9.	---	---			
10.	---	---			
11.	---	---			
12.	---	---			
13.	---	---			
14.	---	---			
15.	---	---			
16.	---	---			
17.	---	---			
18.	---	---			

CLEAR CANCEL BACK NEXT SAVE

Figure 17. Timepoint Number Selection (Apparatus 1/2/5/6)

Method Properties		Dip Interval		DPM	
19.	---	---			
20.	---	---			
21.	---	---			
22.	---	---			
23.	---	---			
24.	---	---			
25.	---	---			
26.	---	---			
27.	---	---			
28.	---	---			
29.	---	---			
30.	---	---			
31.	---	---			
32.	---	---			
33.	---	---			
34.	---	---			
35.	---	---			
36.	---	---			

CLEAR CANCEL BACK SAVE

Figure 18. Timepoint Number Selection (Apparatus 3/7)

Option	Function
Clear	Erases all values entered and returns them to default values.
Save	Saves all method parameters and returns to the Main screen.
Cancel	Exits the timepoint entry screen and returns to the Main screen without saving any parameters.

Timepoint Parameter Entry (Apparatus 1/2/5/6)

This screen specifies the time and the parameters to use for each sample collection.

The screenshot shows the 'Timepoint Properties' dialog box. It contains the following fields and controls:

- Timepoint #**: A text field containing the value '1'.
- Sample Time**: A time input field showing '000:15:00' with a unit label 'hh:mm:ss'.
- Spindle Speed**: A numeric input field showing '50.0' with a unit label 'RPM'.
- Next Media Volume**: A numeric input field showing '900' with a unit label 'mL'.
- Media Change**: A checkbox that is currently unchecked.
- Auto Calculate**: A checkbox that is checked, indicated by a green checkmark icon.
- At the bottom right, there are two buttons: 'CANCEL' (with a red X icon) and 'RETURN' (with a green checkmark icon).

Figure 19. Timepoint Properties Screen (Apparatus 1/2/5/6)

Option	Function
Timepoint #	An indicator of the currently selected timepoint.
Timepoint	The amount of time to wait after the method starts before the sample will be collected.
Spindle Speed	The spindle speed to use just after this sample collection finishes.
Next Media Volume	The nominal volume present in the vessels after each sample cycle. This volume is used to determine the proper sample manifold position if used with an Agilent 708-DS. The Auto Calculate option can be enabled to automatically determine this value.
Media Change	Indicates that you will be performing either a full media change or a media addition after this timepoint.
Auto Calculate	Enable this option to have the firmware automatically determine the Next Media Volume based on the method parameters. If enabled, the Next Media Volume field will become inaccessible.

Timepoint Parameter Entry (Apparatus 3/7)

This screen specifies the time and the parameters to use for each sample collection.

Figure 20. Timepoint Properties Screen (Apparatus 3/7)

Option	Function
Interval	The memory slot of the timepoint being configured.
Dip Interval	Duration of dipping time for each row.
Hold Time	Duration the inner sample tubes / sample holders are stopped at the bottom of the stroke before dipping begins.
Drain Time	Duration the inner sample tubes / sample holders drain over the outer media tubes after being lifted from the media.
Row	Enter the desired row number.
Speed	Enter the desired DPM.
Bath Temp	Enter the desired bath temperature.

Run Method Screen 1

Method Properties	
1.	ACETAMINOPHEN 01/01/2013 12:30 PM
2.	APP 2 METHOD 01/01/2013 12:30 PM
3.	----- --/--/---- --:--
4.	----- --/--/---- --:--
5.	----- --/--/---- --:--
6.	----- --/--/---- --:--
7.	----- --/--/---- --:--
8.	----- --/--/---- --:--
9.	----- --/--/---- --:--
10.	----- --/--/---- --:--




Figure 21. Selecting a Method (Apparatus 1/2/5/6)

The first four screens of Method Properties represent memory slots 1 through 40 (ten methods are displayed per page). The method name and last save time and date are shown for all non-blank methods. Empty slots are indicated by a series of dashes. You can use the following onscreen commands to navigate the screens:

Option	Function
Next	Advances to the next page of available methods.
Back	Returns to the previous page.
Return	Returns to the Main screen without any selections made.

Run Method Screen 2

Figure 22. Method Start Parameters

This screen allows you to enter parameters that are stored with the method run. These parameters are available on a printout if your unit has the optional printer installed.

Option	Function
Operator	25-character field for operator name entry.
Lot	25-character field for lot information.
Batch	25-character field for batch information.
Note	25-character field for additional notes.
Dissolution Apparatus Control	Allows the Agilent 850-DS to control a dissolution apparatus. Note: Enable this option only if the Agilent 850-DS is connected to another Agilent/Varian/VanKel dissolution apparatus. If connected to a third-party dissolution apparatus with resident probes, do not check this box.
Instant	Start a method instantly.
Vessel Temperature	Will start a method once the vessels have achieved a stable temperature as specified in the method to be run within the temperature tolerance settings of the 850-DS.
Bath Temperature	Will start a method once the bath has achieved a stable temperature as specified in the method to be run within the temperature tolerance settings of the 850-DS.

Option	Function
Time Delayed	Will start the method at the time and date specified. Note: Time-delayed method starts are not allowed with the Synchronous Start option.
Synchronized Start	Allows two 850-DS collectors to execute the same method simultaneously. Each 850-DS must have the same method loaded in the same memory location. Methods may be transferred between 850-DS units using an SD card (See "850-DS Method Transfer (to another 850-DS)" on page 54.). A 9-9 pin null modem cable (5075-0449) connects one instrument to the other using the AUX RS-232 port.

NOTE

Time-delayed method starts are not permitted with the Synchronized Start option.

Run Method Screen 3




Figure 23. Run Method Screen 3

Option	Function
RPM	Current spindle speed.
°C	Current bath temperature.
Elapsed Time	Current elapsed time since the start of the method.
Next Sample	Time to the next sample collection time.
Remaining Time	Time required to complete the current method.
View Samples	Graphical indication of the method progress. The shaded circles indicate filled test tubes or vials.
View Timepoints	Provides a view of the timepoints in the current running method.
Stop	Press STOP to display the dialog screen confirming the choice to stop or continue the method.
Eject	Press Eject to access the sample tray. The sample tray may be ejected during an active program as long as the system is not taking a sample. Press RETURN (on the screen displayed) to retract the tray prior to the next sample time point; otherwise, the sample will not be taken.

Run Method Screen 4




Figure 24. Run Method Screen 4

Run Method Screen 5

This screen displays at the conclusion of a method when no errors have been encountered. The total elapsed time, the current time, and the current date are displayed.




Figure 25. Run Method Screen 5

850-DS Method Transfer (to another 850-DS)

NOTE

You must use a 2GB or smaller SD card formatted as FAT16 or FAT32 to complete this procedure.

- 1 Insert an SD card into the SD card slot on the back of the 850-DS.
- 2 Select the **Reports** tab from the 850-DS main screen.




Figure 26. Reports Screen

- 3 Press **SD Card** in the "Print Method..." row. Allow the method to save to the SD card.
- 4 Eject the SD card.
- 5 On 850-DS where you want to import the method, insert the SD card on the back of the unit.

- 6 Select the **Methods** tab from the 850-DS main screen.
- 7 Select an **open location** where the imported method should be stored (press **number**).
- 8 Press **IMPORT**. Allow method to upload from SD card.

Method Properties

Name:

Vessel Temp: °C
Bath Temp: °C
Spindle Speed: RPM
Media Volume: mL
Rotation Start Delay: s

Profile Print Interval:

☐ Apparatus 1
☒ Apparatus 2
☐ Apparatus 5
☐ Apparatus 6

☐ Initial Temp
☐ Initial Sample
☐ Final Temp

☒ Enable Final Spin

RPM:
Duration:

IMPORT **CLEAR** **CANCEL** **NEXT** **SAVE**

Figure 27. Import Method

Alarms

Press **Alarms** from the Main screen to configure alarms on the 850-DS.

Figure 28. Alarm Screen

Option	Function
Timer Alarm	This alarm acts in a manner similar to a laboratory countdown timer. Once a value is entered and you press RETURN , the alarm will start a countdown to zero. You do not need to remain on this screen for the alarm to remain active. Once the timer value reaches zero, an alarm will sound and a dialog box will display to indicate that the alarm has expired.
Pause Duration Alarm	This option sets the duration that you are allowed while performing media changes and additions. The duration of a pause is recorded and printed on the optional printer.
Delayed Start	Beeps prior to starting a method when you have selected Delayed Start .
Moving	Beeps prior to any tray movements of the Agilent 850-DS.
Sample Complete	Beeps at the end of each sample timepoint.
Program Complete	Beeps at the conclusion of a method.

Calibration

Press **Calibration** from the Main screen to calibrate the 850-DS.

Calibration

Volume Calibration

mL to dispense 

Row Number 

Actual mL

mL / stroke 

Full Range Volume Dispense

Test Tube Rack ONLY. Requires test tubes in rows 1 through 7. Will pump 0.1, 0.5, 1, 2, 4, 10 and 15 mL into the test tubes.

Notifications

Calibration Due Date

PM Due Date

Figure 29. Calibration Screen

Volume Calibration

The Agilent 850-DS uses a syringe pump to deliver accurate volumes from 0.1 mL to 14.0 mL. A volume calibration function is available to compensate for pump delivery speed, tubing length, and media viscosity. Volume accuracy is a function of pump speed and aspiration dwell time. To improve volume accuracy, try using slower pump speeds. The pump speed and aspiration dwell time are set on the second Instrument Settings screen. See "Instrument Settings Screen 2" on page 74..

NOTE

Agilent recommends calibrating using a volume of 10 mL. If that does not provide acceptable performance, you may perform a calibration at the typical sample collection volume. Although the Agilent 850-DS can dispense volumes as low as 0.1 mL, it is not recommended to calibrate with volumes less than 2 mL.

Option		Description
Volume Calibration	mL to Dispense	Enter volume (in mL) to be collected in sample tubes / vials during the volume calibration procedure.
	Row Number	Specify the row where the sample volume should be collected.
	Actual mL	Enter the average sample volume collected (in mL) during the volume calibration procedure.
	mL / Stroke	Automatically calculated value based on volume calibration parameters. This value is stored until the next volume calibration is performed. No user entry is required/permitted.
	Prime	Primes (fills) the sample tubing lines prior to executing the volume calibration procedure. This ensures the lines are sufficiently filled and an accurate sample volume is collected.
	Dispense	Dispenses the programmed volume (in mL) into the specified row. Pressing Dispense displays an onscreen dialog that guides the user through the volume calibration process.
		Volume Calibration Prior to performing this volume calibration, ensure that: 1. Sample cannulas are sufficiently submerged. 2. Clean, dry test tubes have been weighed and placed in the specified row. 3. The sample lines have been primed. 4. The volume to be collected does not exceed the sample container capacity. CANCEL RETURN
		Volume Calibration To complete the volume calibration: 1. Record the gross weight of each sample. 2. Calculate the net weight of each sample. 3. Calculate each sample volume using the density of the media. 4. Calculate the average volume and enter it in the Actual mL field. Operator RETURN
Stop		Stops the volume calibration procedure.

Figure 30. Volume Calibration Pre-check

Figure 31. Volume Calibration Completion

Full Range Volume Dispense

This function allows the user to verify the full range accuracy of the pump and is only supported for 15 mL test tubes. Begin with a tray containing seven rows of clean, dry test tubes whose initial weight has been taken. Pressing **Run** will start a pumping cycle where 0.1 mL, 0.5 mL, 1.0 mL, 2.0 mL, 4.0 mL, 10 mL, and 15 mL will be dispensed into the corresponding first seven rows. When complete, gravimetric analysis may be done as above to determine the actual volumes delivered.

Option		Description
Full Range Volume Dispense	Run	Starts the Full Range Volume Dispense cycle.
	Stop	Stops the cycle immediately.
Notifications	Calibration Due Date	Enter the date that instrument calibration is due. On the due date, a reminder notification (wrench icon) will display. The wrench icon will remain visible on the Main screen as long as the current date is past the calibration date. Enter a future date to remove the reminder icon.
	Preventative Maintenance Due Date	Enter the date that preventative maintenance is due. On the due date, a reminder notification (wrench icon) will display. The wrench icon will remain visible on the Main screen as long as the current date is past the PM date. Enter a future date to remove the reminder icon.
History		This screen displays the details of the previous 20 volume calibrations performed on the instrument including date, time, pump speed, and operator. The complete log of calibrations may be printed pressing PRINT on this screen.





Figure 32. Volume Calibration History

Cleaning

Press **Cleaning** on the Main screen to access the cleaning cycle on the 850-DS.



The screenshot shows the 'Cleaning' screen with a blue header. It contains several input fields and buttons. Under 'Cleaning Cycle', 'Volume' is set to 5.0 mL and '# Cycles to Run' is set to 1. There are 'RUN' and 'STOP' buttons. Under 'Source', 'Rinse Port' is selected with a radio button, and 'Sample Cannula' is unselected. Under 'Reservoirs', 'Rinse Media Replacement' and 'Rinse Waste' are both unselected. There is a '# Cycles to Run' field set to 1. A 'RETURN' button with a green checkmark is at the bottom right.

Figure 33. Cleaning Parameters Screen

NOTE

A validated cleaning method should be performed after each dissolution test. This enables the 850-DS to perform at a high level for an extended period of time while reducing ongoing maintenance costs. Cleaning methods are formulation dependant. If you need help determining appropriate cleaning solution and cycles, contact the Dissolution Hotline at dissolution.hotline@agilent.com.

The Agilent 850-DS can perform a manual or an automated flush of all of the tubing. This screen allows a manual cleaning cycle to be started. If Enable Cleaning Cycle is selected in a method, the parameters entered on this screen are used to control the cleaning cycle.

Option		Description
Cleaning Cycle	Volume	Enter the volume that will be pumped through the tubing.
	# Cycles to Run	This controls how many times the above volume will be repeated.
	Run	Starts a cleaning cycle with the entered parameters.
	Stop	Stops a cleaning cycle.
Source	Rinse Port	This selects the rinse port as the source of the rinse media. Connect an appropriate container to the Agilent 850-DS rinse port with a length of 1/4 inch tubing. At the start of the cleaning cycle, an internal pump fills the rinse reservoir on the tray assembly. The needles are lowered into the rinse media and the specified volume is pumped from the needles through the sample tubing and then through the return tubing. All lines will then be purged with air to eliminate as much liquid from the tubing as possible.
	Sample Cannula	Selects the sample cannula as the inlet source. At the start of a cleaning cycle, the syringe pump moves the volume specified from the sample cannula through the return tubing. A smaller fixed volume of 0.7 mL will be pumped through the needles into the waste row. All lines will then be purged with air to eliminate as much liquid from the tubing as possible. It is the responsibility of the user to properly submerge the sample cannulas in the appropriate rinse solution prior to this cycle.
Reservoirs	Rinse Media Replacement	Enable this option to rinse the internal media replacement reservoir after each test with rinse solution. The reservoir will be filled and emptied based on the number of cycles defined.
	Rinse Waste	Enable this option to rinse the internal waste reservoir after each test with rinse solution. The reservoir will be filled and emptied based on the number of cycles defined.
	# Cycles to Run	This value defines how many times the appropriate reservoirs will be filled and emptied with rinse solution at the conclusion of each test. Rinsing of these reservoirs serves to prepare the system for the next test as well as flush the internal volume sensors to prevent long term damage.

CAUTION

Ensure that the rinse connection is not capable of siphoning back into the instrument. The inlet tubing must never fall below the level of rinse in the container.

The media bottles should be located at the same bench level as the 850-DS or lower.

Diagnostics

Press **Diagnostics** on the Main screen to access diagnostics on the 850-DS.




Figure 34. Diagnostics Screen 1

This screen allows manual control of all the sample handling parameters. Tray movement, needle position, and syringe pump control are available. The tray image on the left side will change to indicate which tray type is currently configured. As the tray is moved to valid rows, the appropriate row in the tray image will be filled to provide feedback for row position.

Option		Description
Row Position	Row Indicator	Pressing the numeric row indicator will display the numeric entry screen. Enter the desired row number and press RETURN .
	Up Arrow	Increments the current row number by one.
	Down Arrow	Decrements the current row number by one.
	Move	Moves the tray to the position indicated in the row indicator field.
Needle Position	Up Arrow	Raises the needles.
	Down Arrow	Lowers the needles.
Tray Position	Eject	Moves the tray to the eject position.
	Home	Performs a home of the needles, tray, and pump in order.
	Waste	Moves the tray to the waste row of the waste block.
	Media	Moves the tray to the media-replacement row of the waste block.
	Rinse	Moves the tray to the rinse-row of the waste block.
	Current Position	An indicator that shows the current position of the tray. If the tray has been moved to a particular row number, this indicator will be blank and the tray position will be indicated on the tray image.
Pump Control	Volume	Press the numeric indicator and enter the volume to be dispensed.
	Fill	Aspirates the requested volume into the syringe pump.
	Dispense	Pumps the contents of the syringe pump to the selected location.
	To: Needles	The requested volume is aspirated from the sample cannula and pumped out to the needles.
	To: Sample Return	The requested volume is aspirated from the sample cannula and pumped out to the return lines.

Diagnostics Screen 2




Figure 35. Diagnostics Screen 2

Option	Description	
Rinse Pump	Pump In	Activates the pump to fill the rinse row.
	Pump Out	Activates the pump to empty the rinse row.
	Stop	Stops the pump.
Media Replacement Pump	Pump In	Activates the pump to fill the media-replacement row.
	Pump Out	Activates the pump to empty the media-replacement row.
	Stop	Stops the pump.
Waste Pump	Run	Activates the waste pump (empties the waste row).
	Stop	Stops the pump.
Self Test Results	Run	Runs the self-test routine.
	Print	Prints the results of the self-test routine.
Printer Test	Print	Conducts a test to ensure the printer is functioning properly.

Diagnostics Screen 3




Figure 36. Diagnostics Screen 3

This screen allows manual control of the optional Filter Module (if installed). The filter may be moved to one of four positions and the filter clamp may be opened and closed.

Option		Description
Filter Transport	Load	Moves the filter transport plate to the load position. This is the location of the filters that are placed in the inlet stack.
	Media Repl	Moves the filter transport plate to the media replacement position. This position allows the filter clamp plates to clamp on themselves. This position is not directly under the new filters in order to minimize any carryover or contamination between a used filter and the new ones in the stack. This position is used when a media replacement is being performed.
	Filter	Moves the filter transport plate to the filter position. This position allows the filter clamp plates to seal on the filter and pump samples through it.
	Eject	Moves the filter transport plate to the eject position.
	# Filters to Eject	Pressing the # Filters to Eject indicator will display the numeric entry screen. Enter the number and press RETURN .
	Cycle Filters	Will cycle the filter transport plate between the load and eject position for # Filters to Eject times. This may be used to empty the inlet stack or verify operation.
Filter Clamp	Open	Opens the filter clamp plates. The clamp mechanism will stop automatically when it reaches its open position.
	Close	Closes the filter clamp plates. Caution: Care should be used with this function as damage may occur if there are any obstructions in the clamping path.
	Stop	Turns off the clamp motor. The plates remain in their current position.
	Clamp Force Control	This parameter adjusts the amount of power the motor receives from 90 to 100%. Values below this result in poor clamping and are thus not allowed.

Diagnostic Screen 4 (Apparatus 1/2/5/6)




Figure 37. Diagnostics Screen 4 (Apparatus 1/2/5/6)

Option		Description
Spindle	Setpoint RPM	Pressing the RPM indicator will display the numeric entry screen. Enter the desired RPM and press RETURN .
	Run	Starts the spindles.
	Stop	Stops the spindles.
Dosage Drop Module (DDM)	Position	Pressing the Position indicator will display the numeric entry screen. Enter the desired number and press RETURN .
	Up Arrow	Increments the DDM position.
	Down Arrow	Decrements the DDM position.
	Single	Actuates the individual DDM selected by the position indicator.
	Open All	Actuates all DDM modules simultaneously.
Bath	Setpoint	Pressing the Setpoint indicator will display the numeric entry screen. Enter the desired temperature and press RETURN . The dissolution apparatus bath setpoint will be updated to the value entered.
	Actual	Displays the current actual bath temperature.
Vessel Temperatures		If a manifold is installed and lowered, the current vessel temperatures are displayed here.

Option		Description
Drive	Up	Raises the drive unit
	Down	Lowens the drive unit.
Manifold	Up	Raises the manifold.
	Down	Lowens the manifold.

Diagnostic Screen 4 (Apparatus 3/7)





Figure 38. Diagnostics Screen 4 (Apparatus 3/7)

Option		Description
Sample Head	Home	Returns the drive unit to the home position.
	Row	Enter the desired row number.
	Move	Moves the drive unit to the entered row.
Sample	Speed	Enter the desired speed in DPM.
	Dip Hold Time	Stops the inner sample tubes / sample holders at the bottom of the stroke for the length of time entered.
	Run	Begins manual dipping.
	Stop	Stops manual dipping.
Bath	Setpoint	Pressing the Setpoint indicator will display the numeric entry screen. Enter the desired temperature and press Return. The dissolution apparatus bath setpoint will be updated to the value entered.
	Actual	Displays the current actual bath temperature.

Instrument Settings

The following screens show the parameters that influence the setup and operation of the Agilent 850-DS. These parameters apply to all methods and modes of operation.



The screenshot displays the 'Instrument Settings' screen. At the top, a blue header bar contains the title 'Instrument Settings'. Below this, the screen is organized into several sections. The 'Collector' section at the top contains two columns of fields: 'Serial Number' (DL00000D000), 'LAN MAC Address' (a dashed line), 'Firmware Version' (1.00), 'LCD Version' (1.00), 'Release Date' (01 JAN 2000), 'Collector Name' (LAB001), 'Collector Comm ID' (01), 'Prime Loss Volume' (4.0 mL), and 'Number Active Vessels' (6). Below the 'Collector' section are three main settings areas. The 'Dissolution Apparatus' section has two radio buttons: 'Apparatus 1/2/5/6' (selected) and 'Apparatus 3/7', with a 'Comm ID' field set to 01. The 'Sample Format' section has four radio buttons: 'Test Tube' (selected), 'Vial Blocks', 'Well Plate', and '100 Vial Tray'. The 'Vessel Assignments' section features a graphic of a test tube with blue liquid. At the bottom right, there are two buttons: 'NEXT' with a green arrow icon and 'RETURN' with a green checkmark icon.

Collector	
Serial Number	DL00000D000
LAN MAC Address	-----
Firmware Version	1.00
LCD Version	1.00
Release Date	01 JAN 2000
Collector Name	LAB001
Collector Comm ID	01
Prime Loss Volume	4.0 mL
Number Active Vessels	6

Dissolution Apparatus
☒ Apparatus 1/2/5/6
☐ Apparatus 3/7
Comm ID: 01

Sample Format
☒ Test Tube
☐ Vial Blocks
☐ Well Plate
☐ 100 Vial Tray

Vessel Assignments


NEXT **RETURN**

Figure 39. Instrument Settings Screen 1

Option		Description
Collector	Serial Number	Manufacturer-assigned serial number.
	LAN MAC Address	MAC address that will be used for LAN communications.
	Firmware Version	Version of the firmware installed on the instrument.
	Release Date	Release date of the firmware.
	Collector Name	A 25-character field to describe the instrument.
	Collector Comm ID	A 2-digit numeric field that is used when a PC or other instrument is communicating with the Agilent 850-DS. Note: If multiple instruments are connected together, each unit must have a unique Comm ID.
	Prime Loss Volume	The volume in mL that is pumped at the beginning of each sample point. This volume represents the volume in the tubing between the sample collector and the needles. Pumping this volume to the return lines first ensures that the volume aspirated at the sample time is what is actually delivered to the sample tubes / vials instead of the volume remaining in the lines after prime.
Dissolution Apparatus	Number Active Vessels	Enter the number of active vessel positions on the dissolution apparatus. This value is important if the vessel temperature start is used.
	Apparatus 1/2/5/6	Sets the operating mode to expect an Apparatus 1/2/5/6 instrument to be connected to the Agilent 850-DS.
	Apparatus 3/7	Sets the operating mode to expect an Apparatus 3/7 instrument to be connected to the Agilent 850-DS.
	Comm ID	Communications ID of the attached Agilent dissolution apparatus.
Sample Format	Test Tubes	Indicates a rack holding an array of 8 x 12 test tubes or vials is installed.
	Vial Blocks	Indicates a rack holding two 9 x 6 blocks with 2 mL vials is installed.
	Well Plate	Indicates a rack holding two 96-position well plates is installed.
	100-vial Tray	Indicates a rack holding an array of 10 x 10 (offset spacing) for 2 mL vials is installed.

	1	2	3	4	5	6	7	8
Sample	☒	☒	☒	☒	☒	☒	☒	☒
Control	☐	☐	☐	☐	☐	☐	☐	☐
Blank	☐	☐	☐	☐	☐	☐	☐	☐
Not Used	☐	☐	☐	☐	☐	☐	☐	☐

RETURN

Figure 40. Vessel Assignment Screen

Vessel Assignments

This screen allows you to enter a category for each vessel's usage. These categories are used for reporting only.

Instrument Settings Screen 2

Instrument Settings

Tolerances

Temperature Tolerance

Speed Tolerance

Leak Detect Options

☐ Continue Method

☐ Stop Immediately

Timepoint Auto-Increment

Minutes

LCD Settings

☐ ☐

☐ ☐

☐ ☐

Pump Speed

☐ 12 mL/min

☐ 10 mL/min

☐ 8 mL/min

☐ 6 mL/min

Align Screen

Date / Time

Aspiration Dwell Time

seconds

BACK RETURN

Figure 41. Instrument Settings Screen 2

Option		Description
Tolerances	Temperature Tolerance	Displays the numeric entry screen. Enter the desired tolerance and press Return. Temperatures that are outside of this tolerance to their respective setpoints will raise an error condition. These errors will be reported through dialog screens and on the optional printer.
	Speed Tolerance	Displays the numeric entry screen. Enter the desired tolerance and press Return. Spindle speeds that are outside of this tolerance to their respective setpoints will raise an error condition. These errors will be reported through dialog screens and on the optional printer.
Leak Detect Options	Continue Method	If a leak is detected, a notification will be shown on the Main screen, but otherwise all running methods will continue uninterrupted.
	Stop Immediately	If a leak is detected a notification will be shown and any currently running method will be aborted.
Timepoint Auto-increment		This value controls the auto populate feature of method timepoint entry. After the first timepoint is entered, all subsequent timepoints will have this value added to them.
LCD Settings	Brightness Up	Increases the brightness of the LCD backlight.
	Brightness Down	Decreases the brightness of the LCD backlight.
	Speaker Up	Increases the volume of the internal speaker.
	Speaker Down	Decreases the volume of the internal speaker.
Pump Speed		Sets the speed of the syringe pump to the selected rate. Use slower speeds for higher viscosity media or media containing surfactant. Higher pump speeds may be used for less viscous media. The pump speed may also be reduced to improve volume accuracy.
Aspiration Dwell Time		Sets the time (in seconds) that the syringe pauses after pulling in liquid prior to dispensing. This setting may need to be increased to maintain volume accuracy for certain types of dissolution media. The Dwell Time is automatically adjusted (3-second minimum) when filtration is enabled.




Figure 42. Touch Screen Calibration Screen

Option	Description
Align Screen	Starts the touchscreen alignment routine. Follow the instructions on the screen to calibrate the touchscreen. If the screen is not touched for 10 seconds, the Agilent 850-DS will return to the Instrument Settings screen.

Option	Description
Date / Time	Pressing Data / Time will display the Date and Time setting screen.




Figure 43. Date and Time Change Screen

Option	Description
Operator	30-character field for the operator name. Any change of the date or time requires an operator's name. The name is used in the audit log as described below.
Time	Time Field Pressing the time will display the numeric entry screen. Enter the desired time and press Return.
	Time Format Select either 12- or 24-hour time format to be used in all displays and printouts.
Date	Date Field Pressing the date will display the numeric entry screen. Enter the desired date and press Return.
	Date Format Select one of the three date formats to be used in all displays and printouts.
Audit Log	The audit log shows the last ten time or date changes and who made them. The up and down arrows scroll through the list.
Return	Returns to the Instrument Settings screen.

Manufacturing Configuration

NOTE

The Mfg. Config. menu can only be accessed by a certified Agilent representative. The options enabled or disabled on this screen control functions visible to the user when operating the 850-DS.

The Manufacturing Configuration screen features a blue header bar with the title "Manufacturing Configuration". Below the header, there is a "Serial Number" label followed by a text input field. The screen is divided into three main sections: "Instrument Options" on the left, "Test Mode" on the right, and "Reset Administrator Access Password" at the bottom. The "Instrument Options" section contains three checkboxes labeled "Printer", "Filter Changer", and "Username / Password Access". The "Test Mode" section contains a yellow warning triangle icon and a checkbox. The "Reset Administrator Access Password" section contains a "Reset" button with a circular arrow icon. A "RETURN" button with a green checkmark icon is located in the bottom right corner of the screen.

Figure 44. Mfg Config Screen

Reports

Press **Reports** from the Main screen to print or export reports on the 850-DS. Reports can be printed on the optional printer, exported through the auxiliary serial port (AUX RS-232), or saved to an SD card inserted in the back SD card slot.




Figure 45. Reports

Option	Function
Last Run Results	Prints the last run results. Select the destination by pressing the appropriate button.
Print Method	Displays the select method screen. Once a method has been selected, prints the method parameters. Select the destination by pressing the appropriate button.
Print Instrument Settings	Prints the instrument settings. Select the destination by pressing the appropriate button.
Enable Printer	A toggle function that enables or disables the printer. If the printer is not enabled, all printing functions are suppressed.
Return	Returns to the Main screen.

NOTE

In order to print reports to the printer, the optional printer must be enabled. All other functions are available with the Enable Printer function deselected.

User Access

Press **User Access** from the Main screen to assign the system access level, modify functionality of the User Levels and, if enabled, regulate individual user access. The 850-DS firmware supports a system-based level of control or access to specific functions based on individual users. This distinction is governed by the Username/Password Access feature (see User Access screen on [page 81](#)).

If the Username/Password Access feature is enabled, individual users are assigned a Name, Password, and User Level and must log in to access the 850-DS firmware. Specific functionality is based on their User Level which is defined by an Administrator of the system.

If the Username/Password Access feature is disabled, an Administrator may toggle between the three available access levels available on the 850-DS. The current access level is indicated by a green circle. Any user may change the access level to a lower setting (e.g., Advanced User to User) without entering a password. Any changes to a higher level (e.g., Advanced User to Administrator) require a password.




Figure 46. User Access

Option	Function
Change Password	Allows an administrator to change the password. Displays the password change screen. The new password must be entered correctly twice before it will be changed. The default password is 12345.
User Level Setting	This screen allows an Administrator to define access to specific functions of the 850-DS based on the User Group assigned to the system or individual user. The system settings are stored when Save is pressed.
Users List	This screen allows an Administrator to create, modify and define access for users of the 850-DS. A total of 20 individual users may be stored. Each user is assigned a User Level (Administrator, Advanced User, or User) that determines what functions may be accessed. An individual Name and Password is assigned to each user at the time of creation.

User Level Setting			
Function	1	2	3
Method Settings	✓		
Alarm Settings	✓	✓	✓
System Clean Settings	✓	✓	
Diagnostics	✓	✓	
Instrument Settings	✓	✓	✓
Calibrations Settings	✓	✓	
Run Methods	✓	✓	✓
Report Settings	✓	✓	✓

1. Administrator, 2. Advanced User, 3. User

Figure 47. User Level Setting

Users List			
1.	 Admin	6.	
2.	 User1	7.	
3.	 Operator1	8.	
4.	 Operator2	9.	
5.	 Operator3	10.	


Figure 48. Users List

Dialog Screens (Blue)




Password Reset

Displays after the Administrator Access password has been reset to factory default.




Memory Access

Displays when memory is being accessed. Typically it will be displayed when saving or retrieving methods.




Waiting on Heating - Vessel Temperature Start

Displays when a method is running and the 850-DS is waiting on the dissolution apparatus for temperatures to reach a setpoint or become stable. The current temperature, the target temperature, and the vessel temperatures are displayed.




Waiting on Heating - Bath Temperature Start

Displays when a method is running and the 850-DS is waiting on the dissolution apparatus for temperatures to reach a setpoint or become stable. The elapsed time, current temperature, and the bath temperature are displayed.




Time Delayed Start

Displays when a time delayed method start has been selected. The expected start time is displayed so you know how long it will be until the method begins.




Reading Temperatures

Displays during any temperature-read events on an Agilent dissolution apparatus. Temperature-read events typically last about 40 seconds.




Method Stop Confirmation

Displays after you have pressed STOP while a method is running. It allows you to confirm your intentions to stop or continue running the method.




Change Media

Notifies you that it is time to add or change media during a method. The elapsed time is shown in the dialog. If the Pause Duration Alarm ("**Pause Duration Alarm**" on page 56) time entered in the Alarms screen is exceeded, then an error will be generated and noted in the method log.




Access Control

Used when changing access levels or when locking the display. It serves as a generic password entry screen. You must enter access codes correctly twice before any changes will be accepted.



Cleaning Cycle Complete


Displays when the cleaning cycle is complete.



The User Login screen has a blue header bar with the text "User Login". Below the header, there are two input fields: "Name" and "Password". At the bottom right, there are two buttons: "CANCEL" with a red X icon and "LOGIN" with a green checkmark icon.

User Login

Displays when Lock is pressed from the 850-DS Main screen if the Username/Password Access feature is enabled. The individual user must enter their Name and Password to gain access to their predefined functions.




The Configure User Authority screen has a blue header bar with the text "Configure User Authority". Below the header, there are two input fields: "Name" and "Password". Below these fields, there are three user level icons: a single person (blue), a group of three people (green), and a group of three people with a plus sign (red). Below the icons are three radio buttons. At the bottom left is a "DELETE" button with a trash can icon. At the bottom right are "CANCEL" (red X) and "SAVE" (green checkmark) buttons.

Configure User Authority

Displays when a number is pressed from the Users List screen. This screen is used to create, modify or delete users. Each user is assigned a Name, Password, and User Level from this screen.


Warning Screens (Yellow)



The Apparatus 1/6 Auto Head warning screen has a yellow header bar with the text "Apparatus 1/6". The main text area contains the instructions: "Please raise head and attach dosage(s). Press RETURN to continue." At the bottom right is a "RETURN" button with a green checkmark icon.

Apparatus 1 / 6 Auto Head


Displayed for USP Apparatus 1 or 6 that have an automated drive-unit lift. You are prompted to attach the appropriate dosage form. The head is automatically lowered after you press RETURN.



The Apparatus 1/6 Manual Head warning screen has a yellow header bar with the text "Apparatus 1/6". The main text area contains the instructions: "Please raise head and attach dosage(s). Lower head and press RETURN to start test." At the bottom right is a "RETURN" button with a green checkmark icon.


Apparatus 1 / 6 Manual Head

Displayed for USP Apparatus 1 or 6 methods on units that have a manual drive-unit lift. You are prompted to attach the appropriate dosage form and lower the drive unit manually. The method will resume after you press RETURN.




Manual Dosage Drop

Displayed for USP Apparatus 2 or 5 methods on units that do not have the dosage delivery modules (DDMs) installed. After you insert the dosage and press RETURN, the method will resume.




Head Raised

Generic message on units that have a manual head lift. At any time the head is not fully lowered, this prompt will display.




Sample Tray Full

Used on methods that have more timepoints than will fit on one sample rack. After you replace the rack with a fresh one and press RETURN, the method will continue. If the tray is not retracted prior to the next time point, sample collection will be skipped.




Sample Tray Ejected

Displayed if the tray is ejected during an program. Press **RETURN** to retract the tray prior to the next time point. If the tray is not retracted prior to the next time point, sample collection will be skipped.




709-DS Preheat Special Condition

Displays when connected to an Agilent 709-DS running an Apparatus 1/6 method. In order to achieve reliable start conditions, finned baskets must be used during pre-heat. Press **RETURN** to start the pre-heat cycle.




Alarm Time Expired

Notifies you that the alarm time has expired.




Clear Method Confirmation

Confirmation that all data from a method will be erased and replaced with default values. This prevents you from accidentally erasing methods.



Service Due


Displays when either the Calibration Date Preventative Maintenance Date ("**Calibration**" on page 57) has passed. This dialog appears once. After this screen has been cleared, the service due wrench icon appears in the lower right corner of the Main screen. The service due icon will remain visible until the corresponding date has been updated in the Calibration Screen.



User Delete Confirmation


Displays when an attempt is made to delete a user from the Users List using the Configure User Authority screen. Only an Administrator has the ability to delete users from the 850-DS.

Error Screens (Red)




Method Error

This screen appears at the conclusion of a method when errors have been encountered. The total elapsed time, the current time, and the current date are displayed.



Invalid Data

Displayed when data has been entered that is outside of the acceptable limits. Press **RETURN** to re-enter the data. The prompts in the data entry screens will inform you of the valid data ranges.




Method Not Supported

Displayed when the loaded method requires hardware or options that are not installed. For example, if a method requires filtration and no Filter Module is installed.




Printer Out of Paper

This dialog is displayed when the printer is out of paper. The paper status is only checked at the beginning of a print request, not continuously during a print. See "Printer Maintenance" on page 93..




Timepoints Too Close

This dialog is displayed when entering timepoints and the minimum transfer time is not met. Increase the timepoint spacing or reduce the sampling volume parameters to meet the minimum transfer time parameter.




Not Valid Tray Position

This dialog is only displayed when in Diagnostics. It is displayed if the tray is not in a valid position when a pumping action has been requested (e.g. dispensing a volume through the needles, but the tray is in the home or eject position). Moving the tray to a valid row will resolve this issue.




Clock Audit Log

Displayed if a change to either the time or date has been requested and no operator has been entered. Operator names are required for clock changes and for saving into the audit log.




Filter Clamp Plates Open

This dialog is displayed on units with a Filter Module installed. If a pumping action has been requested and the filter clamp plates are in the open position, this error will appear. Usually, this is encountered during diagnostics. Navigate to the Filter Module section of the Diagnostics screen (See "Diagnostics Screen 2" on page 66.) and close the clamp plates. Once closed, normal pumping is allowed.




Filter Plate Warning

Displayed if a value > 2 is entered for the number of samples to filter prior to filter plate replacement. Saturation may occur and volume inaccuracy or system performance could deteriorate under certain conditions. Proper filter validation should be performed to ensure acceptable performance.



Filter Transport Conflict

This dialog is displayed on units with a Filter Module installed. If in the Diagnostics Screen, a filter transport move is requested and the clamp plates are closed, this error will be displayed. The transport can't move while the clamp plates are closed. Navigate to the Filter Module section of the Diagnostics screen (See "Diagnostics Screen 2" on page 66.) and open the clamp plates. Once opened, normal filter transport movement is allowed.





5

Maintenance and Troubleshooting

Preventive Maintenance 92

Printer Maintenance 93

Obtaining Warranty and Other Services 99

Preventive Maintenance

Preventive maintenance intervals may vary depending on the frequency of instrument usage.

Clean System

The Agilent 850-DS can perform an automated flush of all the tubing. Refer to **“Cleaning”** on page 61 for instructions on how to run the cleaning cycle. A validated cleaning method should be performed after each dissolution test. This enables the 850-DS to perform at a high level for an extended period of time while reducing ongoing maintenance costs.

Periodic Maintenance

During a visual inspection of unit, check for:

- Loose or missing hardware
- Signs of wear on pulleys and belts
- Tubing blockages
- Operation of needles
- Signs of wear in the syringes

Printer Maintenance

Removing a Thermal Paper Roll



Figure 49. Near-empty Paper Roll

NOTE

A red line on the side of the paper roll indicates the paper supply is almost exhausted. Replacement is recommended when the red line is visible. If the paper roll is not changed and the paper supply is exhausted, the 850-DS gets an out-of-paper signal from the printer, and displays an error on the screen.

- 1 Place fingers underneath the light-colored handle and pull forward to pop out the front of the printer.



Figure 50. Accessing the Printer Handle

Opening the chamber reveals the paper roll or empty core inside.




Figure 51. Open Printer Chamber

- 2 Internal side levers hold the roll in place. Since it is empty or nearly empty, use your fingers to pull out the small roll or core.

NOTE

The empty core is disposable. Each new roll includes a cardboard core. New rolls may be ordered using part number 12-9988 (10 rolls) or 12-9989 (100 rolls).




Figure 52. Removing the Thermal Paper Roll

Inserting a Thermal Paper Roll

NOTE

Only a 58 mm wide x 25 meter (maximum) thermal printer paper roll with a 13 mm core size can be used in the panel-mount printer.

- 1 With the paper chamber open, insert the paper roll until the levers snap the roll into place. Load the paper so the leading edge of the paper feeds from over the roll and not from under it.




Figure 53. Inserting a Thermal New Paper Roll

- 2 Lead the paper out for a few inches and close the top lid of the chamber.




Figure 54. Leading Paper and Closing the Lid

- 3 Keeping the paper centered, close the bottom lid of the chamber and snap it into place.




Figure 55. Closing the Printer Chamber

- 4 To verify that the paper is not skewed or jammed, press the **upper-right button** on the printer to feed some paper. If paper does not feed, re-open the chamber and re-center the paper.



Figure 56. Feeding Paper

Obtaining Warranty and Other Services

To place a service order (warranty or other services), please contact your local Customer Care Center. Contact information can be found at www.agilent.com under your country using the Contact Us link. Place your service request using the displayed phone number or E-mail address.

In This Book

- Chapter 1 Safety
- Chapter 2 Introduction
- Chapter 3 Setting Up the 850-DS
- Chapter 4 Operating the 850-DS
- Chapter 5 Maintenance and Troubleshooting

www.agilent.com

© Agilent Technologies, Inc. 2018

70-9074 Rev D March 2018

