



Agilent 7000 Series Triple Quadrupole GC/MS System

Data Sheet



The Agilent 7000 Series GC/MS Triple Quadrupole System with MassHunter software offers outstanding sensitivity and selectivity. The GC and MS components have been designed to deliver consistent sensitivity and accurate quantification of trace-level target compounds in very complex matrices. The Agilent 7000 Series system must be combined with the high-performance 7890A gas chromatograph.

Agilent 7000 Series Triple Quadrupole GC/MS System

Triple Quadrupole Mass Spectrometer

Mode (standard)	EI (High Sensitivity Extraction Source)
Modes (optional)	PCI and NCI
Ion source material	Noncoated, proprietary inert source
Ion source temperature	106 to 350 °C
Filaments	Dual filaments for EI
Electron energy	10 to 300 eV
Mass range (<i>m/z</i>)	10 to 1,050
Resolution	Unit mass adjustable by tune, 0.7 to 2.5 daltons
Dynamic range (electronic)	> 10 ⁶
Scan rate (electronic)	Up to 6,250 u/s
MRM speed (transitions/sec)	500
Minimum MRM dwell	1 msec
Mass filters (2)	Proprietary monolithic hyperbolic gold-coated quadrupole
Mass axis stability	< ± 0.10 u over 24 hours (10-40 °C)
Quadrupole temperature	106 to 200 °C
Collision cell	Linear hexapole
Collision cell gas	Nitrogen with helium quench gas for reduction of metastable helium
Collision energy	Selectable up to 60 eV
Detector	Triple-Axis HED-EM with extended-life EM and dynamically ramped-iris
Tuning	Autotune or manual
Total gas flow	Up to 8 mL/min GC carrier plus another 5 mL/min of methane for CI operation, plus an additional 1–2 mL/min of N ₂ and He for the collision cell gases
Pumping system	Dual stage turbomolecular pump
Software	Agilent MassHunter acquisition, data handling (quant/qual) and reporting
Simultaneous MS and GC	Can collect 2 GC detector signals while acquiring MS data



Agilent Technologies

Gas Chromatograph (Agilent 7890A GC)

For more specifications on GCs refer to the GC Data Sheet

Injector	Split/splitless, Multi-mode inlet, PTV and others available
Autosampler	7693 ALS, CombiPAL, G1888A ¹ Headspace Sampler
Oven temperature	Ambient 4 to 450 °C
Oven ramps/plateaus	20/21. Negative ramps are allowed.
Electronic pneumatic control (EPC)	Auto pressure regulation for split/splitless, septum purge
Carrier gas control modes	Constant pressure and flow modes; pressure and flow programmable
Pneumatic splitter	Capillary Flow Technology devices for effluent splitting, backflushing, and column switching

Installation Checkout Specifications

EI MRM sensitivity	100 fg of octafluoronaphthalene (OFN) will produce > 500:1 RMS S/N for the transition of m/z 272→241 or 272→222 using autotune parameters
PCI MRM sensitivity	100 fg BZP will produce > 50:1 RMS S/N for the transition of m/z 183→105 (using methane)

Other Reference Specifications²

EI MRM area precision	10 sequential injections of 100 fg OFN will give an area precision of < 8% RSD
EI scan sensitivity	1pg OFN will produce > 300:1 RMS S/N for m/z 272 scanning from m/z 50 to 300
NCI SIM sensitivity	10 fg OFN will produce > 200:1 RMS S/N for m/z 272 (using methane)

Physical Requirements

Dimensions (Triple Quad MS)	35 cm (w) × 86 cm (d) × 47 cm (h) Additional space should be added for the data system and printer.
Weight (Triple Quad MS)	59 kg
Dimensions (mechanical pump)	18 cm (w) × 35 cm (d) × 28 cm (h)
Weight (mechanical pump)	21.5 kg
Dimensions (7890A GC)	58 cm (w) × 54 cm (d) × 57 cm (h)
Weight (7890A GC)	45 kg

1. G1888A programmed through module keyboard and operated with Remote Start.

2. EI scan and NCI SIM function are demonstrated at installation; reference specifications are not confirmed.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

Dwell (msec)	10	5	1
# Transitions	18	33	99
Cycle Time	198	198	198
MRM/sec	91	167	500
n = 1	37453	39629	36246
2	37854	38444	34524
3	36464	39402	35007
4	38547	37270	33895
5	38372	38376	36955
6	36349	37881	33315
7	39439	32900	34861
8	37654	37787	37986
9	37588	38620	34307
10	37140	36425	34700
Avg Area	37686	37673	35180
SD Area	940	1925	1447
% RSD	2.50%	5.10%	4.10%

Typical peak area precision at short dwell times and high MRM/sec conditions (sample: 10 pg propylamide)

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2009
Printed in the USA
October 29, 2009
5989-9625EN