



# MAC26

*Autosampler Device*



## Microtrac MAC26 *Features*

- Holds up to 26 organic or inorganic samples in the carousel
- Constructed out of surgical-grade stainless steel
- FLEX software provides fully programmable measurement schedules, pre-defined SOP's, and manual control allowing users to define a "home" or "dump and rinse" vial
- Runs dry powder or liquid dispersions
- "Tap and Rinse"

Size:	30.48W X 30.48D X 31.75H cm (12W X 12D X 12.5H in)
Number of Samples:	Up to 26
Sample Volume:	Up to 1.8 ml wet or dry
Sample Vial:	Glass, compatible with most organic and inorganic solvents
Construction:	All external parts are constructed of stainless steel or painted with acid resistant polyurethane paint for corrosion protection. Sample carousel is constructed with stainless steel.
Power:	115 or 230 V AC, 50/60 Hz.

*For more information on the unique benefits of the MAC26, a software demo, or to request a quote, please scan this QR code with your smartphone or tablet or visit [microtrac.com/mac26](http://microtrac.com/mac26)*



## Systems/Accessories

S3500



BLUEWAVE



SDC



*Complex Particles. Clear Data.*

**SL-PS-06 Rev. D**