


New Reagents for Dako CoverStainer.

Choose the H&E staining intensity you want.




New flexible staining intensities to fit y


Dako H&E Staining




Skin




Lymph nodes




Blood vessel




Kidney




Intestine




Colon adenocarcinoma



Intestine



Skin



Kidney

Proven consistency


Consistently stain up to 3,000 slides with the same reagents for accurate H&E results that you can trust, slide after slide. The Dako Reagent Management System consistently circulates reagents and optimizes reagent consumption, which prevents precipitation. Dako Ready-to-Use Reagents streamline your workflow and come in closed containers to ensure that your reagents last longer. A unique rack design, only available for Dako CoverStainer, eliminates cross contamination risk and minimizes reagent carryover. H&E consistency made easy.

Product	Description	Order No.
Dako Hematoxylin	Ready-to-Use	CS700
Dako Eosin	Ready-to-Use	CS701
Dako Bluing Buffer	Ready-to-Use	CS702




our preferences.


Dako Progressive Staining




Skin




Lymph nodes




Blood vessel




Kidney




Intestine




Colon adenocarcinoma



Intestine



Skin



Kidney


Exceptional staining quality

Obtain crisp and clear staining with the Dako Progressive Staining protocol. Dako Harris Hematoxylin and Dako Modified Eosin Y show great morphology detail across different tissue types even from a low magnification, allowing quick and easy assessment of slides and shorter time to diagnosis.


Product	Description	Order No.
Dako Harris Hematoxylin	Ready-to-Use	CS709
Dako Modified Eosin Y	Ready-to-Use	CS711
Dako Bluing Buffer	Ready-to-Use	CS702




Dako Regressive Staining




Skin




Lymph nodes




Blood vessel




Kidney




Intestine




Colon adenocarcinoma



Intestine



Skin



Kidney

Excellent differentiation

Dako Gill's 3 provides a color intense hematoxylin staining showing well-defined nucleus and blue goblet cells. Used with Dako Differentiation Solution and Dako Eosin Y Phloxine B, your staining will show beautiful contrast with blue nucleic appearance and very distinct eosinophilic components varying from pink to orange, giving clear and bright details with an excellent differentiation of cell structures.

Product	Description	Order No.
Dako Gill's 3 Hematoxylin	Ready-to-Use	CS708
Dako Eosin Y Phloxine B	Ready-to-Use	CS710
Dako Differentiation Solution	Ready-to-Use	CS707
Dako Bluing Buffer	Ready-to-Use	CS702




Every lab is different.


Different pathologists have different H&E staining intensity preferences. Two new sets of reagents for Dako CoverStainer provide you with the staining intensity that your lab prefers, combined with the unmatched proven consistency and high-quality staining of Dako CoverStainer.

Three validated and customizable protocols ensure that your lab's individual needs are met, with ease.


H&E is the most widely used diagnostic test in the pathology lab, so consistently stained slides at the intensity you prefer are a necessity for fast and easy reading.



Slide 1



Slide 3,000



Optimize your workflow with Dako CoverStainer and provide completed slides to the pathologist more quickly and continuously throughout the day.

Dako CoverStainer.

The efficient, safe and accurate choice for H&E

Efficient

Dako Coverstainer fully automates every step of the H&E process from baking to drying, reducing manual processes and freeing up staff to work on other things. A high throughput of up to 240 slides per hour maximizes productivity, and continuous batch loading/unloading will even out your workload and contribute to an easy and efficient operation.

Safe

Dako CoverStainer enables you to improve operator safety by reducing exposure to hazardous reagents and with quiet operation of only 38 decibels, providing a working environment that is free from unhealthy fumes and noise. It's also safe for the environment, with less waste produced and disposed of in closed containers.

Accurate

The Dako CoverStainer solution enables your lab staff to stain with consistent high quality and provide the right diagnosis every single time.

Trusted Answers. Together.



Dako

Agilent Pathology Solutions

www.agilent.com

Represented in more
than 100 countries

Australia
+61 1800 802 402

Austria
+43 1 408 43 34 0

Belgium
+32 16/93 00 30

Brazil
+55 11 50708300

Canada
+1 800 387 8257

China
+80 08 20 3278

Denmark
+45 44 85 97 56

Finland
+358 9 348 73 950

France
+33 1 64 53 61 44

Germany
+49 40 69 69 470

Ireland
+353 1 479 0568

Italy
+39 02 58 078 1

Japan
+81 3 5232 9970

Korea
+82 80 004 5090

The Netherlands
+31 20 42 11 100

Norway
+47 23 14 05 40

Poland
+48 58 661 1879

Spain
+34 93 344 57 77

Sweden
+46 8 556 20 600

Switzerland
+41 41 760 11 66

United Kingdom
+44 (0)1 353 66 99 11

United States of America
+1 805 566 6655